PATH NOSTICS MOVING CARE FORWARD The Development of an Integrated Molecular and Phenotypic Method for the Management of Urinary Tract Infections

Dave Baunoch, PHD

Urinary Tract Infections





Cost of Urinary Tract Infections (UTI)

Cost to the Healthcare System

Responsible for **≈10.5** MILLION office visits/year¹ UTI complications result in DAYS longer for each hospital stay²

Urinary Tract Infections account for

50%

of all Medicare admissions⁵ **Problem: Urinary Tract** Infections cost the US healthcare system

\$13**B**

a year. Without efficacious diagnostic tests, ineffective antibiotics will continue to lead to more resistant microorganisms.

Cost to Humanity

Up to 1/3 of infections illustrate resistance to an antibiotic³

According to the CDC, antibiotic **resistance** gives rise to at least





- Sheng, W. et al. 2005. Impact of Nosocomial Infections on Medical Costs, Hospital Stay, and Outcome in Hospitalized Patients. J 2. Formos Med Assoc. 104 (2005), 318-326
- Schmiemann, G. et al. 2012. Resistance profiles of urinary tract infections in general practice -- an observational study.BMC urology, 12, 1 (2012), 33.
- http://www.cdc.gov/drugresistance/threat-report-2013/
- 5. Strausbaugh, L.J. et al. The Burden of Infection in Long-Term Care 2016 Infect Con and Hosp Epid 2000:21:674-679





Impact of Urosepsis

50% Of the sepsis cases in Nursing Homes originate from urinary tract infections 60%

Of patients with sepsis caused by urinary tract infections die² \$2.8 B

Total cost of sepsis due to urinary tract infections in the US annually³ Problem: Culture has significant limitations including the inability to detect gram-positive organisms. 57.1 %

of patients with gram-positive urosepsis die.

Impact of an incorrect Diagnosis and Treatment



67.8%

Of patients prescribed the wrong antibiotic on initial treatment died (as compared to 28.7% of patients receiving the right antibiotic)⁴



Initial treatment delay increases the odds of admission to ICU up to **1.98**

And decreases the overall survival rate⁵





The Problem with Routine Culture and Sensitivity





Problem Statement #1

URINE IS NOT STERILE AND CONTAINS A MICROBIOME

- Urine has historically been viewed to be sterile with the presence of bacteria in the urine seen as an indicator of infection.
- Recent work has clearly demonstrated that the urinary tract contains a microbiome that is characterized by a preponderance of *Lactobacillales* in women and *Corynebacterium*in in men.¹ Often times, the bladder contains additional uropathogenic organisms (including E.coli) that can coexist in a persistent state of asymptomatic bacteriuria.²



Traditional Method

Step 1: Identify Organism: 24-48 hours Step 2: Susceptibility: 24-48 hours

Characteristics of bacterial colonies can help in the process of identification











Problem Statement #2

URINE CULTURE MISSES THE VAST MAJORITY OF POTENTIAL UROPATHOGENS

- **Current Methodology** Developed in the 1950's, the standard method involves applying 1ul of urine onto Blood and MacConkey Agar plates and incubating them at 35 degrees centigrade for 24 hours in the presence of oxygen.¹
- What is the Fundamental Problem with culture as a detection method? From the composition of the agar, to the pH, gas ratio's, and time of incubation culture is a methodology that has been biased for the detection of a subset of pathogens primarily *E. Coli.*²
 - The biased results developed using this methodology often creates findings that are not consistent with the clinical symptoms
 - It is unable to detect slow growing organisms including fastidious and anaerobic organisms as well as most gram-positive organisms
 - Time consuming process that can take up to 72 hours to complete
 - Inherent methodology issues limit the number of organisms reported to no more then 2 with 3 or more considered indications of contamination

Uropathogens detected Standard urine culture (SUC) versus enhanced quantitative urine culture (EQUC)

Organism	SUC	EQUC
Actinobaculum schaalii	0	5
Aerococcus urinae	1	11
Alloscardovia omnicolens	0	6
Candida albicans	0	2
Candida parapsilosis	0	4
Citrobacter koseri	0	1
Corynebacterium riegelii	0	3
Corynebacterium urealyticum	0	1
Enterobacter aerogenes	1	3
Enterococcus faecalis	1	12
Escherichia coli	22	24
Klebsiella pneumoniae	3	8
Proteus mirabilis	0	1
Pseudomonas aeruginosa	1	1
Serratia marcescens	0	1
Staphylococcus aureus	1	5
Streptococcus agalactiae	1	10
Streptococcus anginosus	1	25



Blood agar, 1 µL, 24 hours, aerobic



Adapted from June 9, 2016 ASM Microbe presentation by Alan J. Wolfe, PhD: "Urine Is Not Sterile: Why We Should Care."

Blood agar, 100 µL, 48 hours, 5% CO2

PATH NOSTICS

Problem Statement #3

THE ORGANISMS MISSED BY URINE CULTURE ARE CLINICALLY SIGNIFICANT

The Loyola Study- Pivotal study demonstrating limitations of standard culture in identifying the pathogens associated with UTI's.

- Followed 150 patients who were split into two groups based on whether they believed they were symptomatic for UTI.¹
- Compared the results obtained when they used standard culture to an enhanced version which had modified growth conditions including an increased incubation time.
- In the group who believed they were symptomatic standard culture detected only 57% of the uropathogens where the enhanced methodology detected 91%.



Why is this important?

IN SYMPTOMATIC PATIENTS, STANDARD CULTURE MISSED 2/3 OF ALL POSITIVE PATIENTS

- Outcome
 - Culture positive: 38.5%
 - Failed to see symptoms improve with standard of care based upon culture results
 - Symptoms worsen
 - Culture negative 67%
 - Failed to see symptoms improve
 - Symptoms worsen





Problem Statement #4

URINE CULTURE MISSES POLYMICROBIAL INFECTIONS

- Wolfe and Brubaker have proposed moving from an E.coli-centric view of urinary tract infections.¹
 - With an increasing number of studies demonstrating that most urinary tract infections have multiple urinary pathogens present in the same sample, we should begin to shift our thinking away from a monocentric view of urinary tract infections.
- Studies by Price et al, 2016, and Wolfe et al, 2012 have demonstrated that a significant number of cases identified as having only E.coli present actually have multiple species present.^{2,3}
 - In the Price et al, 2016 paper 35 of 43 patients had additional species
 - In 25 of 35 cases the additional species were uropathogens
- A simple truth: It is becoming increasingly clear that a significant number of urinary tract infections are polymicrobial in nature. Because of the polymicrobial nature of infections, efficacy of treatment is dropping significantly.

1. Price, T.K., Hilt, E.E., Dune, T.J. et al. Int Urogynecol J (2018) 29: 205. https://doi.org/10.1007/s00192-017-3528-8



^{2.} Wolfe AJ, Toh E, Shibata N, Rong R, Kenton K, Fitzgerald M, Mueller ER, Schreckenberger P, Dong Q, Nelson DE, Brubaker L. 2012. Evidence of uncultivated bacteria in the adult female bladder. J. Clin. Microbiol. 50:1376–1383.

Price, T.K., Dune, T., Hilt, E.E., Thomas-White, K.J., Kliethermes, S., Brincat, C., Brubaker, L., Wolfe, A.J., Mueller, E.R. and Schreckenberger, P.C. 2016. The Clinical Urine Culture: Enhanced Techniques Improve Detection of Clinically Relevant Microorganisms. Journal of clinical microbiology. 54, 5 (May 2016), 1216–22

Problem Statement #5

INTERACTION OF BACTERIA IN POLYMICROBIAL INFECTIONS RESULTS IN CHANGES TO ANTIBIOTIC SENSITIVITY

The sharing of metabolic products provides polymicrobial infections an advantage¹

- In the presence of antibiotics the sharing of metabolic products plays a protective role increasing resistance and virulence
- Testing isolated bacteria may over or under estimate the degree of antibiotic resistance



Gram-Positive Uropathogens, Polymicrobial Urinary Tract Infection, and the Emerging Microbiota of the Urinary Tract

> KIMBERLY A. KLINE¹ and AMANDA L. LEWIS² Singapore Centre on Environmental Life Sciences: Engineering, School of Biological Sciences.

Naryang Technological University, Sugapore 637551; "Department of Molecular Microbiology, Washington University School of Medicine, Sr. Loeis, MD 63110

Increased human pathogenic potential of Escherichia coli from polymicrobial urinary tract infections in comparison to isolates from monomicrobial culture samples

Gemma Croxall,¹ Vivienne Weston,² Susan Joseph,¹ Georgina Manning,¹ Phil Cheetham¹ and Alan McNally¹

Pathogen Research Group, Nottingham Trent University, Clifton Lane, Nottingham NG11 BNS, UK Nottingham University Hospitals, Nottingham, UK

The current dispracise standard processions sufficient by the Health Protection Agency for uningtract infections (UTIII) in crinical aboratories does not report bacteria isolated from samples containing these or more upmortain and the same have (UTII) surported and untrasted, particularly in defined patients, where polymicinala UTII mamples are especially present. This subst apports the presence of the major upmoltgravity and upmoltants. The hypothese isolations of mestance to polymicinala UTII samples are especially highlights that Exclusivity and upmolt the interview of the sample same satisfically more invasive (P<2000) in a vie explained and interview polymicinala UTII samples are satisficially more invasive (P<2000) in a vie explained and interview polymicinala UTII samples are satisficially more invasive (P<2000) in a vie explained and interview polymicinala UTII samples are satisficially more invasive (P<2000) in a vie explained and interview polymicinala UTII samples are satisficially more invasive provident UTII samples in more in the satisfic damon more interview polymicinala UTII samples are satisficially more invasive provident UTII samples in more in more in the satisfic damon more interview polymore in the satisfic damon more interview polymore polymore





THE LENGTHY TURNAROUND TIME FOR STANDARD CULTURE AND SENSITIVITY SUPPORTS EMPIRIC TREATMENT AND ANTIBIOTIC RESISTANCE



2018-PUUS Western IRB Number: 20181661 Not yet published.

Summary of the Problem

ROUTINE CULTURE AND SENSITIVITY

- **1.** Urine is not sterile and contains a Microbiome
- 2. Urine culture misses that vast majority of potential uropathogens
- 3. The organisms missed by urine culture are clinically significant
- 4. Urine culture misses polymicrobial infections
- 5. Interaction of bacteria in polymicrobial infections results in changes to antibiotic sensitivity
- 6. The lengthy turnaround time for standard culture and sensitivity supports empiric treatment and antibiotic resistance



The Development of an Integrated Molecular and Phenotypic Method for the Management of Urinary Tract Infections







A unique rapid molecular test for both pathogen identification and antibiotic sensitivity providing personalized therapy options that work the first time.



ORGANISMS DETECTED:

BACTERIAL/YEAST ORGANISMS

- Acinetobacter baumannii
- Actinobaculum schaalii
- Aerococcus urinae
- Alloscardovia omnicolens
- Candida albicans
- Candida glabrata
- Candida parapsilosis
- Citrobacter freundii
- Citrobacter koseri

BACTERIAL GROUPS

- Coagulase neg. staphylococci*
- Viridans group streptococci**

KEY:

- Guidance Basic and Guidance Comprehensive
- Guidance Comprehensive only
- Add-on tests
- Corynebacterium riegelii
- Enterobacter aerogenes
- Enterococcus faecalis
- Escherichia coli
- Klebsiella oxytoca
- Klebsiella pneumoniae
- Morganella morganii
- Mycoplasma hominis

SEXUALLY

Mycoplasma genitalium

Chlamydia trachomatis

Neisseria gonorrhoeae

Trichomonas vaginalis

TRANSMITTED ORGANISMS

- Mycobacterium tuberculosis
- Pantoea agglomerans
- Proteus mirabilis
- Providencia stuartii
- Pseudomonas aeruginosa
- Serratia marcescens
- Staphylococcus aureus
- Streptococcus agalactiae
- Ureaplasma urealyticum



Emerging Organisms

	Gram Negative											
Citrobacter freundii	Known to cause several nosocomial infections of the urinary tract (5). Wiping posterior to anterior following a bowel movement and normal intercourse can both bring about C. freundii UTIs.											
Citrobacter koseri	Can cause UTI's more frequently in individuals with urinary diversions.											
Enterobacter aerogenes*	E. aerogenes is acquired by unintentional bacteria transmission in hospital settings, and E. aerogenes infections mainly occur when host immunity defenses are already suppressed.											
Escherichia coli	Uropathogenic E. coli (UPEC) is one of the leading root causes of UTIs.											
Klebsiella oxytoca	Cause UTI's in individuals with catheters.											
Klebsiella pneumoniae	UTIs caused by multidrug-resistant K. pneumoniae isolates are a significant public health issue.											
Morganella morganii	UTIs are typically connected with complicated features like indwelling catheters.											
Pantoea agglomerans	Is associated with catheter-related bacteremia.											
Proteus mirabilis	Can cause symptomatic infections of the urinary tract, including cystitis and pyelonephritis, and can appear in cases of asymptomatic bacteriuria.											
Providencia stuartii	Are among the most common cause of catheter-associated UTIs, particularly in the elderly with long-term indwelling urinary catheters.											
Pseudomonas aeruginosa	Is a common cause of UTIs, and patients with impaired immunity possess a greater risk for colonization by this organism.											
Serratia marcescens	Has now been implicated as an etiological agent in a wide range of infections, including UTIs.											

Gram Positive											
Actinobaculum schaalii	Is a frequent cause for urinary tract infections amongst the elderly, in children, and in those with underlying urological conditions.										
Aerococcus urinae	Can cause simple and complicated UTIs, bacteremia, and endocarditis in elderly adults who have multimorbidity, chronic urinary retention, or indwelling catheters.										
Alloscardovia omnicolens	Is hard to identify and one of the organisms identified by Enhanced Quantitative Urine Culture in women complaining of UTI symptoms.										
Corynebacterium riegelii	Have been reported to cause opportunistic infections in both immunocompromised as well as immunocompetent patients.										
Enterococcus faecalis	Can trigger endocarditis and septicemia, UTI, meningitis, and other infections in people.										
Staphylococcus aureus	Has been identified in children with urinary tract abnormalities and vesico-ureteral reflux.										
Streptococcus agalactiae	Has been isolated in patients with UTI symptoms.										





Real-time PCR workflow





2-10mL urine sample



Samples undergo extraction





DNA isolated and loaded



Real-time PCR "CHIP"

- Ability to run high numbers of samples simultaneously, with large numbers of targets.
- Isolated DNA and Master Mix is plated onto a nanofluidic plate with 3,072 throughholes.
- Each chip can hold 21-24 samples (18 patient samples and 3 controls) 56 Targets in duplicate
- Estimated Assay Run Time: 2 hours





Software analysis



A. Raw data is gathered through the software for analysis.

- Amplification plots are displayed for individual targets, samples, and subarrays.
- CT's, Amplification scores, and Cq Confidence values are given for each subarray and its 64 wells.







GENOTYPE ANSWERS ONLY PART OF THE COMPLEX PROBLEM OF • Guidance tests for the ANTIBIOTIC RESISTENCE

- Antibiotic Class Antibiotic Class Gene Gene ErmA + Erm B VIM КРС Macrolide resistance ermC mefA IMP-2 group TFM IMP-1 group OXA-23 CTX-M group 1 IMP-16 SHV VEB IMP-7 Carbapenem resistance OXA-1 **OXA-72** CTX-M group 2 **OXA-40** Extended-Spectrum-Betalactamase CTX-M group 9 **OXA-58** CTX-M group 8/25 **OXA-48** PER-1 NDM PER-2 blaOXA-48 GES TetM Tetracycline TetB blaNDM-1 QnrA aaC6-aph3 Quinolone and fluoroquinolone resistance Aminoglycoside QnrB anti-la-aph2 Trimethoprim/Sulfamet drf(A1, A5), sul (1,2) Methicillin resistance mecA hoxazole ampC, FOX, ACC vanA1 AmpC resistance vanA2 Vancomycin resistance DHA, MOX/CMY, BIL/LAT/CMY vanB
- presence of 38 genes know to be associated with resistance to certain antibiotics
- Does Not Provide the Complete Answer – Why?
 - Limited number of resistance genes that can be identified via molecular assay
 - Gene resistance continuously change

Pathnostics Moving Care Forward

 Resistance gene may not be active.

Developing a Methodology for Measuring Pooled Sensitivity







It is becoming increasingly clear that a significant number of urinary tract infections are polymicrobial in nature. Because of the polymicrobial nature of infections, efficacy of treatment is dropping significantly.



Gram-Positive Uropathogens, Polymicrobial Urinary Tract Infection, and the Emerging Microbiota of the Urinary Tract

KIMBERLY A. KLINE¹ and AMANDA L. LEWIS²

¹Singapore Centre on Environmental Life Sciences Engineering, School of Biological Sciences, Nanyang Technological University, Singapore 637551; ²Department of Molecular Microbiology, Washington University School of Medicine, St. Louis, MO 63110 Increased human pathogenic potential of *Escherichia coli* from polymicrobial urinary tract infections in comparison to isolates from monomicrobial culture samples

Gemma Croxall, 1 Vivienne Weston, 2 Susan Joseph, 1 Georgina Manning, 1 Phil Cheetham 1 and Alan McNally 1

¹Pathogen Research Group, Nottingham Trent University, Clifton Lane, Nottingham NG11 8NS, UK ²Nottingham University Hospitals, Nottingham, UK

The current diagnostic standard procedure outlined by the Health Protection Agency for urinary tract infections (UTIs) in clinical laboratories does not report bacteria isolated from samples containing three or more different bacterial species. As a result many UTIs go unreported and untreated, particularly in elderly patients, where polymicrobial UTI samples are especially prevalent. This study reports the presence of the major uropathogenic species in mixed culture unine samples from elderly patients, and of resistance to font-line antibiotics, with potentially increased levels of resistance to ciprofloxacin and trimethoprim. Most importantly, the study highlights that *Escherichia coli* present in polymicrobial UTI samples are statistically more invasive (P<0.001) in *in vitro* epithelial cell infection assays than those isolated from monomicrobial culture asmples. In summary, the results of this study suggest that the current diagnostic standard procedure for polymicrobial UTI samples may pose an increased risk to human health.

Polymicrobial Infections are in the Community



PATHNO

MOVING CARE FORWARD

Resistance Higher For Polymicrobial Infections





Pathnostics



Polymicrobial Interactions Change MIC Levels

Increased human pathogenic potential of Escherichia coli from polymicrobial urinary tract infections in comparison to isolates from monomicrobial culture samples

Gemma Croxall,¹ Vivienne Weston,² Susan Joseph,¹ Georgina Manning,¹ Phil Cheetham¹ and Alan McNally¹

¹Pathogen Research Group, Notlingham Trent University, Clifton Lane, Notlingham NG11 8NS, UK ²Notlingham University Hospitala, Notlingham, UK

The current diagnostic standard procedure outlined by the Health Protection Agency for urinary tract infections (UTIs) in clinical laboratories does not report bacteria isolated from samples containing these or more different bacterial species. As a result many UTIs go unreported and untreated, particularly in elderly patients, where polymicrobial UTI samples are especially prevalent. This study reports the presence of the major unopathogenic species in mixed culture unive samples from elderly patients, and of resistance to front-ine antibiotics, with potentially increased levels of resistance to ciprofloxacin and trimethoprim. Most importantly, the study highlights that Escherichia coli present in polymicrobial UTI samples are statistically more invasive (P<0.001) in in vitro epithelial cell inflection assays than those isolated from monomicrobial outure samples. In summary, the results of this study suggest that the current diagnostic standard procedure for polymicrobial UTI samples may pose mis present in polymicrobial UTI samples may pose an increased risk to human health.

Table 1. Prevalence of antibiotic resistance in UTI isolates

'Poly' refers to strains of polymicrobial infection origin, whilst 'mono' refers to strains of monomicrobial infection origin. -, Particular species/ antibiotic combination not tested.

Antibiotic concentration (µg ml ⁻¹)	Escherich	Escherichia coli (%)		Eacherichia coli (%)		coccus is (%)	Protesta (9	minabilis 6)	Staphyl	ососсыя я (%)	Paradomonas aeraginosa (%)		
	Poly (n=129)	Mone (n=21)	Poly (n=110)	Mone (s=4)	Poly (s=56)	Mone (s=1)	Poly (s=18)	Mone (n=4)	Poly (s=46)	Mone (s=5)			
Gentamicin (2)	12.4	4.76	-	-	5.3	0	4.5	0	-	-			
Cefotaxime (1)	17.8	14.29	-	-	33.9	100.0	-	-	-	-			
Ceftazidime (1)	18.6	9.52	-	-	35.7	100.0	-	-	-	-			
Meropenem (2)	0	•	15.5	•	0	•	-	-	-	-			
Piperacillin- tazobactarn (16)	6.2	4.76	45	25.0	23.2	100.0	-	-	-	-			
Co-amosiclav (32)	5.4	•	2.7	0	16.1	100.0	0	0	-	-			
Trimethopeim (2)	44.2	28.57	-	-	89.0	•	22.2	0	-	-			
Ciprofloxacin (4)	23.3	9.52	28.2	•	0	•	55.5	50.0	-	-			
Cefradine (32)	28.7	19.05	-	-	55.3	100.0	50.0	25.0	-	-			
Nitrofurantoin (32)	17.1	9.52	10.9	•			11.0	0	-	-			
Amoxicillin (32)	45.0	42.86	7.3	•	37.5	100.0	66.6	75.0	-	-			
Gentamicin (10)	-	-	-	-	-	-	-	-	2.17	0			
Piperacillin- tazobactara (85)	-	-	-	-	-	-	-	-	0	0			
Ceftazidime (30)	-	-	-	-	-	-	-	-	2.17	0			
Meropenem (10)	-	-	-	-	-	-	-	-	2.17	0			
Ciprofloxacin (1)	-	-	-	-	-	-	-	-	4.35	0			

Mutualism: Cross Feeding







Antibiotic Sensitivity Testing

Upon completion of UTI Molecular Testing, results are used to set up Antibiotic Sensitivity Testing for each sample.



30

GUIDANCE" U	TI EVE: RECURRENT, PERSISTENT, OR CON	Ocusetions, Please call 716-966-1221 or vielt pathnottics.com MPLICATED UTI REPORT + TB/STI
Patient: First Last	DOB: 10-16-1959	Case#: PUXRI9-010846-TS
POOLED PHENOTYPIC SENSITIVITY	DETECTED (\$):	
Ciprofloxacin (PO/IV) Levofloxacin (PO/IV)	 Nitrofurantoin (PO) Piperacillin/Tazobactam (IV) 	 Tetracycline (PO)
POOLED PHENOTYPIC RESISTANCE	DETECTED (R):	
Ampicillin (PO/IV) Ampicillin/Sulbactam (IV) Cefeclor (PO) Cefeclol (V) Cefepime (IV)	- Cefoxitin (IV) - Ceftazidime (IV) - Ceftriaxone (IM/IV) - Fosfomycin (PO)	- Gentamicin (IM/IV) - Meropenem (IV) - Sulfamethoxazole/Trimethoprim (PO/IV) - Vancomycin (IV)
RESISTANCE GENE(S) DETECTED (RGD):	
Ampicillin Resistance ESBL Resistance	Macrolide Resistance Methicillin Resistance	 Vancomycin Resistance
RESISTANCE GENE(S) TESTED - NO	T DETECTED:	
Broad Spectrum*	 Carbapenem Resistance 	 Quino/Fluoroquinolone
ORGANISM(S) TESTED - NOT DETE	CTED:	
BACTERIA:		
Acinetobacter baumannii	 Klebsiella pneumoniae 	 Providencia stuartii
Aerococcus urinae Citrobacter koseri Coagulase Negative Staph Group** Corynebacterium riegelii Klebsiella aerogenes	 Morganella morganii Mycobacterium tuberculosis Mycoplasma genitalium Mycoplasma hominis Pantoea apolomerans 	 Pseudomonas aeruginosa Serratia marcescens Staphylococcus aureus Streptococcus agalectiae Streptococcus ayoarnes
 Klebsiella oxytoca 	Proteus mirabilis	 Ureaplasma urealyticum
VIRUS:		
• Adenovirus • BK Virus • Human Herpes Virus-1,-2 (HSV -1,-2)	Human Herpes Virus-5 (CMV) Human Herpes Virus-6 (HHV-6)	 Human Herpes Virus-7 (HHV-7) JC Virus
STI:		
Chlamydia trachomatis	 Neisseria gonorrhoeae 	 Trichomonas vaginalis
YEAST:		
Candida albianas	- Condido alabasta	- Candida o acanellasia

References:

- Braid Spectrum Redittert genes identified confers redittence ecross multiple classes of entibiotics including periodilin derivatives, cephelosporins, monobactams and cerbapenens ** Congulare Negative Staphylococcus includes: Staphylococcus epidemidis, Staphylococcus heenolyticus, Staphylococcus iugdunenesis, and Staphylococcus
- **Asprophyticu**
- Eliti, Pacifire For estandad-spactrum bate-lactamese (Elitik), which are enzymen that confer realizance to most bate-lactam antibiotic, including periolities, establishing only, and the monobactam astronom. Infectione with Elitik, producing arguinant have been astroched with poor automes which are also plane planetare. Bioreconces an establishing cardination are been astroched with poor automes which are also plane planetare. Bioreconces and planetare. Bioreconces partners where which are also plane planetare.

Discluinver: This text was developed and ht performance characteristics determined by Pathnostics. It has not been cleared or approved by the UR Food and Drug Administration. The FDA has determined that such clearance or approvals is not necessary. This text is used for clinical purposes. It should not be regarded as investigational or for research. This laboratory is certified under the Chickel Laboratory Improvement Amendments of 1989 (CLIA-98) as qualified to perform high complexity clinical testing. Unite speciment resolved geneter than 5 days poor collection may give unveloble cells/nc, counts due to overgrowth of microorganism(c). Nativalways and Clockal Interference Microbie and reliance pares are detected through multiples PCR Pathogens are reported as: to pathogenic through exception of the pathogenic transmission of the pathogenic tr

Electronically Signed By: Bryan Chow, CLS, CGMBS, MB (ASCP) at 01/29/2020 09:31 am

CONFIDENTIAL HEALTHCARE INFORMATION Testing performed by Pathnostics | Hedical Director: Haher Badir, HD, FCAP CLIA# 05D2024468 | 17661 Coven, Irvine, CA 92614 | Phone: 714-866-821

First Last PAGE 2 OF 2

GUIDANCE UTI

Phonec MRN#

PATHNOSTICS GUIDANCE*COMPREHENSIVE: RECURRENT, PERSISTENT, OR COMPLICATED UTI REPORT + TB/STI Physician: Doctor

Facility: Facility

Phone:

Fax

Patient: First Last

DOB: 10-16-1959 Gender: F

Case#: PUXR19-010846-TS Collection Method: Voided Date Collected: 06-25-2019 Date Received: 06-26-2019 Date Reported: 01-29-2020

(?) Questions. Please call 714-966-1221 or visit pathnostics.com

RESULTS: PATHOGENIC DNA DETECTED

ORGANISM(S) TESTED - DETECTED: (See last page for Organism(x) Tested - Not Detected)

- Citrobacter freundli >100.000 cells/ml.
- Enterococcus faecalis >100,000 cells/ml.
- Viridans Group Strep >100,000 cells/mL Actinotignum schaalii 50,000-99,999 cells/mL
- Escherichia coli >100,000 cells/ml.

Alloscardovia omnicolens 50,000-99,999 cells/mL

LEGEND S = Pooled Phenotypic Serebikty Detected R = Pooled Phenotypic Resistance Detected R8D = Resistance Gene(x) Detected	Levelbrach	Tetracycline	Clareflexaeln	Nitrefurantoin	Piperacijin / Taxobactam	Forfomycin	Sulfamethozazolo / Trimethopolm	Gentamicin	Heropenem	Amstellts	Cetaclor	Ceffeitazone	Cettaddine	Nancomyclin	Amptellin / Sulbactam	Catapima	Catazolin	Cefestitin
Formulations	Payne		POW	P+	•	P0	Paylo	IM/TV	w	Poyte	PO	M/W	w	w	N	w	w	N
Pooled Phenotypic Sensitivity	\$	\$	\$	\$	5	R		R							n	R		
Resistance Gene(s) Detected					RGD					RGD	RGD	RGD	RGD	RGD	RGD	RGD	RGD	RGD
MIC Results (ug/mL)	•	2	•	32	16/4													

Citrobacter freundil	\checkmark	\checkmark	\checkmark	~	~		~	\checkmark	\checkmark			~	~			\checkmark		
Enterococcus faecalis	~	<	~	~	~	<			~	~				<	<			
Escherichia culi	~	~	~	~	~		~	\checkmark	~	~	~	~	~		\checkmark	\checkmark	\checkmark	~
Viridans Group Strep	~	~			~				~	~		~	~	~	~	~	~	
Actinotignum schaall	~	~	~	~				~		~		~		<				
Alfoscardovia omnicolens	~	~	~	\checkmark	~	~		\checkmark		\checkmark			\checkmark	~				

Check marks indicate situations for which antibiotic use is either FDA-approved or off label use for antibiotics is illustrated in peer review Iterature. References available upon request.

CONFIDENTIAL HEALTHCARE INFORMATION Testing performed by Pathnostics | Hedical Director: Haher Badir, HD, FCAP CLA# 05D3034468 | 17661 Cowen, Irvine, CA 93614 | Phone: 714-666-9321





PATHNOSTICS MOVING CARE FORWARD

RESULTS: PATHOGENIC DNA DETECTED

ORGANISM(S) TESTED - DETECTED: (See last page for Organism(s) Tested - Not Detected)

Citrobacter freundii >100,000 cells/mL
 Enterococcus faecalis >100,000 cells/mL

- Viridans Group Strep >100,000 cells/mL
- Actinotignum schaalii 50,000-99,999 cells/mL

Escherichia coli >100,000 celis/mL

Alloscardovia omnicolens 50,000-99,999 cells/mL

LEGEND							£											
S =Pooled Phenotypic Sensitivity Detected R =Pooled Phenotypic Resistance Detected RGD =Resistance Gene(s) Detected	Levofloxacin	Tetracycline	Ciproflexacin	Nitrofurantoin	Piperacillin / Tazobactam	Fastomycln	Sulfamethoxazole / Trimethopri	Gentamicin	Meropenem	Ampleillin	Cefacior	Ceftriaxone	Coftazidime	Vancomycin	Amplellin / Sulbactam	Cefepime	Cefazolin	Cefexitin
Formulations	PQ/1V	PO	PQ'IV	PO	IV	РО	PQ/IV	IM/IV	N	PQIV	PO	IM/IV	IV	IV	IV	N	N	IV
Pooled Phenotypic Sensitivity	s	s	s	s	s	R	R	R	R	R	R	R	R	R	R	R	R	R
Resistance Gene(s) Detected					RGD					RGD	RGD	RGD	RGD	RGD	RGD	RGD	RGD	RGD
MIC Results (ug/mL)	1	2	1	32	16/4													

Organism(s) Tested - Detected: 🗸 = Check marks are supportive data and are NOT patient specific.

Citrobacter freundii	~	\checkmark	~	~	~		\checkmark	\checkmark	\checkmark			\checkmark	\checkmark			\checkmark		
Enterococcus faecalis	\checkmark	\checkmark	<	\checkmark	<	~			\checkmark	~				~	<			
Escherichia coli	~	~	<	~	~		\checkmark	~	\checkmark	~	\checkmark	~	\checkmark		~	\checkmark	\checkmark	\checkmark
Viridans Group Strep ****	\checkmark	~			~				\checkmark	\checkmark		~	\checkmark	\checkmark	~	\checkmark	\checkmark	
Actinotignum schaalii	\checkmark	~	~	~				\checkmark		\checkmark		\checkmark		\checkmark				
Alloscardovia omnicolens	~	\checkmark	\checkmark	~	~	\checkmark		\checkmark		\checkmark			\checkmark	\checkmark				



= Check marks indicate situations for which antibiotic use is either FDA-approved or off label use for antibiotics is illustrated in peer review literature. References available upon request.



Clinical Trials Results





Current Clinical Studies

34

Studies	Internal Validation	Retrospective UTI Urology Study (2018-RUUS)	Prospective UTI Urology Study (2018-PUUS)	Retrospective VPA Home Care UTI Study (2019 R- VPA-HUS)
Purpose	Technical Validity	Clinical Validity	Clinical Utility	Clinical Utility
Cohort Size	NA	500	2,518	110,000
Population	NA	60 and older	Mean age of 73	Mean age 70
Sites	Internal Test Validation	Urology Clinic (1)	37 Urology Offices in 7 states with 75 physicians	Home Care Patients
Study Design	NA	Retrospective analysis comparing Guidance UTI to culture	Compare outcomes of patients treated based upon Culture or Guidance UTI results	Comparison of culture to Guidance UTI in diagnosing and managing patients with UTI
Principal Investigator	Michael Opel, PhD	Kirk Wojno, MD	Kirk Wojno, MD	Kirk Wojno, MD
Status	Completed	Completed	Paper Submitted	Paper Submitted

NOSTICS

Bacterial Detection

	Guidance Positive n (%)	Guidance Negative n (%)	Total n (%)	Agreement n (%)	p-value
Culture Positive	1,018 (40.5%)	80 (3.2%)	1,098 (43.7%)	92.7% (91.4,94.1%)	<0.001
Culture Negative	557 (22.2%)	856 <mark>(34.1%)</mark>	1,413 (56.3%)	60.6% (58.0,63.1%)	
Total	1,575 (62.7%)	936 <mark>(37.3%)</mark>	2,511 (100.0%)	74.6% (72.9,76.3%)	





Compare Ability to Detect Polymicrobial Infections

	Guidance Polymicrobial (<u>≥</u> 2 bacteria)	Guidance Positive (Monomicrobial)	Guidance Negative	Total	Agreement
Culture Polymicrobial (≥2 bacteria)	141 (5.6%)	25 (1.0%)	2 (0.1%)	168 (6.7%)	83.9%
Culture Positive (Monomicrobial)	425 (16.9%)	427 (17.0%)	78 <mark>(3.1%)</mark>	930 (37.0%)	45.9%
Culture Negative	268 (10.7%)	289 <mark>(11.5%)</mark>	856 (34.1%)	1,413 (56.3%)	60.6%
Total	834 (33.2%)	741 (29.5%)	936 (37.3%)	2,511 (100.0%)	56.7%





Bacterial Histogram-Prospective Study



Number of Prospective Patients Testing Positive

Retrospective Paper – In Press

Infectious Diseases

Multiplex PCR Based Urinary Tract Infection (UTI) Analysis Compared to Traditional Urine Culture in Identifying Significant Pathogens in Symptomatic Patients

Kirk J. Wojno, David Baunoch, Natalie Luke, Michael Opel, Howard Korman, Colieen Kelly, S. Mohammad A. Jafri, Patrick Keating, Dylan Haæiton, Stephany Hindu, Bridget Makhloouf, David Wenzler, Mansour Sabry, Frank Burks, Miguel Penaranda, David E. Smith, Andrew Korman, and Larry Sirls

OBJECTIVE To evaluate whether multiplex PCR-based molecular testing is noninferior to urine culture for detection of bacterial infections in symptomatic patients. METHO DS Retrospective record review of 582 consecutive elderly patients presenting with symptoms of lower urinary tract infection (UTI) was conducted. All patients had traditional urine cultures and PCR molecular testing run in parallel. RESULTS A total of 582 patients (mean age 77; range 60-95) with symptoms of lower UTI had both urine cultures and diagnostic PCR between March and July 2018. PCR detected unpathogens in 326 patients (56%, 326/582), while urine culture detected pathogens in 217 patients (37%, 217/582). PCR and culture agreed in 74% of cases (431/582): both were positive in 34% of cases (196/582) and both were negative in 40% of cases (235/582). However, PCR and culture disagreed in 26% of cases (151/582): PCR was positive while culture was negative in 22% of cases (130/582), and culture was positive while PCR was negative in 4% of cases (21/582). Polymicrobial infections were reported in 175 patients (30%, 175/582), with PCR reporting 166 and culture reporting 39. Further, polymicrobial infections were identified in 67 patients (12%, 67/582) in which culture results were negative. Agreement between PCR and urine culture for positive cultures was 90%, exceeding the noninferiority threshold of 85% (95% conflict of interest 85.7%-93.6%). CONCLUSION Multiplex PCR is noninferior to urine culture for detection and identification of bacteria. Further investigation may show that the accuracy and speed of PCR to diagnose UTI can significantly improve patient ourcomes. UROLOGY 136: 119-126, 2020 @ 2019 Elsevier Inc.

The additional write culture is commonly egolicit on the gold smalled for detection and identification of pathogens. However, evidence has here accumulating to appopriate of molecular methods such as RCR. With a natimicabilit esistence becoming both more comnon and complex, effective transmers of (univery marinfection) UTIs is even more dependent on the accurate identification of pathogens. Some organism can be

PCR scale cas be detained in a day or less, while calues can expart 2 or more days. Previous article have reported PCR to have both superior sentitivity and specificity, and have recommended PCR for pixel identification of pathogene in sepsi-¹³ and for diagnosis of general infections and accually manustited disease.¹⁴ meaning infections² tubercalosis, and generalization directions³. Few studies have compared multiples PCR with urine

factidious, and therefore difficult to grow in culture. Further,

Peaced a Diedenmen Dr Larry 54s received hourants for provention main a ALM-2019 and SLEU-2018. From the Competencies Undage A Distance of Michigan Heckbare Professionals.

From the Competitudes Orange-A Database of Michgan Healthoure Proportionals, Royal Oak, Mij the Pathnonics, Irvine, CA; and the Kelly Statistical Coundring, Calilevil CA.

Addres compandenze su Ma Wojne, M.D., Computenzier Undege - A Division of Mailigan Hashinaw Professionali, 31157 Woodward Ave, Royal Cali, M14007. David Bauerds, Ph.D., Patterative, 7064 Canes, in Neuro, CA, SSA-9, E solid Instrumentary and California and Systematic ress. Submitted: Amazer 32, 2019. Account for historication. Conder 10, 2019.

© 2019 Pathnostics. Published by Elsevier Inc.

This is an open access article under the CC BYNCND license. (http://creativecommons.org/licenses/byno.nd/4.0/)

against 20,12 and the fourth against 9 lacteria.¹¹ Polymicrobial infections may occur in as many of 39% of UTIs^{16,15} and can display enhanced virulence and https://doi.org/10.1016/j.urdog.2019.10.018 0030.4295

culture for diagnosis of UTIs and acute cystitis. Although

several studies have compared performance of PCR with

urine culture for detection of a single pathogen, only 4

have tested multiplex PCR: one against 15 bacteria,10 a

second against 14 bacteria together with 6 fungi,11 a third





Genotype v Phenotype





Resistance Gene and Susceptibility Agreement

ABR and Genotype Agree		ABR and Genotype Disagree	
Bacteria Sensitive & No R. Genes	Bacteria Resistant & R. Genes Present	Bacteria Resistant & R. Genes Not Present	Bacteria Sensitive & R. Genes Present
42%	16%	24%*	18%**
58%		42%	

*There are either more genes to detect or organism interactions impacting resistance patterns

**Genes may be present in low quantity, inactive, or nonfunctional





Bacteria Interactions











Figure 5: Correlations between organisms, excluding E. coli. The strength of the correlation is represented by the width of the edge connecting the genes. Only correlations greater than 0.1 shown.



Organism Interactions Impact Susceptibility Results



Pair increases resistance.

Would expect the antibiotic resistance levels to equal *E. coli. K. pneumonia* causes resistance to drop.

PATHNOSTICS

MOVING CARE FORWARD



Organism Interactions: Cephalosporin

Antibiotic	Bacteria that increase odds of	Bacteria that decrease odds	Interactions
	resistance (odds-ratio)	of resistance (odds-ratio)	
Cefaclor	E. faecalis ▲ (2.3)Pseudomonas aeruginosa (8.6)	Klebsiella pneumonia (0.69)	 E. coli & E. faecalis and Klebsiella & E. faecalis together <u>decrease</u> the odds of resistance E. coli & Klebsiella together increase the odds of resistance
Cefazolin	Actinobaculum schaalli (1.3) Pseudomonas aeruginosa (2.2)		 Pseudomonas & E. faecalis together <u>decrease</u> the odds of resistance E. coli & Klebsiella together <u>increase</u> the odds of resistance
Cefepime	E. faecalis ▲ (9.4)E. coli (1.7)		• E. coli & Klebsiella and E. coli & Proteus together <u>increase</u> the odds of resistance
Cefoxitin			 Pseudomonas & E. faecalis together <u>decrease</u> the odds of resistance E. coli & Klebsiella together <u>increase</u> the odds of resistance
Ceftazidime	E. faecalis \blacktriangle (7.3)	Klebsiella pneumonia (0.53)	 CNS & Klebsiella together <u>decrease</u> the odds of resistance E. coli & Klebsiella and E. coli & Pseudomonas together <u>increase</u> the odds of resistance
Ceftriaxone	E. faecalis ▲ (6.9)	Klebsiella pneumonia (0.65) Proteus mirabilis (0.41) VGS (0.73)	 Pseudomonas & CNS together<u>decrease</u> the odds of resistance E. coli & Klebsiella together <u>increase</u> the odds of resistance





Conclusions

- Multiplex PCR for UTI demonstrates
 - improved detection of organisms in urine
 - improved detection of polymicrobial urinary tract infections
- Pooled Antibiotic Sensitivity and identification of Antibiotic Resistance Genes may better clinical utility than culture
 - Reduces the need for empiric therapy with quicker turn around time
 - Identifies effective antibiotics for the patient's specific group of organisms causing the UTI which improves time to symptom resolution







Thermo Fisher Scientific and its affiliates are not endorsing, recommending, or promoting any use or application of Thermo Fisher Scientific products presented by third parties during this seminar. Information and materials presented or provided by third parties are provided as-is and without warranty of any kind, including regarding intellectual property rights and reported results. Parties presenting images, text and material represent they have the rights to do so.



