



Fisher Scientific



Protector® Laboratory Fume Hoods

HIGH PERFORMANCE, ENERGY SAVINGS



*Protecting your
laboratory environment*

LABCONCO®



Protector® XStream® Laboratory Hoods

The most energy efficient Labconco fume hood ever

The patented* Protector XStream Laboratory Hood was engineered to be the best containing fume hood possible. Testing shows the Protector XStream easily meets containment per SEFA-1† low velocity hood standards when subjected to the ASHRAE 110‡ test protocol with results of less than 0.05 ppm leak rate when tested at 4.0 lpm at OSHA-recognized 60 fpm face velocity.

During independent testing**, the Protector XStream Hood was challenged well beyond the SEFA-1† standards. With a face velocity of 40 fpm and sash fully open, the Protector XStream was subjected to 50 fpm cross drafts, NIH‡ protocol, and tracer gas measurements in the chest of the mannequin. In all scenarios, the Protector XStream allowed **0.00 ppm** average

level of tracer gas outside the fume hood. Although your safety officer or industrial hygienist will determine the actual face velocity setting for your laboratory, **the ability of the Protector Hood to contain under these adverse conditions sets a new standard of safety.**

Safety is foremost, but energy savings is equally impressive. Although face velocity is a factor, it's the **volumetric rate (CFM)** that determines the energy consumption of a fume hood. Operating a 6' Protector XStream Hood at 60 fpm face velocity, with the sash in its fully open position, requires only **690 CFM**. Regardless of your desired operating face velocity, the Protector XStream yields the **lowest required CFM**.

Energy savings translates to dollar savings. The Protector XStream Hood provides an excellent economic payback when compared to traditional by-pass hoods operated at 80 or 100 fpm. For example, a 6' Protector XStream Hood with sash fully open when operated at 60 fpm face velocity consumes a mere 690 CFM as previously mentioned. Compared to a traditional by-pass hood operated at 100 fpm (1250 CFM), the Protector XStream Hood offers significant savings, which adds up to annual dollar savings per year of \$3920.††


For even greater savings, the Protector XStream may be factory-prepared to accommodate a variable air volume system without the need for by-pass modifications.

Using the concepts of fluid dynamics, Labconco engineers designed the Protector XStream Laboratory Hood to produce **horizontal airflow, which reduces the tendencies for turbulence.** The innovative and aerodynamic designs of the sash

handle, air foil, upper dilution air supply and rear downflow baffle work in concert to produce horizontal airflow patterns that significantly reduce concentrations of chemical contaminants throughout the work area, particularly near the oper-


ator's breathing zone and at the work surface. Depending on sash position, tendencies for air turbulence, vortexing and **"the roll"** frequently observed during traditional fume hood smoke tests are virtually eliminated.

Traditional By-Pass Hood Design



Smoke tests on traditional hoods show the tendency for contaminants generated in the interior to roll forward producing high concentrations of contaminants behind the sash in close proximity to the user's breathing zone.

Protector XStream Hood Design



In contrast, smoke tests on Protector XStream Hoods show contaminants removed in a single pass and a remarkable lack of turbulence. Horizontal laminar air flowing toward the baffle forces contaminants to the rear interior, away from the user. The upper dilution air supply sweeps the upper interior to eliminate stagnant pockets of air and to prevent contaminants from concentrating behind the sash.

*U.S. Patent No. 6,461,233


**Independent testing by AccuTec Services, Inc., Lee's Summit, MO, National Environmental Balancing Bureau (NEBB)-Certified, Professional Engineer

†See back cover for a list of regulations, standards and registered trademarks.

††See page 5 for energy savings details.



Protector® XStream® Laboratory Hoods




Requires
Ductwork and
Blower

Upper Dilution Air Supply

The sash interior is constantly bathed with room air from the dilution supply above the work area to eliminate chemical fumes along the sash plane, near the critical breathing zone. A small percentage (5-10%) of the required air volume is introduced through the dilution air supply to ensure maximum containment.

No additional blowers are required.



Rear Downflow Dual Baffle System

The slots in the primary baffle direct inflow air in non-turbulent streams from the hood face into the baffle in a single pass. The secondary baffle, located between the primary baffle and the back wall, counteracts the upward air streams that create roll in traditional hoods by forcing the air movement downward before exhausting.


No moving components are used.

The best containing Labconco fume hood ever

Clean-Sweep™ Sash Handle

The sash handle includes Clean-Sweep slots to bleed air into the hood chamber and direct chemical fume concentrations away from the user's breathing zone.


The slim-line radiused sash handle sweeps airflow into the hood with minimal turbulence.



Opti-Zone™ Baffle


The Opti-Zone Baffle decreases the typical face velocity variations found with other baffles. The unique slot pattern and sizes increase velocities in the middle and at the work surface of the hood where it is needed while slowing velocities at the corners. This uniformity lowers the required average face velocity necessary for containment.

Tapered slots decrease resistance to air entering the baffle.



Eco-Foil™ Air Foil

The Eco-Foil reduces energy consumption by 7-10% compared to flat air foils while its aerodynamic curve allows air to sweep the work surface for maximum containment. Clean-Sweep openings pull inflow air from under the air foil forcing air into non-turbulent air streams. The curve is comfortable for arms resting on it while encouraging users to keep fume-generating items well within the hood's interior.





Protector® XStream® Laboratory Hoods



5' Protector XStream Laboratory Hood 10-369-105 is shown with SpillStopper Work Surface 10-369-182, Protector Standard Storage Cabinet 16-305-81 and Protector Acid Storage Cabinet 16-305-87.

All models feature:

- By-pass airflow design with variable air volume compatibility.
- Eco-Foil Air Foil with aerodynamic Clean-Sweep™ airflow openings.*
- Cord-Keeper™ Slots on left and right side of air foil.
- Upper Dilution Air Supply.*
- Glacier white powder-coated steel exterior.
- Rear Downflow Dual Baffle System.*
- Chemical-resistant, fiberglass-reinforced, composite panel liner and baffles with flame spread index less than 25 per ASTM E84**. Baffles are removable for cleaning.
- Opti-Zone™ Baffle with tapered slots.
- Tempered safety glass vertical-rising sash.
- Powder-coated sash handle with aerodynamic Clean-Sweep™ openings.*
- 37.5" (95.3 cm) high sightline from the work surface to the header panel.
- Removable front and side panels, and front and interior service access panels for access to plumbing and electrical wiring.

*U.S. Patent No. 6,461,233

**See back cover for a list of regulations, standards and registered trademarks.



Heights of switches, electrical receptacle and service fixtures meet requirements of Americans with Disabilities Act (ADA).



Exclusive Feature

- Pre-wired T8 fluorescent lighting with vapor-proof design and ADA-compliant light and blower switches.
- Sash stop located at 18" (45.7 cm) sash opening position.
- Powder-coated stainless steel, 12.81" (32.5 cm) ID exhaust connection(s).

All models conform to the following regulations and standards:**

- CFR 29, Part 1910
- SEFA 1-2010
- NFPA 45-2011
- ASTM E84-09C
- ASHRAE 110-1995
- ANSI Z9.5-2011
- UL 61010-1
- CAN/CSA C22.2 No. 61010-1
- UL 1805
- CE Conformity Marking (230 volt models)
- SEFA 8-2010, Cabinet Surface Finish Tests

Fixtured models feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing for gas and lower left side with copper tubing for cold water. Components for converting either or both fixtures to air and vacuum are provided. Inlet tubing is not provided.
- One pre-wired GFCI electrical duplex receptacle on lower right side and, on 8' models only, one additional pre-wired GFCI electrical duplex receptacle on the lower left side.

Required accessories not included:

- Remote Blower.
- Ductwork.
- Work Surface. See page 15.
- Base Cabinet or Stand.

Optional accessories for on-site installation include:

- Service Fixture Kits.
- Electrical Duplex Kits.
- Guardian Airflow Monitor Kits.
- Ceiling Enclosure and Rear Finish Panel Kits.
- Distillation Grid Kits.
- Sash Stop Kits.

Contact Fisher for ordering information on accessories.



Ordering Information, Airflow Data & Energy Savings

PROTECTOR® XSTREAM® LABORATORY HOODS

Nominal Width	Fisher Catalog Number	Electrical Requirements	Service Fixtures	Electrical Duplex	Shipping Weight lbs
4 Feet	10-369-100	100-115V, 50/60Hz, 20A	0	0	440
4 Feet	10-369-101	100-115V, 50/60Hz, 20A	2	1	440
4 Feet	10-369-102	208-230V, 50/60Hz, 5A	0	0	440
4 Feet	10-369-103	208-230V, 50/60Hz, 5A	2	0	440
5 Feet	10-369-104	100-115V, 50/60Hz, 20A	0	0	525
5 Feet	10-369-105	100-115V, 50/60Hz, 20A	2	1	525
5 Feet	10-369-106	208-230V, 50/60Hz, 5A	0	0	525
5 Feet	10-369-107	208-230V, 50/60Hz, 5A	2	0	525
6 Feet	10-369-108	100-115V, 50/60Hz, 20A	0	0	600
6 Feet	10-369-109	100-115V, 50/60Hz, 20A	2	1	600
6 Feet	10-369-110	208-230V, 50/60Hz, 5A	0	0	600
6 Feet	10-369-111	208-230V, 50/60Hz, 5A	2	0	600
8 Feet	10-369-112	100-115V, 50/60Hz, 20A	0	0	770
8 Feet	10-369-113	100-115V, 50/60Hz, 20A	2	1	770
8 Feet	10-369-114	208-230V, 50/60Hz, 5A	0	0	770
8 Feet	10-369-115	208-230V, 50/60Hz, 5A	2	0	770

Total Exhaust CFM and Static Pressure @ 28" Sash Opening (100% Open)

Face Velocity (fpm)			Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)					
Sash @ Full Open (28")	4' Hood		5' Hood		6' Hood		8' Hood	
	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.
	100	705 0.26	930 0.32	1150 0.41	1600 0.29			
	80	565 0.17	745 0.20	920 0.26	1280 0.19			
	60	425 0.09	560 0.12	690 0.15	960 0.10			

Total Exhaust CFM and Static Pressure @ 18" Sash Opening (62.5% Open)

Face Velocity (fpm)			Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)						
Sash @ 62.5% Open (18")	4' Hood		5' Hood		6' Hood		8' Hood		
	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.	
	100	440	0.10	580	0.12	720	0.16	1000	0.11
	80	350	0.06	465	0.08	575	0.10	800	0.07
	60	265	0.04	350	0.05	430	0.06	600	0.04



Energy Savings & Dimensional Data

PROTECTOR® XSTREAM® LABORATORY HOODS

The Protector XStream Laboratory Hood shows significant savings over its lifetime when compared to a typical fume hood. The CFM usage and related energy costs associated with exhausting tempered air from the laboratory to the outside are provided below. Maximum savings are achieved using a Protector XStream Laboratory Hood

operating at 60 fpm with a variable air volume system. Protector XStream Laboratory Hoods are compatible for use with variable air volume (VAV) systems. Please Contact Fisher for ordering information on factory preparing Protector XStream Laboratory Hoods to a specific VAV controller cutout.

Energy Savings Dollars Compared to a Typical Fume Hood

	CFM	Dollars/Year	Dollars/Lifetime	Lifetime Dollar Savings Compared to Typical Hood
6' Typical Hood @ 100 fpm, full open sash (28"), constant volume†	1250	\$8,750	\$131,250	0
6' XStream at 100 fpm, full open sash (28"), constant volume†	1150	\$8,050	\$120,750	\$10,500
6' XStream at 60 fpm, full open sash (28"), constant volume†	690	\$4,830	\$72,450	\$58,800
6' XStream at 60 fpm, 62.5% open sash (18"), constant volume†	430	\$3,010	\$45,150	\$86,100
6' XStream at 60 fpm, variable air volume††	N/A	\$1,883	\$28,245	\$103,005


†Based on average annual dollars per CFM of \$7.00, fume hood operating 24 hours a day and 5 days per week (6240 hours per year).

Average annual dollar per CFM range from \$5.00 to \$12.00 depending on geographic location. Lifetime calculations are based on 15 years.


††Based on 8 hours per day with 18" sash opening and 60 fpm face velocity, and remaining time with sash closed.

Closed sash air volume is based on ANSI Z9.5 minimum of 150 Air Changes per hour (ACH), and \$0.0000187/ft³ air.


TOP 8' Hoods




TOP 4', 5' & 6' Hoods



	A	B	C
4' Hood	48.0" (121.9 cm)	38.1" (96.8 cm)	24.0" (61.0 cm)
5' Hood	60.0" (152.4 cm)	50.1" (127.3 cm)	30.0" (76.2 cm)
6' Hood	72.0" (182.9 cm)	62.1" (157.7 cm)	36.0" (91.4 cm)
8' Hood	96.0" (243.8 cm)	86.1" (218.7 cm)	—



FRONT



SIDE

BIM Revit® and detailed AutoCAD® drawings available at www.labconco.com. See back cover for trademark information.



Protector® Premier® Laboratory Hoods

Protector® Premier® Laboratory Hoods incorporate a sleek interior with a molded one-piece fiberglass liner, the signature feature of Labconco's leading line of fume hoods since 1961. The one-piece liner of specially-formulated, fiberglass-reinforced polyester offers corrosion and fire resistance and easy clean up. Without seams, the interior has fewer points of deterioration for longer life.

Like the Protector XStream Hoods, Protector Premier Hoods incorporate many containment-enhancing features including

Clean-Sweep™ technology and the Eco-Foil™ air foil. Testing confirms the Protector Premier Hood meets the SEFA-1 * standard of a low velocity, high performance hood and may be operated as low as 60 fpm. These hoods are offered for use with a remotely-located blower or with a built-in blower — the only high performance hood with built-in blower available anywhere.

Features and benefits unique to Protector® Premier® Laboratory Hoods are described below. Additional features are detailed on page 8-12.



One-piece molded fiberglass liner offers superior corrosion and chemical resistance, durability and light reflectivity. Its seamless and smooth, radiused corners make cleaning easy and results in less deterioration for longer life. On models for use with remote blowers, as shown above, the exhaust connections are also seamless, molded fiberglass.



Models with built-in blower are available. Built-in blower is belt-driven with molded thermoplastic housing and coated aluminum impeller that is non-sparking and corrosion-resistant. The blower is available with standard or explosion-proof motor.

*See back cover for a list of regulations, standards and registered trademarks.



Protector® Premier® Laboratory Hoods

FOR USE WITH REMOTE BLOWER



4' Protector Premier Laboratory Hood 10369117 is shown with SpillStopper Work Surface 10369147 and Protector Standard Storage Cabinet 1630580.

Fixtured models feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing for gas and lower left side with copper tubing for cold water. Components for converting either or both fixtures to air and vacuum are provided. Inlet tubing is not provided.
- One pre-wired GFCI electrical duplex receptacle on lower right side and, on 8' and larger models only, one additional pre-wired GFCI electrical duplex receptacle on the lower left side.

Required accessories not included:

- Remote Blower.
- Ductwork.
- Work Surface. See page 15.
- Base Cabinet or Stand.

Optional accessories for on-site installation include:

- Service Fixture Kits.
- Electrical Duplex Kits.
- Guardian Airflow Monitor Kits.
- Distillation Grid Kits.
- Sash Stop Kits.
- Ceiling Enclosure and Rear Finish Panel Kits.

All models feature:

- By-pass airflow design.
- Eco-Foil™ Air Foil with aerodynamic Clean-Sweep™ airflow openings.*
- Cord-Keeper™ Slots on left and right side of air foil.
- Glacier white powder-coated steel exterior.
- One-piece molded fiberglass liner and pre-set baffle(s) with flame spread less than 25 per ASTM E-84.**
- Tempered safety glass vertical-rising sash with powder-coated sash handle.
- 37.5" (95.3 cm) high sightline from the work surface to header panel.
- Removable front and side panels, and front service access panels for access to plumbing and electrical wiring.
- Pre-wired T8 fluorescent lighting with vapor-proof design and ADA-compliant light and blower switches.
- Molded fiberglass 12.81" ID exhaust connection(s).

All models conform to the following regulations and standards**:

- CFR 29, Part 1910
- NFPA 45-2011
- ASHRAE 110-1995
- UL 61010-1
- UL 1805
- SEFA 1-2010
- ASTM E84-09C
- ANSI Z9.5-2011
- CAN/CSA C22.2 No. 61010-1
- CE Conformity Marking (230 volt models)
- SEFA 8-2010, Cabinet Surface Finish Tests



*U.S. Patent No. 6,461,233

**See back cover for a list of regulations, standards and registered trademarks. Heights of switches, electrical receptacle and service fixtures meet requirements of ADA.



Exclusive Feature



Dimensional Data & Airflow Data

PROTECTOR® PREMIER® LABORATORY HOODS FOR USE WITH REMOTE BLOWER


Total Exhaust CFM and Static Pressure @ 28" Sash Opening (100% Open)

Face Velocity (fpm)	Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)							
	4' Hood		5' Hood		6' Hood		8' Hood	
Sash @ Full Open (28")	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.
100	725	0.22	955	0.31	1180	0.41	1640	0.28
80	580	0.14	765	0.20	945	0.26	1310	0.18
60	435	0.08	575	0.11	710	0.15	985	0.10


Total Exhaust CFM and Static Pressure @ 18" Sash Opening (62.5% Open)

Face Velocity (fpm)	Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)							
	4' Hood		5' Hood		6' Hood		8' Hood	
Sash @ 62.5% Open (18")	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.
100	450	0.09	595	0.12	735	0.16	1025	0.11
80	365	0.06	480	0.08	590	0.10	820	0.07
60	270	0.03	360	0.04	440	0.06	615	0.04


TOP 8' Hoods




TOP 4', 5' & 6' Hoods



	A	B	C
4' Hood	48.0" (121.9 cm)	38.1" (96.8 cm)	24.0" (61.0 cm)
5' Hood	60.0" (152.4 cm)	50.1" (127.3 cm)	30.0" (76.2 cm)
6' Hood	72.0" (182.9 cm)	62.1" (157.7 cm)	36.0" (91.4 cm)
8' Hood	96.0" (243.8 cm)	86.1" (218.7 cm)	—



FRONT



SIDE



Protector® Premier® Laboratory Hoods

WITH BUILT-IN BLOWER



Requires
Ductwork

5' Protector Premier Laboratory Hood 10-369-130 is shown with SpillStopper Work Surface 10-369-154, Protector Standard Storage Cabinet 16-305-81 and Protector Acid Storage Cabinet 16-305-87.

All models feature:

- By-pass airflow design.
- Built-in belt-driven, corrosion-resistant exhaust blower with adjustable sheave, molded thermoplastic housing and non-sparking, coated aluminum impeller.
- Eco-Foil™ Air Foil with aerodynamic Clean-Sweep™ airflow openings.*
- Cord-Keeper™ Slots on left and right side of air foil.
 - Glacier white powder-coated steel exterior.
- One-piece molded fiberglass liner and pre-set baffle(s) with flame spread less than 25 per ASTM E-84.**
- Tempered safety glass vertical-rising sash with powder-coated aluminum sash handle.
- 37.5" (95.3 cm) high sightline from the work surface and header panel.
- Removable front and side panels, and front service access panels for access to plumbing and electrical wiring.
- 10.8" ID exhaust connection (4' models); 12.8" ID exhaust connection (5' and 6' models).

All models conform to the following regulations and standards**:

- CFR 29, Part 1910
- NFPA 45-2011
- ASHRAE 110-1995
- UL 61010-1
- UL 1805
- SEFA 1-2010
- ASTM E84-09C
- ANSI Z9.5-2011
- CAN/CSA C22.2 No. 61010.1
- CE Conformity Marking (230 volt models)
- SEFA 8-2010, Cabinet Surface Finish Tests

Standard models feature:

- Pre-wired T8 fluorescent lighting with vapor-proof design, and ADA-compliant light and blower switches.

Explosion-proof models feature:

- Explosion-proof blower and incandescent light fixture (bulb not included).
- Furnished without switches, electrical receptacles and wiring.

Fixtured models feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing for gas and lower left side with copper tubing for cold water. Components for converting either or both fixtures to air and vacuum are provided. Inlet tubing is not provided.
- One pre-wired GFCI electrical duplex receptacle on lower right side.

Required Accessories not included:

- Ductwork.
- Work Surface. See page 15.
- Base Cabinet or Stand.

Optional accessories for on-site installation include:

- Service Fixture Kits.
- Electrical Duplex Kits.
- Guardian Airflow Monitor Kits.
- Distillation Grid Kits.
- Sash Stop Kits.
- Ceiling Enclosure and Rear Finish Panel Kits.

Built-in Blower Maximum External Static Pressure @ 100 fpm and with Sash Full Open (28")

Hood Width	CFM	S.P.	Nominal Ductwork Diameter	Equivalent Resistance†
4 Feet	725	0.17"	10"	75
5 Feet	955	0.12"	12"	75
6 Feet	1180	0.17"	12"	70

*U.S. Patent No. 6,461,233

**See back cover for a list of regulations, standards and registered trademarks.

†Equivalent resistance in feet of straight duct.

Heights of switches, electrical receptacle and service fixtures meet requirements of ADA.




Exclusive Feature

BIM Revit® and detailed AutoCAD® drawings available at www.labconco.com. See back cover for trademark information.





Dimensional Data

PROTECTOR® PREMIER® LABORATORY HOODS WITH BUILT-IN BLOWER



	A	B	C
4' Hood	48.0" (121.9 cm)	38.1" (96.8 cm)	10.8" ID (27.5 cm)
5' Hood	60.0" (152.4 cm)	50.1" (127.3 cm)	12.8" ID (32.5 cm)
6' Hood	72.0" (182.9 cm)	62.1" (157.7 cm)	12.8" ID (32.5 cm)





Ordering Information

PROTECTOR® PREMIER® LABORATORY HOODS

Nominal Width	Fisher Catalog Number	Electrical Requirements	Service Fixtures	Electrical Duplex	Built-in Blower	Shipping Weight lbs
4 Feet	10-369-116	100-115V, 50/60Hz, 10A	0	0	No	325
4 Feet	10-369-117	100-115V, 50/60Hz, 20A	2	1	No	335
4 Feet	10-369-118	208-230V, 50/60Hz, 5A	0	0	No	325
4 Feet	10-369-119	208-230V, 50/60Hz, 20A	2	0	No	335
4 Feet	10-369-120	100-115V, 50/60Hz, 10A	0	0	Yes	365
4 Feet	10-369-121	100-115V, 50/60Hz, 20A	2	1	Yes	375
4 Feet	10-369-122	208-230V, 50/60Hz, 5A	0	0	Yes	370
4 Feet	10-369-123	208-230V, 50/60Hz, 20A	2	0	Yes	380
4 Feet	10-369-124	100-115V, 50/60Hz, 10A	0	0	Yes (EP)	375
5 Feet	10-369-125	100-115V, 50/60Hz, 10A	0	0	No	410
5 Feet	10-369-126	100-115V, 50/60Hz, 20A	2	1	No	420
5 Feet	10-369-127	208-230V, 50/60Hz, 5A	0	0	No	410
5 Feet	10-369-128	208-230V, 50/60Hz, 20A	2	0	No	420
5 Feet	10-369-129	100-115V, 50/60Hz, 10A	0	0	Yes	450
5 Feet	10-369-130	100-115V, 50/60Hz, 20A	2	1	Yes	460
5 Feet	10-369-131	208-230V, 50/60Hz, 5A	0	0	Yes	455
5 Feet	10-369-132	208-230V, 50/60Hz, 20A	2	0	Yes	465
5 Feet	10-369-133	100-115V, 50/60Hz, 10A	0	0	Yes (EP)	460
6 Feet	10-369-134	100-115V, 50/60Hz, 10A	0	0	No	485
6 Feet	10-369-135	100-115V, 50/60Hz, 20A	2	1	No	495
6 Feet	10-369-136	208-230V, 50/60Hz, 5A	0	0	No	485
6 Feet	10-369-137	208-230V, 50/60Hz, 20A	2	0	No	495
6 Feet	10-369-138	100-115V, 50/60Hz, 10A	0	0	Yes	525
6 Feet	10-369-139	100-115V, 50/60Hz, 20A	2	1	Yes	535
6 Feet	10-369-140	208-230V, 50/60Hz, 5A	0	0	Yes	530
6 Feet	10-369-141	208-230V, 50/60Hz, 20A	2	0	Yes	540
6 Feet	10-369-142	100-115V, 50/60Hz, 10A	0	0	Yes (EP)	535
8 Feet	10-369-143	100-115V, 50/60Hz, 10A	0	0	No	650
8 Feet	10-369-144	100-115V, 50/60Hz, 20A	2	1	No	660
8 Feet	10-369-145	208-230V, 50/60Hz, 5A	0	0	No	650
8 Feet	10-369-146	208-230V, 50/60Hz, 20A	2	0	No	660






Dimensional Data

SPILLSTOPPER™ WORK SURFACES

SpillStopper Work Surfaces are molded from a special formulation of corrosion-resistant epoxy resins. Dished and contoured to conform to the interior liner of Protector Laboratory Hoods. The front edge has a large radius to aerodynamically direct airflow into the hood. SpillStopper Work

Surfaces may be ordered with a pre-cut 6" x 3" oval cupsink cutout. Cupsink **(16-306-15)** and Cupsink Cover are sold separately unless otherwise noted.



16-306-15 Oval Polypropylene Cupsink
Mounts in work surface with cupsink cutout, 3.0" x 6.0" (7.6 x 15.2 cm). 1.5" National Pipe Straight Mechanical (NPSM) thread. Shipping weight 4 lbs. (2 kg)



Ordering Information

SPILLSTOPPER™ WORK SURFACES FOR **PROTECTOR XSTREAM** LABORATORY HOODS

Hood Width	Fisher Catalog Number	Left Cupsink Cutout	Right Cupsink Cutout	Shipping Weight
4 Feet	10-369-175	None	None	120
4 Feet	10-369-177	None	Rear	120
4 Feet	10-369-178*	Side	None	120
4 Feet	10-369-179*	None	Side	120
4 Feet	10-369-180	Rear	Rear	120
4 Feet	10-369-181*	Side	Side	120
4 Feet	10-369-176*	Rear	None	120
5 Feet	10-369-184*	None	None	160
5 Feet	10-369-184	None	Rear	160
5 Feet	10-369-185*	Side	None	160
5 Feet	10-369-186*	None	Side	160
5 Feet	10-369-187	Rear	Rear	160
5 Feet	10-369-188*	Side	Side	160
5 Feet	10-369-183*	Rear	None	160
6 Feet	10-369-189	None	None	220
6 Feet	10-369-191	None	Rear	220
6 Feet	10-369-192*	Side	None	220
6 Feet	10-369-193*	None	Side	220
6 Feet	10-369-194	Rear	Rear	220
6 Feet	10-369-195*	Side	Side	220
6 Feet	10-369-190*	Rear	None	220
8 Feet	10-369-196	None	None	250
8 Feet	10-369-198	None	Rear	250
8 Feet	10-369-199*	Side	None	250
8 Feet	10-369-200*	None	Side	250
8 Feet	10-369-201	Rear	Rear	250
8 Feet	10-369-202*	Side	Side	250
8 Feet	10-369-197*	Rear	None	250

*Not compatible with Protector Solvent Storage Cabinets

BIM Revit® and detailed AutoCAD® drawings available at www.labconco.com. See back cover for trademark information.



Ordering Information

SPILLSTOPPER™ WORK SURFACES FOR PROTECTOR PREMIER LABORATORY HOODS

Hood Width	Fisher Catalog Number	Left Cupsink Cutout	Right Cupsink Cutout	Shipping Weight
4 Feet	10-369-147	None	None	110
4 Feet	10-369-149	None	Rear	110
4 Feet	10-369-150*	Side	None	110
4 Feet	10-369-151*	None	Side	110
4 Feet	10-369-152	Rear	Rear	110
4 Feet	10-369-153*	Side	Side	110
4 Feet	10-369-148*	Rear	None	110
5 Feet	10-369-154	None	None	150
5 Feet	10-369-156	None	Rear	150
5 Feet	10-369-157*	Side	None	150
5 Feet	10-369-158*	None	Side	150
5 Feet	10-369-159	Rear	Rear	150
5 Feet	10-369-160*	Side	Side	150
5 Feet	10-369-155*	Rear	None	150
6 Feet	10-369-161	None	None	205
6 Feet	10-369-163	None	Rear	205
6 Feet	10-369-164*	Side	None	205
6 Feet	10-369-165*	None	Side	205
6 Feet	10-369-166	Rear	Rear	205
6 Feet	10-369-167*	Side	Side	205
6 Feet	10-369-162*	Rear	None	205
8 Feet	10-369-168	None	None	240
8 Feet	10-369-170	None	Rear	240
8 Feet	10-369-171*	Side	None	240
8 Feet	10-369-172*	None	Side	240
8 Feet	10-369-173	Rear	Rear	240
8 Feet	10-369-174*	Side	Side	240
8 Feet	10-369-169*	Rear	None	240

*Not compatible with Protector Solvent Storage Cabinets



Standards & Registered Trademarks

Standards

Key aspects of standards and codes as they relate to laboratory ventilation are summarized below.

ASHRAE 110-1995 Method of Testing Performance of Laboratory Fume Hoods (ANSI Approved)

Evaluates fume hood's containment characteristics.

- Three part test: Smoke generation, face velocity profile, tracer gas release @ 4 liters per minute.
- Rated As Manufactured (AM), As Installed (AI) and As Used (AU).

American Society of Heating, Refrigerating and Air-Conditioning Engineers

1791 Tullie Circle NE
Atlanta, GA 30329
(404) 636-8400
www.ashrae.org

ANSI Z9.5-2011 Standard—Laboratory Ventilation

Covers entire laboratory ventilation system.

- Vertical stack discharge @ 2000-3000 fpm.
- New and remodeled hoods shall have a monitoring device.
- Ductless hoods should only be used with non-hazardous materials.

American Industrial Hygiene Association

2700 Prosperity Avenue, Suite 250
Fairfax, VA 22031
(703) 849-8888
www.aiha.org

Federal Register 29 CFR Part 1910

Occupational exposure to hazardous chemicals in laboratories

National Research Council Recommendations Concerning Chemical Hygiene in Laboratories (Non-mandatory) from "Prudent Practices."

- Fume hoods should have a continuous monitoring device.
- Face velocities should be between 60-100 linear feet per minute (lfpm).
- Average 2.5 linear feet of hood space per person.

Occupational Safety & Health Administration U.S. Department of Labor

200 Constitution Avenue, NW
Washington, DC 20210
(800) 321-6742
www.osha.gov

ASTM E84-09C Standard Test Method for Surface Burning Characteristics of Building Materials

Determines the relative burning behavior of the material by observing the flame spread along the specimen.

- Measures the flame spread and smoke development.
- Material is exposed to flaming fire for 10 minutes and the results measured and recorded.
- Results are compared to the indexes of mineral fiber cement board (flame spread and smoke development of zero) and red oak flooring (smoke development of 100).

ASTM International

100 Barr Harbor Drive
P.O. Box C700
West Conshohocken, PA 19428-2959
(610) 832-9585
www.astm.org

NFPA 45: Standard on Fire Protection for Laboratories Using Chemicals, 2011 edition

- Laboratory hoods should not be relied on for explosion protection.
- Fume hood exhaust air should not be recirculated.
- Services should be external to the hood.
- Materials of construction should have flame spread of 25 or less.

National Fire Protection Association

1 Batterymarch Park
Quincy, MA 02169-7471
(800) 344-3555 or (617) 770-3000
www.nfpa.org

NIH - Section 15991 Onsite Testing for Constant Volume Hoods - June 2006

- Follows ASHRAE test methods except for the following:
 1. 6 L tracer gas release rate instead of 4 L.
 2. Hood is loaded with boxes and cans.
 3. Rapid walk-by test.

National Institutes of Health

9000 Rockville Pike
Bethesda, MD 20892
(301) 496-4000
www.nih.gov

SEFA 1-2010 Laboratory Fume Hoods Recommended Practices

- High performance fume hood definition: hood with sash fully open and operating at 60 fpm contains at 4.0 AM 0.05
- Covers design, installation, testing, maintenance and safe use of laboratory fume hoods

SEFA 8-2010 Recommended Practices For Metal Laboratory Grade Furniture, Casework, Shelving and Tables, 8.0 Cabinet Surface Finish Tests

Defines test methods for evaluating the finish of laboratory furniture.

- Laboratory grade paint finishes shall withstand chemical exposure, hot water, and impact from a one-pound ball dropped from 12".
- Paint coating shall sufficiently adhere to the substrate.
- Paint shall be resistant to scratches.

Scientific Equipment & Furniture Association

1205 Franklin Avenue, Suite 320
Garden City, NY 11530
(516) 294-5424
www.sefalabs.com

UL 61010-1 Electrical Equipment for Laboratory Use

Specifies the general safety requirements for electrical equipment.

- Based on International Electrotechnical Commission (IEC) Publication 61010-1 with differences noted for U.S. use.
- Tests for protection against electrical shock, mechanical hazards, spread of fire, radiation, liberated gases, explosion and implosion.
- Tests for resistance to shock, vibration, impact, heat, moisture and liquids.

Underwriters Laboratories Inc.

333 Pfingsten Road
Northbrook, IL 60062-2096
(847) 272-8800
www.ul.com

CAN/CSA Standard C22.2 No. 1010.1 Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use

Specifies general safety requirements for electrical equipment.

- Design and methods of construction should provide adequate protection to the operator and the surrounding area against shock or burn, mechanical hazards, excessive temperature, spread of fire from the equipment, gas liberation, explosion or implosion.

Canadian Standards Association

5060 Spectrum Way, Suite 100
Mississauga, Ontario
L4W 5N6, CANADA
(800) 463-6727 or (416) 747-4044
www.csa.ca

ETL listing

ETL Testing Laboratories is a Nationally Recognized Testing Laboratory (NRTL). The ETL mark signifies that a product conforms to the following:

- UL Standard 61010-1 in the U.S.
- CAN/CSA Standard C22.2 No. 61010.1 in Canada.
- Products that bear the ETL mark are subjected to a comprehensive safety program that includes testing, listing, labeling and quarterly follow-up inspections.

Intertek Group

www.intertek.com

CE Marking

Indicates an electrical apparatus conformity to all safety and other directives/specifications presently required by the Council of European Communities.

- Electrical safety.
- Electromagnetic emissions testing — interference signals being output by the product.
- Electromagnetic immunity testing — the product does not respond to outside electromagnetic interference signals.

European Union

www.europa.eu

Registered Trademarks

AutoCAD® is a registered trademark of Autodesk.

ANSI® is a registered trademark of American National Standards Institute.

Revit® is a registered trademark of Autodesk.

SEFA® is a registered trademark of Scientific Equipment and Furniture Association.

UL® is a registered trademark of UL, LLC.



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com