



# Solana<sup>®</sup> Trichomonas ASSAY



## Stop compromising. Start treating.

Traditional methods for diagnosing Trichomonas require tradeoffs between cost and accuracy. Wet mount may be rapid and inexpensive, but it lacks in sensitivity (51% to 65%).<sup>1</sup> Similarly, rapid antigen tests are fast and easy, but are not as sensitive as their molecular counterparts.<sup>2</sup> The CDC recommends testing for Trichomonas using highly sensitive Nucleic Acid Amplification Tests, like Solana, to avoid underdiagnosing and undertreating the disease.

### Prompt initiation of patient management and antimicrobial therapy

- Results in just 30 minutes means you can help your patients get accurate, same-day treatment
- High volume throughput keeps pace with your busy laboratory

### True, Molecular Accuracy

- Confidence in your diagnosis
- More sensitive than wet mount, rapids, and culture

### Ease and Affordability

- Clinician collected vaginal swab and female urine samples
- Objective results
- Simple procedure
- Competitive pricing and flexible instrument placement options


The Accurate. Sustainable.  
Molecular **Trichomonas** Solution.

1. <http://www.cdc.gov/std/tg2015/trichomoniasis.htm>  
2. Garber, GE. The laboratory diagnosis of trichomonas vaginalis. *The Canadian Journal of Infectious Diseases and Medical Microbiology*. 2015;16(1):35-38

Solana Trichomonas Assay run with the Solana bench top instrument a mere 9.4" x 9.4" x 5.9"



# Procedure



# Clinical Performance

## Performance Characteristics of the Solana Trichomonas Assay with Clinician Collected Vaginal Swabs by Symptom Status compared to the Composite Reference Method Wet Mount and Culture

Site Number	Symptom Status	N	TP	FP	TN	FN	Prev %	Sensitivity % (95% CI)	Specificity % (95% CI)
Combined	Asymptomatic	501	50	5	446	0	10.0	100 (92.9 to 100)	98.9 (97.4 to 99.5)
	Symptomatic	542	69	7	465	1	12.9	98.6 (92.3 to 99.7)	98.5 (97.0 to 99.3)
	All	1043	119	12*	911	1*	11.5	99.2 (95.4 to 99.9)	98.7 (97.7 to 99.3)

\*Of the one thousand forty-three (1043) specimens evaluated a total of thirteen (13) specimens were discordant. Of the twelve (12) discordant (Solana Positive/Composite Reference Method Negative) specimens, four (4) were positive by an FDA-cleared *Trichomonas vaginalis* molecular assay. The one (1) discordant (Negative Solana/Composite Reference Method Positive) specimen, it was negative by an FDA-cleared *Trichomonas vaginalis* molecular assay.

## Performance Characteristics of the Solana Trichomonas Assay with Female Urine Specimens by Symptom Status compared to the Composite Reference Method

Site Number	Symptom Status	N	TP	FP	TN	FN	Prev %	Sensitivity % (95% CI)	Specificity % (95% CI)
Combined	Asymptomatic	501	49	7	444	1	10.0	98.0 (89.5 to 99.6)	98.4 (96.8 to 99.2)
	Symptomatic	543	65	10	463	5	12.9	92.9 (84.3 to 96.9)	97.9 (96.2 to 98.8)
	All	1044	114	17	907	6	11.5	95.0 (89.5 to 97.7)	98.2 (97.1 to 98.8)



Solana Trichomonas Assay – Swab – 48 Test Kit: Cat. #23043102  
 Solana Trichomonas Assay – Urine – 48 Test Kit: Cat. #23043103  
 Quidel Molecular Trichomonas Control Set – 80-Reaction Kit: Cat. #23043135  
 Solana Trichomonas Startup Kit – Swab: Cat. #23043105  
 Solana Trichomonas Startup Kit – Urine: Cat. #23043104

SSM304203EN00 (03/17)

© 2015 Thermo Fisher Scientific Inc. All rights reserved.  
 Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

### In the United States:

For customer service, call 1-800-640-0640  
 To fax an order, use 1-800-290-0290  
 To order online: [www.fisherhealthcare.com](http://www.fisherhealthcare.com)

