

LAB

R E P O R T E R

Science Innovations and Discoveries

NO. 2, 2021

Reducing the Environmental Impact of Disposable Face Masks

Also Featured:

Ability to Adapt, Prioritize Are Keys to **Boosting Productivity**

Addressing **Classroom Safety** Concerns During the Pandemic

The Basics of **Cleanroom Upkeep**

Pharmaceutical Applications of **Artificial Intelligence**

This issue of
Lab Reporter
is sponsored by 

 **fisher scientific**
part of Thermo Fisher Scientific

CONTENTS



20

20
COVER STORY
Reducing the Environmental Impact of Disposable Face Masks



32



4

4
Ability to Adapt, Prioritize Are Keys to Boosting Productivity

12
Addressing Classroom Safety Concerns During the Pandemic



12

32
The Basics of Cleanroom Upkeep

36
Pharmaceutical Applications of Artificial Intelligence



36

16
Check out more great content from the sponsor of this issue of *Lab Reporter*:
DuPont Protection Technologies

Note: To check pricing, please sign in to your account.

Supplier and Other Trademarks: Air-Tite, Ansell, AnsellGUARDIAN, Avantor, BDD, BioTek, Boekel, BRAND, BrandTech, Cidecon, Comfort-Air, Conflitk, Cytiva, Decon, Dentec, DuPont, DWK, FlaskScrubber, Greiner Bio-One, Guardian, Heidolph, Hei-VAP, IRADECON, IsoClean, J.T.Baker, Kinematica, KNF Neuberger, Labconco, LABOPORT, Laxco, Macron, METTLER TOLEDO, MICROLON, Milli-Q, MilliporeSigma, Nasco, OHAUS, ProShield, SafetyHeat, Sera-Xtracta, SmartHeat, SmartHousing, SmartLink, SmartPresence, SmartRate, Sonoco, SteamScrubber, Superdex, Synergy, TechNIGlove, ThermoSafe, Transferpette, Tyvek, Wells Lamont, Whirl-Pak, Youngstown Barrel and Drum

Trademarks of Thermo Fisher Scientific: Acros, AcroSeal, Alfa Aesar, Applied Biosystems, CryoExtra, Cryomed, CryoPlus, Fisher Scientific, Fisherbrand, Invitrogen, *Lab Reporter*, SuperScript, Taqman, Thermo Scientific, VILO

For trademark ownership, see fishersci.com/trademarks.

CONTACT US: United States: 1-800-766-7000 fishersci.com • Canada: 1-800-234-7437 fishersci.ca

LAB REPORTER

Science Innovations and Discoveries

Lab Reporter provides quick and easy access to today's cutting-edge products and trusted solutions for all of your scientific research and applications.

SUPPLIER ARTICLES



6
Ergonomic Glove Design
Ansell



8
Air-Sensitive Chemicals
Acros Organics/Alfa Aesar



15
Glassware Safety
DWK Life Sciences



16
Tyvek Protective Apparel
DuPont Protection Technologies



24
Cryogenic Storage
Thermo Scientific



28
Nitrile Allergies
TechNiGlove



31
Bead Mill Basics
Fisherbrand



38
Balance Static
METTLER TOLEDO



42
Pipetting Accuracy
BrandTech

SUPPLIER PRODUCT GUIDE

Air-Tite Nonsterile Syringes 26	Dentec Comfort-Air NxMD Respirators..... 34	Laxco SLI-6000 Fluorescence Microscopes 41
Ansell AnsellGUARDIAN Webinar..... 10	DuPont Protection Technologies Tyvek Apparel . 30	MilliporeSigma Milli-Q IX Pure Water Systems... 10
Applied Biosystems PCR Products..... 19	Fisherbrand Guarantee 47	Nasco Whirl-Pak Sampling Bags..... 44
Avantor J.T.Baker and Macron Fine Chemicals 35	Fisherbrand Sonicators..... 18	OHAUS Guardian Hotplate Stirrers..... 18
BD Biosciences Reagents 48	Greiner Bio-One ELISA Plates..... 26	SONOCO ThermoSafe EOS Shippers..... 34
BioTek Microplate Washers and Readers..... 45	Heidolph Hei-VAP Rotary Evaporators..... 47	Wells Lamont Glove Liners 40
Boekel Block Heaters..... 44	Kinematica Homogenizers 40	Youngstown Barrel and Drum
Cytiva Pathogen Kits and Columns..... 27	KNF LABOPORT Vacuum Pumps 46	Spill Containment Products..... 46
Decon Labs Disinfectants..... 11	Labconco Glassware Washers 14	

This publication was printed on recycled paper.



Visit fishersci.com/labreporter
or fishersci.ca/labreporter to
subscribe for free.

fisher scientific
part of Thermo Fisher Scientific



Ability to Adapt, Prioritize Are Keys to Boosting Productivity

By Kevin Ritchart

With limited staffing, a finite number of working hours, and a wide array of tasks to complete, lab managers are always on the lookout for methods that will make their facility more efficient. The challenge is finding ways to ramp up productivity without sacrificing the health and safety of employees or the accuracy of results.

From helping employees learn to manage their time more effectively to updating lab layouts to create a more functional working space, there are a variety of things lab managers can do to have a positive effect on productivity.

Adapt and Overcome

At its core, time management is the decision-making process that structures, protects, and adjusts a person's time as it relates to changing environmental conditions. Years of research have shown that the three key skills separating success from failure in time management are awareness, arrangement, and adaptation.

Awareness is maintaining a realistic view of time by recognizing it as a limited resource. Arrangement is organizing your tasks and goals to use your time effectively. Adaptation is monitoring your use of time while performing tasks, and adjusting to interruptions and changing priorities.

Prioritization and Productivity

When it comes to optimizing productivity, it's vital to focus on tasks that are aligned with your overall objectives and leave the less pressing matters for later. One method that can help in this pursuit is to employ the Eisenhower Matrix of Prioritization.

This approach helps you to prioritize your tasks by classifying them into four categories: urgent and important, important but not urgent, urgent but not important, and neither urgent nor important.

Using this approach can help reduce the number of unimportant tasks that are completed as a means of appearing busy, when in fact those tasks are having the opposite effect on your overall productivity.

Staffing and Sterilization

Like any area where people congregate and interact, labs had to adapt their operations and designs because of the COVID-19 pandemic.

Some of the more successful strategies for encouraging social distancing and enforcing new health and safety guidelines include staggering shifts to have fewer employees in the lab at one time, implementing new protocols for traffic flow, and enforcing the use of additional personal protective equipment.

But despite these additional steps, sanitizing remains the most effective means of combating infection, both inside and outside the lab. Life sciences labs, many of which already employ strict sterilization protocols, have adopted enhanced cleaning measures used by the hospitality and medical industries.

The equipment that's used most often is identified by the lab manager and cleaned on a regular basis to help reduce the risk of contamination. By employing this approach, the need to shut down equipment for longer periods of time to perform deep cleaning is minimized, which can help increase productivity.

Room to Move

While they've been given greater focus during the COVID-19 pandemic, expansion of lab space and enhanced sanitization are here to stay when it comes to optimizing safety in the lab.

Some companies may seek out new, larger facilities or reallocate existing space to accommodate HVAC enhancements, store additional cleaning supplies and personal protective equipment, and add bench space to give employees ample room to perform their tasks.

While these strategies for improving productivity and adapting lab layouts can be useful, it's important to note that there's no one-size-fits-all approach to enhancing productivity. It may take some time to develop an approach that fits your unique needs.

This content was inspired, in part, by "Designing Labs for Productivity," Lab Manager, February 28, 2020; and "Time Management Is About More Than Life Hacks," Harvard Business Review, January 29, 2020.

Kevin Ritchart is a Thermo Fisher Scientific staff writer.

Better by Design

How Ergonomically Designed Hand Protection Delivers a Safer, More Productive Workforce

The true economic cost of workplace injury is more than lost wages or time. Ongoing medical expenses, lowered production, and increased worker compensation premiums are the more obvious and measurable aspects, but there are hidden costs as well.

Although harder to assess, increased workplace injury rates contribute to heightened stress or anxiety among workers and can potentially lower both morale and productivity. Identifying and managing hazard risks before any issues arise requires the recognition of less obvious relationships and choosing personal protective equipment (PPE) that is ergonomically designed for specific conditions and tasks.



Common Working Conditions

In many industrial environments, repetitive manual tasks are unavoidable. Workers perform labor-intensive activities like lifting or lowering, pushing or pulling, and holding or restraining tools and other items. Such tasks require repetitive hand and arm movement and can lead to hand fatigue that puts more stress on other parts of the body, contributing to injuries.

These actions can become hazardous through repetition, the use of sustained pressure or force, maintaining prolonged or awkward postures, and ongoing exposure to vibration. These stress the

body and potentially lead to a variety of musculoskeletal disorders (MSDs). MSDs typically result from gradual wear and tear (triggered by repetition) or are brought on by sudden strenuous activity or unexpected movement.

The hands and arms are particularly susceptible to a range of conditions brought on by performing manual tasks, including:

- Muscle, ligament, or tendon sprains and strains
- Joint and bone injury or degeneration in the shoulder, elbow, or wrist
- Nerve injuries and compression
- Muscular or vascular disorders

Most of these can be acute or become enduring chronic conditions that threaten long-term health and productivity.

Prevalent but Preventable

MSDs are categorized as “body stressing” injuries and diseases. Workplace-related body stressing injury and disease is expensive and widespread. In the European Union, musculoskeletal disorders are the most frequently reported work-related health problem.¹ In Australia, over one third of the total number of worker injury cases and total economic costs are associated with body stressing or manual handling injury.²

In many cases, operations and safety managers include gloves in programs to mitigate the risk of hand or arm injury. However, issues arise when the chosen “solution” fails to address the immediate requirements. And a poorly selected glove can introduce its own risks.

For example, gloves that restrict movement of the hand or fingers require more muscle effort to perform the required tasks. This increases the risk of strain, which can lead to hand fatigue or persistent and painful conditions such as carpal tunnel syndrome.

Establishing Ergonomics

While safety personnel may be familiar with the concept of ergonomics, they may fail to consider the potential impact of PPE choices or incorporate ergonomic best practices and design principles. Primarily referencing the interaction between the worker’s musculoskeletal system and workspace, ergonomic design aims to minimize exposure to risk factors for MSD while increasing efficiency and comfort.

For gloves, this approach negates the oft-cited “comfort versus protection” argument — anecdotal evidence suggests that workers will simply ignore hand protection solutions that impede function or are otherwise uncomfortable.



If recurring or repetitive tasks are known to place strain on the muscles, nerves, and tendons in one’s hands, performing those activities while wearing thick, rigid, ill-fitting, slippery, or otherwise uncomfortable gloves can exacerbate the problem. Safety and operations managers should select a protective glove designed for the specific hazard types and the functions being performed, considering:

- **Fit:** gloves that are too small compromise movement; those that are too large restrict dexterity
- **Grip:** the amount of grip significantly affects the amount of muscle effort required to securely handle, hold, or manipulate objects

Content provided by:



- **Construction:** materials and comfort dictate overall wearability
- **Application-specific issues:** wet conditions, contact with abrasive materials, or the use of vibrating machinery present differing needs; assess each application to determine the optimal glove choice

Selection will vary based on these and other factors, but the process should include an assessment of the effect on manual performance to help workers operate both safely and efficiently.

Advances in glove technology deliver superior hand protection while providing much-needed support for musculoskeletal health, so workers no longer have to choose between comfort and protection. By supplying ergonomically designed protective gloves, employers, managers, and operators can help minimize the risk of MSD injury, support safety compliance, and maintain productivity levels.

Ansell ERGOFORM is a new technology used to design hand protection that supports musculoskeletal health during repetitive tasks. We measure the toll of occupational activities and apply cutting-edge technologies to produce our gloves. ERGOFORM certified products have been scientifically proven to deliver measurable improvements in worker comfort, fit, and productivity while reducing risk factors for ergonomic injury.



1. European Agency for Safety and Health at Work, Estimating the cost of work-related accidents and ill-health: An analysis of European data sources, p12

2. Safe Work Australia, The Cost of Work-related Injury and Illness for Australian Employers, Workers and the Community: 2012–13, p31



Air-Sensitive Chemistry: Practical and Safety Considerations

Air-sensitive chemicals are commonly used in academic and industrial chemistry laboratories. Examples of air-sensitive reagents include organometallic compounds like organo-lithium, -magnesium, -zinc, and -aluminium; hydrides; borane complexes; and alkaline metals and several transition metals in their (0) oxidation state. Some of these chemicals are also pyrophoric and spontaneously ignite under standard conditions when they come into contact with air.

Many incredibly useful synthetic reactions are air sensitive, including Grignard reactions, reductions using hydrides, metalations and transmetalations, carbolithiation, and many more. These reactions are widely used to synthesize fine chemicals, drugs, polymers, and many other products.


The exposure of these reactions to air, oxygen, or moisture can:

- Favor side-reactions that produce undesired products
- Decompose reagents so that no reaction occurs
- Cause fires, explosions, or other hazardous conditions

Handling Air-Sensitive Reactions

Air sensitivity can be separated into two categories: catalytically and stoichiometrically driven.

Figure 1



A typical example of catalytic air sensitivity, reported in Figure 1a, is the

ignition of residual organic solvents or H_2 gas adsorbed on the surface of a reduced Pd (0) catalyst, because of its exothermic oxidation in the presence of atmospheric oxygen. Stoichiometric air sensitivity is very common for many organometallic species, such as the formation of butane gas by reaction of butyl lithium with water, illustrated in Figure 1b, with strong release of heat that can easily cause a fire.

Handling catalytic air sensitivity requires specialized laboratory apparatus: vacuum lines, Schlenk lines, inert-atmosphere glove boxes, and special reactors. Most of the synthetically useful compounds, however, belong to the stoichiometric category and are frequently used reagents by the majority of chemists.

Extra-dry solvents are another class of air- and moisture-sensitive products, but they do not present intrinsic hazards. In their case, inappropriate handling and storage can cause the product to degrade over time. This is hard to detect but can interfere with experiments, affecting results and often requiring time-consuming troubleshooting.

Laboratory Safety Considerations

While the contact of an extra-dry solvent with moisture may not itself be a hazard, using a water-contaminated solvent in air-sensitive reactions can lead to very dangerous situations. Extra-dry solvents would for this reason be stored with the same care as pyrophoric compounds or solutions and other hazardous substances.

These considerations are very important because laboratory safety is a high priority for industrial and academic researchers alike. During the past fifteen years, there have been multiple

high-profile accidents, several of which have been summarily attributed to the inexperience of students working in academic labs. Accidents can rarely be attributed to a single cause, and further investigation of these incidents reveals several contributing factors.

There are no comprehensive datasets about the types and frequency of accidents and near-misses, which makes meaningful laboratory safety research difficult, according to Dana Menard and John F. Trant's review published in *Nature Chemistry* in 2020¹. Available data, however, show a high incidence of accidents involving reactive gases and pyrophoric materials. Butyllithium, for example, is known to cause problems in academic organic chemistry labs where correct management of air sensitivity can often be an issue.

The presence of more experienced research personnel in industrial and government laboratories may make these labs somewhat safer. However, a study by Schroder and associates² found that standardized risk assessments were not routinely used. Industrial and government labs reported higher compliance rates than the average in academic environments, but the average remains worryingly low at 43% and 36% respectively.

Safely Handling Air-Sensitive Chemicals

The use of air-sensitive reagents in chemical research is not new. A significant amount of literature is available to provide foundations for safe handling and operation³. An appropriate risk assessment and procedure must be based on the following two practices:

Content provided by:

ThermoFisher
SCIENTIFIC

- The use of clean and dry glassware and equipment
- The use of specialized packaging, syringes, and inert dry gases

While the first practice is straightforward, it must be taken seriously. Reagents can react violently in the presence of even tiny amounts of water, and minor air moisture condensation caused by temperature differences between labware and the environment can be enough to cause a fire. The second practice requires a different type of attention. Many air-sensitive reagents and ultra-dry solvents and solutions are specially packaged to help simplify their handling.

AcroSeal Packaging

The Acros Organics and Alfa Aesar brands offer AcroSeal packaging, the industry-leading solution. Its multi-layer septum and innovative anti-tampering cap (with ample surface area) provide convenient and safe storage and

dispensing while limiting exposure to the atmosphere. Syringes pass easily through the septum and create only a small hole that self-heals quickly under normal circumstances.

To dispense chemicals from AcroSeal packaging, use a syringe with an 18- to 21-gauge needle and a dry inert gas like nitrogen or argon. First, pressurize the bottle by injecting the gas and then withdraw the desired amount of liquid. Alternatively, use a double or double-tipped needle — one needle to withdraw the liquid and the other to add the inert gas from a gas line or balloon.

Best Practices

Some protocols recommend using glass syringes for pyrophoric reagents, but recent studies⁴ suggest that less

experienced users may find it easier to substitute a single-use polypropylene Luer lock syringe.

Drying solvents or preparing solutions of air-sensitive reagents can be tedious and time-consuming, and it often presents significant hazards. Ready-made products in specialized safe packaging provide a convenient, cost-effective, and safety-oriented solution.

The AcroSeal technology has undergone several innovations during its 25 years on the market. More than 2,000 Acros Organics and Alfa Aesar products are currently available in AcroSeal packaging.

Visit fishersci.com/acrosealpromo or fishersci.ca/acrosealpromo to save on select AcroSeal products. Visit fishersci.com/acroseal to shop the full line of products with AcroSeal packaging.



References:

1. A.D. Ménard and J.F. Trant, A review and critique of academic lab safety, *Nature Chemistry* 2020, 12, 17–25
2. Schröder, I., Huang, D. Y. Q., Ellis, O., Gibson, J. H. & Wayne, N. L. Laboratory safety attitudes and practices: A comparison of academic, government, and industry researchers. *J. Chem. Health Saf.* 2016, 23, 12–23
3. D.F. Shriver, M.A. Dredzdon, *The manipulation of air sensitive compounds*, John Wiley & Sons, 1986
4. E.S. Von Nehring, V. Dragojlovic, Handling of Air-Sensitive and Moisture-Sensitive Reagents in an undergraduate Chemistry Laboratory: The Importance of the Syringe, *J. Chem Ed.* 2021, 98, 246-249

Focus on Safety Webinar

Ansell's Integrated Approach to Improve Your Business Performance

Date: **Thursday, June 10, 2021**

Time: **1 p.m. ET**

Learn how AnsellGUARDIAN safety services can help you reduce personal protective equipment costs without sacrificing worker safety. Ansell experts can assess your work environment, identify hazards, and provide the safety solutions needed to maintain your operation while limiting hazards and improving worker performance.

Visit fishersci.com/ansellguardian-webinar to learn more.



Reproducibility Begins Here

MilliporeSigma

Milli-Q IX Pure Water Systems from MilliporeSigma

Milli-Q IX water purification systems produce consistently pure water that meets pharmacopeia and ISO requirements and provide safe water storage with continuous monitoring at the dispenser.

- Touchscreen controls
- Optimized water and energy use
- Mercury-free UV lamps
- Automatic e-record archives

Contact your **sales representative** or visit fishersci.com/MilliporeSigma for more information.



Product Description	Production Flow Rate	Cat. No.
Milli-Q IX 7003 Water Purification System	3 L/Hour	ZIX7003T0
Milli-Q IX 7005 Water Purification System	5 L/Hour	ZIX7005T0
Milli-Q IX 7010 Water Purification System	10 L/Hour	ZIX7010T0
Milli-Q IX 7015 Water Purification System	15 L/Hour	ZIX7015T0

Disinfect Your Workspace

Decon Labs Disinfectants

Decon Labs offers a variety of products effective against COVID-19, including ready-to-use disinfectant sprays, concentrated disinfectants, surface wipes, phenol-based disinfectants, quaternary ammonium, and sodium hypochlorite. These products are on EPA List N: disinfectants with efficacy.



Product Name	Minimum Contact Time
Conflikt	1 min.
BDD	10 min.
IRADECON	30 sec.
CiDecon Plus Wipes	5 min.
CiDecon II	10 min.

Product Description	Packaging	Quantity	Cat. No.
Conflikt: Ready-to-Use Quaternary Ammonium Disinfectant Surface Spray	Trigger-Spray Bottle, 16 oz.	12/Case	04-355-33
	Trigger-Spray Bottle, 32 oz.	6/Case	04-355-34
	Jug, 1 gal.	4/Case	04-355-52
BDD: Concentrated Quaternary Ammonium Disinfectant	Jug, 1 gal.	4/Case	18-800-101
IRADECON: Ready-to-Use Sodium Hypochlorite Spray	Trigger-Spray Bottle, 32 oz.	12/Case	04-355-150
CiDecon Plus Wipes: Phenolic-Based Disinfectant Surfaces Wipes	Canister, 180 Wipers	12/Case	04-355-87
CiDecon II: Concentrated Phenolic Disinfectant	Jug, 1 gal.	4/Case	04-355-129



Addressing Classroom Safety Concerns During the Pandemic

By Kylie Wolfe

To keep both students and staff safe, many universities turned to virtual learning this past year. Professors adapted their curriculum to help students get the most from each lesson while changing the way in-person science classes, including chemistry and biology labs, were offered.

Science courses typically require hands-on learning in the form of experiments and demonstrations. Not only do students learn valuable skills, like pipetting and problem solving, they also gain knowledge of lab safety techniques. With the obvious health and safety hurdles resulting from COVID-19, it was necessary for universities to find the right combination of in-person and online education.

Educating from a Distance

At the beginning of the pandemic, some college campuses went completely remote, but the 2020 school year brought about a hybrid model for many. Lab sizes were reduced and social distancing measures were put in place. Yale University installed clear dividers between lab benches for students attending in person and offered live video sessions for students participating online.

“In addition to limiting the number of live participants, we have rearranged our spaces, purchased additional equipment and reorganized our routines in order to minimize foot traffic, limit or eliminate sharing of equipment and maintain social distance at all times,” Christine DiMeglio, organic chemistry lab professor at Yale, told *Yale Daily News*.

Labs were well ventilated and staff disinfected surfaces regularly. They understood that online learning may not always be the best option, especially when it’s hard to recreate hands-on experiments that require a physical lab space. It was important to keep students safe while giving them a chance to learn in person as often as possible.

At the University of Pittsburgh, biology labs were limited to six students at a time. Chemistry labs also had reduced capacities and some physics labs were conducted at home. One approach paired an in-person student with an online student who could watch the experiment live over Zoom, regulating attendance in the classroom. This unique learning environment gave students flexibility while adhering to safety guidelines.

Protecting the Campus Community

Establishing proper protocols was necessary for laboratory courses as well as university research facilities. While some labs had shut down completely last year, others ramped up to help with COVID-19 studies. For Environmental Health and Radiation Safety (EHRS) offices, this was a time to act, managing traditional responsibilities alongside new obligations.

To bring students back to campus safely, EHRS offices played an essential role. At the University of Pennsylvania, the EHRS team helped educate and support students, staff, and the community. They led the on-campus contact tracing effort, reviewed safety plans, and installed signs and sanitization stations.

It was important to keep students safe while giving them a chance to learn in person as often as possible.

“When you get a lot of information and individuals are not used to dealing with certain terminology or approaches, that’s where we step in: To try to facilitate understanding, provide resources, and enable them to do what they need to do in short order and to help somebody get their job done,” said Andrew Maksymowych, associate director of biosafety at the University of Pennsylvania, in an interview with *Penn Today*.

Through consistent communication, COVID-19 safety sweeps, and overall compliance, the team helped the campus reopen. As the pandemic evolves and universities plan for the next school year, safety standards and new ways of holding lectures and labs will already be a tool in their toolbox. The ultimate goal will continue to be finding a balance between safety and education for students and staff.

This content was inspired, in part, by “From PPE to lab safety, supporting the campus community is a full-time job,” Penn Today, December 18, 2020; “COVID-19 brings a wave of change to in-person lab courses,” Yale Daily News, September 11, 2020; and “Students adjust to new online science lab structure,” The Pitt News, January 27, 2021.

Kylie Wolfe is a Thermo Fisher Scientific staff writer.



5,000 solutions. 1 washer.



We hear you, keeping up with cleaning all of your labware is a challenge. Labconco's glassware washers are up for it, with the ability to wash over 5,000 different labware configurations. New software lets you store, customize, and protect programs and data. Labconco. Let us clean, so you don't have to.

Learn more at fishersci.com

Search: [FlaskScrubber](#) or [SteamScrubber](#)

LABCONCO®

The Cost of Risk May Be Defined Only After an Accident Occurs

We all share the common experience of breaking a wine glass or a coffee cup. But in the lab, damage from broken glass can have much greater consequences than spilling your favorite merlot on your new shoes.

For example, if a flask breaks, there is risk of getting cut. And if the flask contains dangerous contents, that risk extends to nearby personnel through splash, vapors, fire risk, equipment damage, and more.

You can better address these risks by understanding the various types of glass and their use, care, and safe disposal practices. It is recommended that you implement a proper glass safety standard operating procedure (SOP) to strengthen lab safety, mitigate risk, and help avoid financial loss. Unfortunately, consideration of financial loss sometimes occurs only after an accident occurs. So how do you get started?

Regulatory Standards

A major difference between consumer and laboratory glassware involves standards, regulations, and practices for manufacture and use. Lab glassware has purity and calibration requirements, should be handled to reduce thermal stress, and requires proper cleaning. In regulated labs, your records of the physical inspection of glassware for proper function and conformance will be reviewed by inspectors.

ASTM International (formerly the American Society of Testing & Materials) is a global organization that helps establish standards for the manufacture of laboratory glassware, including materials, testing, calibration, and other manufacturing procedures. For example, ASTM Method E671 sets limits on acceptable thermal stress levels for annealed glass and offers guidelines for testing, calibration, and maintenance.

Regulatory authorities that follow ASTM standards include the U.S. Environmental Protection Agency (EPA), the U.S. Department of Agriculture (USDA), the U.S. Food and Drug Administration (FDA), the AOAC (formerly the Association of Official Agricultural Chemists), and the U.S. National Institute of Standards and Technology (NIST), among others.

Make sure the manufacturer of your lab glassware complies with these standards for safety, purity, and performance.

Selection of Glass Types in the Lab

In the laboratory, soda lime and borosilicate glass predominate. Soda lime glass is used mostly for jars and containers, and it cannot be autoclaved (sterilized). Borosilicate glass is formulated to withstand thermal shock, can be autoclaved, and resists breakage.

Awareness and understanding of glass differences is critical. Analytical and chromatography glassware should be either Type 33 or 51 borosilicate glass. Never autoclave glass unless the manufacturer indicates it is safe to do so.

Application Considerations

Factors such as chemical attack, sample-glass interactions, thickness, content photosensitivity, washing solvents, and calibration are also important when choosing the right glass type for a procedure.

Life Cycle Management

Life cycle management includes the regular inspection and proper maintenance of reusable glassware. The human eye cannot detect some indications of stress — by the time scratches, scrapes, and damage are visible, it's too late.

Alan Arnold Griffith, an English engineer, studied glass stress and fracture. In the 1920s, he discovered what are known as Griffith Flaws, which are introduced by the mechanical stress of day-to-day laboratory use and improper uses like scratching, flexing, bending, and dropping.

Use a polariscope periodically to inspect your glassware for signs of damage. Inspect and dispose of stressed glass to help prevent injuries and collateral damage.

You may also opt to simply replace your glassware on a scheduled basis. If you use serialized or color-coded glassware, you can track usage, or change colors after a span of two or three years. For maximum risk

prevention, use safety-coated bottles that help contain the contents if the glass breaks.

Ask these questions before an accident occurs:

- What risks are associated with routine daily use or improper cleaning or handling of lab glassware?
- Are we reducing or increasing the creation of Griffith Flaws during use?
- Do our glass usage SOPs address replacement, disposal, or calibrations?

A well-designed and implemented SOP limits your financial loss, improves lab operations, and mitigates safety and regulatory risks.

Now that you understand more about the risks, take action! Try the free DWK Life Sciences on-site or webinar-based Glassware Safety Training. Contact your Fisher Scientific sales representative to schedule a training session.

Visit fishersci.com/glassafety to learn more.

Content provided by:



The Unmatched Versatility of Tyvek Protective Apparel

By Daniel Hammel, North America Marketing Lead, Chemical Industrial & Controlled Environments, DuPont Personal Protection

Both in and out of the laboratory, DuPont Tyvek garments have long been the first choice in worker protection. Much lighter than spunbond meltblown spunbond (SMS) and microporous film (MPF), Tyvek flash-spun fabric provides a barrier to small hazardous particles. This protection is built into the fabric itself, with no films or laminates to abrade or wear away. Abrasion tests show that Tyvek barriers are more durable than MPF fabrics when tested for hydrostatic resistance.



Garments and Accessories

A breathable barrier against hazardous dry particles (as small as one micron), aerosols, and non-hazardous light liquid splash, Tyvek 400 protective apparel can help minimize risk when used in combination with other personal protective equipment (PPE) in hazardous situations.

Many options are available for garments and accessories, including coveralls, lab coats and frocks, sleeves, hoods, shirts, pants, and aprons. Boot and shoe covers are produced with Tyvek 400 or Tyvek 400 FC, the latter adding a friction coating for better skid resistance.

Tyvek is versatile — other styles and innovations have been developed to

provide solutions to growing industrial requirements and end-user needs.

Controlled Environments

Laboratories, cleanrooms, and controlled environment settings — like pharmaceutical, medical device manufacturing, biotech, and electronics manufacturing sites — have high standards for particle and microbiological contamination control.

Tyvek IsoClean garments are part of the DuPont controlled environments product line of single-use apparel and accessories. Tyvek IsoClean products have a long history of cleanroom use and offer an excellent barrier to particles, microorganisms, and non-hazardous liquids.

Find Tyvek IsoClean garments and accessories in clean-processed and sterile options, and in non-sterile bulk packaging. Tyvek Micro-Clean 2-1-2 coveralls are individually packaged as sterile or non-sterile. Sterile Tyvek apparel is gamma sterilized to a Sterility Assurance Level (SAL) of 10^{-6} and is fully traceable.

Available in white and blue coverall options, Tyvek apparel allows you to designate workers who are allowed in restricted areas or can perform certain tasks in your laboratory or controlled environment. Or you can use color coding to easily identify team leaders or supervisors.

Enhanced Chemical Resistance

Tyvek 800 is a liquid-tight garment that combines comfort and durability with

protection against low-concentration, water-based inorganic chemicals and particles as small as one micron.

Tyvek 800 garments even stand up to pressurized jets of liquids and are suitable for a variety of industrial cleaning, food manufacturing, agriculture, petrochemical, wastewater, environmental remediation, laboratory, and maintenance applications.

Tyvek 500 protects against light liquid aerosols and airborne solid particles. Tyvek 600 offers a chemical permeation barrier to low-concentration water-based inorganic chemicals, including infectious agents and bodily fluids.

Tyvek 500, Tyvek 600, and Tyvek 800 are tested and classified according to EN 14126.

Breathability, Form, and Fit

Your body generates heat and moisture vapor whether you're working hard or standing still, but some garments don't



Content provided by:



“breathe.” Trapped heat and moisture can lead to heat stress that not only makes you uncomfortable but also compromises safety. Based on moisture vapor transmission rate (MVTR), Tyvek 400 fabric is 30% more breathable than typical microporous film fabrics.

Tyvek 400 D garments suit applications with limited frontal exposures. These include exposure to tacky resin molds that may stick to (or delaminate) other

garments and sanding and grinding operations that release fine airborne particulates.

These comfortable garments combine the protection and durability of Tyvek fabric on the front with the softness and breathability of DuPont ProShield 10 fabric on the back.

The Need to Be Seen

Low-visibility hazards are common in construction, transportation, waste handling, mining, and other industries, and high-visibility protective apparel lets workers be seen more easily.

With Tyvek 500 HV coveralls, workers benefit from breathable and durable high-visibility garments with the added particle protection of Tyvek. Suitable for the dirty jobs common in many

industries, Tyvek 500 HV coveralls also helps ensure compliance with PPE requirements.

The bright fluorescent orange color of the Tyvek background fabric combined with silver-gray retroreflective material makes these garments highly visible during the day or at night when exposed to a light source. The reflective symmetric “X” on the back of the coverall helps distinguish the worker’s front and back.

Making a Difference

Tyvek offers unmatched versatility, superior overall barrier protection, improved breathability, and better garment durability for laboratory applications, pharmaceutical handling, general manufacturing, and dirty jobs.



Precise Heating and Stirring

OHAUS Guardian Hotplate Stirrers

Durable OHAUS Guardian Hotplate Stirrers deliver intelligent performance and enhanced safety features*:

- SmartHousing is cool to the touch and chemical resistant
- SafetyHeat internal protection system monitors heating and unit function
- SmartRate temperature selection enables fast or slow ramp rate speed
- SmartHeat controls and monitors maximum temperature to prevent overheating
- SmartLink turns off the heating function if the Bluetooth connection is broken
- SmartPresence proximity sensor technology turns off the unit when no one is nearby



Description	Surface Dimensions	Surface Material	Cat. No.
Hotplates			
Guardian 5000 Hotplates	7 × 7 in.	Ceramic	01-922-000
Stirrers			
Guardian 5000 Stirrers	7 × 7 in.	Ceramic	01-922-001
Stirring Hotplates			
	7 × 7 in.	Ceramic	01-922-002
Guardian 5000 Hotplate Stirrers	10 × 10 in.	Ceramic	01-922-003
	Dia. 5.3 in.	Aluminum	01-922-004
	7 × 7 in.	Ceramic	01-922-005
Guardian 7000 Hotplate Stirrers	10 × 10 in.	Ceramic	01-922-006
	Dia. 5.3 in.	Aluminum	01-922-007

*Guardian 5000 units have SmartHousing and SafetyHeat features only; Guardian 7000 units have all six safety features.

fisherbrand

Direct and Indirect Sonication

Fisherbrand Sonic Dismembrators

Fisherbrand Sonic Dismembrators create energy that is transmitted through a titanium probe into a liquid sample to create cavitation (the implosion of micro-bubbles with high shear forces).

Each Fisherbrand Sonic Dismembrator includes a generator, converter, cables, wrench set, and one probe. Other probes and accessories, including the stand and clamp shown here, are sold separately.



Model 505

Model	Applications	Capacity	Power	Cat. No.
50	• Basic Cell Disruption	0.2 to 50 mL	50w	FB50110
120	• Cell Disruption • Protein Extraction • DNA Shearing/ChIP	0.2 to 50 mL	120w	FB120110
505	• Cell Disruption • Nanoparticle Dispersion • Homogenization/Mixing	0.2 to 1000 mL	500w	FB505110
705	• Cell Disruption • Protein Extraction • DNA Shearing/ChIP • Nanoparticle Dispersion • Homogenization/Mixing • Sonochemistry	0.2 to 1000 mL	700w	FB705110

Faster Real-Time PCR Results

Maximize your PCR workflow when you pair Applied Biosystems TaqMan Fast Advanced Master Mix and Array Plates with Invitrogen SuperScript IV VILO Reverse Transcriptase. These real-time PCR reagents work together to give you fast and accurate data and help accelerate your research.

Invitrogen SuperScript IV VILO Master Mix

Choose Invitrogen SuperScript IV VILO Master Mix, a first-strand cDNA synthesis reaction mix, for your two-step RT-qPCR. This product combines optimized buffer conditions with thermostable reverse transcriptase (RT) to maintain linearity across a broad range of RNA inputs.

- Efficient: Up to 100 times higher cDNA yield
- Sensitive: Reduce as many as eight RT-qPCR cycles
- Robust: Works with degraded RNA samples or those that contain inhibitors
- Fast: Ten-minute reaction times



Quantity	Cat. No.
50 Reactions	11-756-050
500 Reactions	11-756-500

Applied Biosystems TaqMan Fast Advanced Master Mix

Designed to match or exceed standard master mix performance, Applied Biosystems TaqMan Fast Advanced Master Mix delivers results in less than 40 minutes.

- Enhanced benchtop stability
- Reduced run times (fast or standard instruments)
- Optimized for multiplexing
- Seamless transitions into your protocols



Quantity	Cat. No.
1 x 1 mL	44-445-56
1 x 5 mL	44-445-57
1 x 50 mL	44-445-58
2 x 5 mL	44-449-63
5 x 5 mL	44-449-64
10 x 5 mL	44-449-65

Applied Biosystems TaqMan qPCR Array Plates

Measure pathway- and disease-related gene expression with Applied Biosystems TaqMan qPCR Array Plates (pre-spotted panels of TaqMan Real-Time PCR Assays). These arrays feature the most popular TaqMan gene expression assays targeting relevant genes. They're ideal for profiling and verification applications that require the analysis of thousands of targets.



**Save 10%
When You
Buy All 3**

When you purchase the following products, you'll save 10% on your order.

- Applied Biosystems TaqMan Fast Advanced Master Mix (5 mL)
- Invitrogen SuperScript IV VILO Master Mix (50 Reactions)
- Applied Biosystems Predesigned TaqMan Array Plates

Visit fishersci.com/fast-taqman to learn more.

Restrictions apply. See back cover for details.

Visit fishersci.com or fishersci.ca to learn more about these products.



Reducing the Environmental Impact of Disposable Face Masks

By Christina P. Hooton

Since the beginning of the pandemic, disposable face masks have become a part of everyday life as we know it. It's estimated that 129 billion face masks are used globally every month.¹ While these are a proven and necessary way to curb the spread of COVID-19, according to data shared by the Centers for Disease Control and Prevention (CDC), they are also adding to the already significant amount of plastic pollution throughout the world.^{2,3}

A Plastic Problem

Although most disposable face masks appear to be made from cloth, they're usually produced from non-woven polymer fabrics made of polypropylene or polycarbonate. The typical disposable blue mask has three layers: an outer layer of non-woven fibers, a filtering middle layer, and a soft inner layer. The filter in the middle is made by forcing melted polymer through tiny nozzles using high-speed blowing gas.

As these disposable masks make their way into the environment, they sometimes pollute waterways or are disposed of in landfills. Like other plastics, they can also break down into microplastics, particles smaller than 5 mm.

Plastics and microplastics are a known threat to the health of marine species, a large part of our food supply. They also enter the food we eat and significantly contribute to climate change through carbon emissions. And, in a rather perverse twist, these plastic particles are also known to spread microbes such as invasive pathogens.⁴

Other Considerations

One solution to this growing problem might be for everyone to wear reusable instead of disposable masks, especially for personal activities like grocery shopping. If a reusable cloth mask offers the right level of protection for a given situation, it must be washed at least daily, and several must be on hand so one is always ready to wear.

However, according to CDC data, a cloth mask alone blocks only 51.4 percent of airborne cough particles compared to properly worn surgical masks that block 77 percent of airborne particles. Wearing both a cloth and a surgical mask proves even more effective, blocking 85.4 percent of particles.⁵ Laboratories and other workplace settings may require the level of protection these disposable face masks provide.

Face Mask Recycling

It may seem obvious to dispose of your face mask along with all other recyclable materials. However, face masks and other forms of single-use personal protective equipment (PPE) are not easy to recycle and may be too expensive for local recyclers to collect. It is not likely that your current waste management solution is offering this recycling service.

It's estimated that 129 billion face masks are used every month globally.¹

The world needs an economical, efficient solution for recycling disposable face masks as safety restrictions from the pandemic persist. Face mask recycling programs present an attractive option especially for organizations like universities and colleges trying to establish and maintain sustainability targets.



continued from page 21

Reducing the Environmental Impact of Disposable Face Masks

TerraCycle, a company known for its ability to recycle challenging materials, offers a mask recycling program through the Fisher Scientific channel. To start, an institution purchases the TerraCycle Disposable Masks Zero Waste Box, a cardboard box clearly marked for face mask disposal and recycling, and places it in their facility.

“It’s a very simple process. The box is shipped direct to the customer. They assemble the box and then they start collecting,” said Madeleine Chadwick, General Management, Leadership Development Program, Thermo Fisher Scientific.

The boxes come in three sizes with the largest option accommodating up to 2,100 masks. Once the box is full, the institution seals and ships it back to TerraCycle where face masks are recycled and processed into new materials. Chadwick said, “The price of the box includes shipping costs and the recycling service.”

Although the program is not intended to recycle conventional medical waste, the boxes are quarantined upon arrival for a minimum of 72 hours to protect the employees that handle the waste. Items are hand-sorted, and any visibly contaminated masks are removed and forwarded to third-party facilities for processing and recycling.

Turning Trash into Treasure

The mask waste is collected until there is enough volume to perform processing, then sorted based on material composition. Metal wires or nosepieces are removed from the masks, grouped together, and processed separately, for example.

The materials are then sent for third-party processing where they are recycled into usable forms like plastic pellets and granules. The polypropylene mixture typically found in face masks is turned into a raw material used in plastic lumber and composite decking applications. That material is then sold to manufacturers for the creation of a variety of new products, including outdoor furniture, shipping pallets, and more.

It’s encouraging to think that some of these recycled face masks could one day become something useful like park benches or picnic tables that people will freely gather on and around once pandemic restrictions are lifted. Fisher Scientific customers alone have already diverted more than 735,000 face masks from landfills using this program.

The TerraCycle program offers organizations a flexible, easy-to-set-up solution. Boxes can be placed in high-traffic areas along with hand sanitizers to reinforce both safety and sustainability. They can be especially useful in areas where a new disposable face mask is required upon every entry, for example outside a research lab. The program was recently expanded to include boxes for gloves, eyewear, and other hard-to-recycle items.

“It’s a very simple process. The box is shipped direct to the customer. They assemble the box and then they start collecting.”

It’s unclear how much longer we’ll be wearing face masks. Even after the pandemic, face mask use might be more prevalent than before, especially during periods like flu season or when traveling. One thing that’s absolutely clear is the need to understand and mitigate their environmental impact, one face mask at a time.

Visit fishersci.com/terracycle-mask-box or fishersci.com/terracycle-glove-box for more information about TerraCycle recycling opportunities.

References:

1. Prata, J.C., Silva, A.L.P., Walker T.R., Duarte, A.C., and Rocha-Santos, T. (2020). COVID-19 Pandemic Repercussions on the Use and Management of Plastics. *Environ. Sci. Technol.* 54(13), 7760–7765.
2. Centers for Disease Control and Prevention. (2020, November 20). *Science Brief: Community Use of Cloth Masks to Control the Spread of SARS-CoV-2*. <https://www.cdc.gov/coronavirus/2019-ncov/more/masking-science-sars-cov2.html>
3. Schnurr R.E.J., Alboiu V., Chaudhary M., Corbett R.A., Quanz M.E., Sankar K., Strain H.S., Thavarajah V., Xanthos D., Walker T.R. (2018). Reducing marine pollution from single-use plastics (SUPs): a review. *Mar. Pollut. Bull.* 137:157–171.
4. Fadare, O. O., & Okofio, E. D. (2020). Covid-19 face masks: A potential source of microplastic fibers in the environment. *Science of The Total Environment*, 737, 140279. doi:10.1016/j.scitotenv.2020.140279
5. Centers for Disease Control and Prevention. (2021, February 19). *Maximizing Fit for Cloth and Medical Procedure Masks to Improve Performance and Reduce SARS-CoV-2 Transmission and Exposure, 2021*. <https://www.cdc.gov/mmwr/volumes/70/wr/mm7007e1.htm>

Christina P. Hooton is a Thermo Fisher Scientific staff writer.

Preparing Samples for Long-Term Cryogenic Storage

How can controlled-rate freezers be used to prepare samples and therapies for long-term cryogenic storage?

The term “cryopreservation” is derived from the Greek word *kryos*, which means cold. Practically, cryopreservation is a method for preserving biological samples in a vitrified state.

Rapid freezing prevents the formation of ice crystals by chilling the sample to the amorphous ice phase while circumnavigating the hexagonal ice phase. Because cryopreservation temperatures are so low, the water inside cells in a sample is vitrified or removed by osmosis and stopping metabolic processes.

Advantages of sample cryopreservation:

- Reduced risk of microbial contamination or cross-contamination
- Less risk of morphological and/or genetic changes
- Techniques are well established and methods and outcomes are well documented

Cryopreservation is effective but not trouble free. Researchers must determine optimal conditions for the viability and recovery of their samples.

What Are the Freezing Phases?

ISBER guidelines state that LN₂ cryopreservation is optimal for long-term specimen storage. These systems offer storage below the glass transition phase temperature of -132°C. In this case, the term “glass” refers to a liquid that is too cold to flow.

Nucleation is the onset of the change from liquid to crystalline. It involves the progressive separation of ice from the unfrozen part of the solution. The solute concentration increases as water separates as ice. This change is also associated with an energy change. The localized increase in temperature is called the latent heat of fusion, which evolves as the solution reaches the equilibrium freezing point.

Best Practices

While storage below the glass transition temperature is important, properly preparing samples for storage at very low temperatures is equally important. When the rate at which a specimen is frozen is controlled, it fosters greater cellular viability when the sample is removed from its cryogenically frozen state.

The International Society for Biological and Environmental Repositories (ISBER)

publishes *The ISBER Best Practices: Recommendations for Biorepositories Fourth Edition*, which offers guidelines for sample preservation. In 2019, the *ISBER Best Practices Addendum: Liquid Nitrogen-Based Cryogenic Storage of Specimens* was also published.

These guidelines comprise the most effective evidence- or consensus-based practices for managing biological and environmental specimen collection, long-term storage, retrieval, and distribution and support the availability of high-quality research specimens.

Thermo Scientific CryoMed Controlled-Rate Freezers

Controlled-rate freezers can help reduce the stress on samples during the freezing and vitrification processes.

CryoMed Controlled-Rate Freezers have six factory-set freezing profiles and users can create up to 14 custom profiles. For samples and therapies, CryoMed users can take advantage of a new pre-set profile that has an end-temperature of -140°C.

This profile, labeled Pre-Set Profile #3, controls sample freezing through the glass transition phase (when amorphous ice can form). It creates a 1°C/minute



Content provided by:

ThermoFisher
SCIENTIFIC

cooling rate from nucleation to -45°C and a $10^{\circ}\text{C}/\text{minute}$ cooling rate to the final temperature of -140°C . It is commonly used for hemopoietic stem cells in volumes from 65 to 100 mL.

Long-Term Sample Storage

Choose from two LN_2 cryogenic freezer platforms for long-term sample storage:

- Thermo Scientific CryoExtra High-Efficiency LN_2 Storage Tanks
- Thermo Scientific CryoPlus LN_2 Storage Devices

To choose the right platform, consider the number of samples being stored and whether liquid or vapor phase liquid nitrogen will be used.

- CryoExtra is preferred for larger numbers of samples stored in vapor phase; four capacity options accommodate between 19,500 and 93,000 vials (2 mL)
- CryoPlus is preferred for fewer samples and accommodates both liquid or vapor phase storage; 2 mL vial capacity ranges from 6,318 to 39,000

Thermo Scientific CryoMed Controlled-Rate Freezers, used in combination with CryoExtra or CryoPlus platforms, provide a workflow solution for laboratories or biorepositories with reliable sample preparation and cryogenic storage.

Visit fishersci.com/CryoMed or fishersci.ca/CryoMed to learn more.



Modern Antibody Analysis and Testing



Greiner Bio-One ELISA Microplates

Why Greiner Bio-One?

- ELISA standard microplate and microplate strip formats
- Consistent and reproducible protein binding
- Constant and uniform binding properties guaranteed



Description	Binding	Well Shape	Mfr. No.	Cat. No.
96-Well MICROLON 200 Plate	Medium	Round Bottom	650001	07-000-131
96-Well MICROLON 600 Plate	High	Round Bottom	650061	07-000-608
96-Well MICROLON 200 Plate	Medium	Flat Bottom	655001	07-000-100
96-Well MICROLON 600 Plate	High	Flat Bottom	655061	07-000-102
96-Well MICROLON 200 Plate, Chimney Well	Medium	Flat Bottom	655080	07-000-099
96-Well MICROLON 600 Plate, Chimney Well	High	Flat Bottom	655081	07-000-101
96-Well MICROLON 200 Strip Plate, 1 x 8 Single Break	Medium	C-Bottom	705070	07-000-779
96-Well MICROLON 600 Strip Plate, 1 x 8 Single Break	High	C-Bottom	705071	07-000-780

BULK UNSTERILE SYRINGES

For Lab and Industry Use

Air-Tite's bulk, unsterile syringes are our most economical non-sterile option. The general-purpose, three-part syringes are designed for applications that require a simple, low-cost transfer vessel. Made to our exacting specifications, these syringes are latex free, non-pyrogenic, and non-toxic. Along with high-quality markings, the design includes siliconized gaskets for a consistently smooth plunger motion and positive plunger stop.



14-817-117

Luer Lock and Luer Slip Syringes available in 1 mL, 3 mL, 5 mL, 10 mL, 20 mL, 30 mL, and 60 mL sizes.

Sterile options are also available.

Contact your Fisher Scientific sales representative today for more info.

Visit fishersci.com/air-tite or fishersci.ca/air-tite to shop.



Air-Tite
Products Co., Inc.

Cytiva Kits and Columns

Cytiva Sera-Xtracta Virus/Pathogen Kits

Isolate total nucleic acid (DNA/RNA) from bacteria and viruses, including adenovirus (type 14), influenza A (H3N2), and COVID-19, with Cytiva Sera-Xtracta Virus/Pathogen Kits. Designed for sensitive detection at low concentrations, this simple and rapid method optimizes your workflow and helps you achieve reproducible results.

Quantity	Cat. No.
96 Extractions	09-920-080
1,000 Extractions	09-920-081



Cytiva Superdex 200 Prepacked Columns

The Cytiva Superdex 200 Increase is a next-generation size exclusion chromatography (SEC) resin used for high-resolution purification and analysis of monoclonal antibodies (mAbs) and proteins. It's available pre-packed in three sizes to meet the needs of various applications. With smaller and more rigid beads than their predecessors, Superdex Increase columns deliver improved resolution in shorter run times.

Description	Cat. No.
Superdex 200 Increase 10/300 GL Prepacked Columns	45-002-570
Superdex 200 Increase 5/150 GL Prepacked Columns	45-002-571
Superdex 200 Increase 3.2/300 Prepacked Columns	45-002-572



Visit fishersci.com/cytiva or fishersci.ca/cytiva to shop.

Nitrile Allergies? Accelerator-Free Nitrile Can Help

by Fritz Maskrey, Vice President, Business Development, TechNiGlove International

When thinking of allergies in the disposable glove world, natural rubber immediately comes to mind.

A simple solution for this problem was introduced in the 1990s using synthetic polymers. This new material would mimic the properties of natural rubber without the protein allergens that were causing some users to develop allergic reactions. Hence the product we know today as “thin-walled disposable nitrile gloves” was born.

Nitrile Glove Production: Simple but Complicated

To create nitrile products, liquid nitrile is refined and coated onto a ceramic form then cured to create an elastic polymer in the desired shape. The elastic polymer, which has been molded into the shape of varying sizes and lengths of gloves, is then chlorinated inside and out and finished based on the specifications of the intended cleanroom environment or general laboratory use.

Since the inception of nitrile, manufacturers have been competing to create a better and more efficient process for making gloves. The competition has led manufacturers to add certain accelerators intended to decrease the time needed for a glove to cure, reduce the cost of the overall material, give the material a softer feel, and generally provide the customer with an affordable high-performance product.

Certain sulfur-based chemicals have been developed to accomplish these goals. These chemicals have incredibly complicated chemical structures and even more complicated scientific names. While this list of chemicals can have many positive effects, one unfortunate side effect can be an allergic reaction (now known as “nitrile allergies”) similar to that of latex.

Current examples of accelerators used in nitrile glove production include:

- N-cyclohexyl-2-benzothiazole sulfenamide

- 2-Dibenzothiazole disulfide
- Tetramethylthiuram disulfide
- Diphenylguanidine
- 2-Mercaptobenzothiazole
- 3-Methylpiperidine

Does This Mean I’m Allergic to Nitrile Gloves?

Since most of these reactions are not to nitrile itself but to chemicals added during manufacturing, they are considered “Type 4” allergies. Also called cell-mediated or delayed allergies, they occur when T-cells are activated and produce inflammation in the affected area. Type 4 allergies are not antibody related but rather are a type of white blood cell response.

One of the newer technologies used in the production of thin-walled disposable gloves is “accelerator-free” nitrile. Production managers, line engineers, and corporate purchasers should be aware of this option. Choosing gloves made



Content provided by:



with this material and these methods can drastically reduce reactions for users with Type 4 allergies.

Unlike latex allergies, which can cause more life-threatening allergic reactions, people who react to the accelerators in nitrile gloves usually experience only contact dermatitis (Type 4 Hypersensitivity).

Reactions to nitrile gloves may include hives, blisters, itching, burning, and sensitivity to the sun. These symptoms usually appear one or two days after exposure. Nitrile glove wearers can also experience non-allergic irritant contact

dermatitis, the symptoms of which include blisters, lesions, ulcers, and dry skin.

Techniglove's Nitrile and Accelerator-Free Nitrile Gloves

Key attributes of our accelerator-free products include:

- Elongation and tensile strength over 300% stronger than traditional nitrile formulations
- Fit, feel, and function almost identical to latex gloves
- Less material with superior strength

Contact your Fisher Scientific sales representative for information about a wide selection of nitrile gloves including Rival (RV400 Series) accelerator-free, low extractable, powder-free nitrile examination gloves.

TechNiGlove International is a leading manufacturer of disposable gloves for contamination-controlled environments. TechNiGlove's nitrile, sterile nitrile, and latex gloves offer consistent quality for use in cleanrooms, sterile non-medical environments, pharmaceutical environments, and industrial settings.





The unmatched versatility of Tyvek in & out of the laboratory



DuPont™ Tyvek® 400 garments are composed of 100% flashspun high-density polyethylene, creating a unique, nonwoven material that provides the ideal balance of protection, durability, and comfort.



Tyvek® IsoClean® garments are offered in clean-processed and sterile options, as well as non-sterile bulk-packaged options. Tyvek® IsoClean® sterile garments are gamma sterilized to a sterility assurance level (SAL) of 10^{-6} and are fully traceable.



Bead Mill Basics

Content provided by:

fisherbrand

In the science community, time is of the essence and results are critical. Efficiency is crucial when lab personnel need to extract different analytes from a variety of samples.

Sample preparation is especially important with problematic materials or specimens and may create a bottleneck in the overall extraction process. Traditionally, samples were individually and manually processed, sometimes ground using a mortar and pestle or subjected to long enzymatic digestions. Bead mill homogenizers help speed this process up and keep your work on schedule.

How Do Bead Mills Work?

Bead mills are one of the most effective methods for processing samples into smaller particle sizes and can be used to process multiple samples simultaneously.

Bead mills are used for dispersion processing — they reduce the size of solid particles and disperse them into the sample. The beads are placed inside a tube with the sample and agitated vigorously. Homogenization occurs with impact or collisions between the sample and the beads (macro-scale), between the beads and the tube (micro-scale), and shear forces created by rapid bead movement (macro-scale).

Choosing the Right Beads

The type and size of the beads are important factors when processing samples. Larger beads can break down large or dense structures while small beads are more effective with cellular components.

Bead shape and density matter as well. Denser beads will be more effective on tougher or harder materials. Irregularly shaped beads create greater mechanical shear forces and work well with tough or fibrous materials.

The variation in beads makes bead mills suitable for multiple applications. Bead mills can be used to dissociate tissue and recover live cells or extract small molecules like nucleic acids and proteins. Bead mill tubes do not require additional liquid, so they can be used for milling or dry grinding solid samples.

Advantages of Bead Mill Sample Prep

Bead mill homogenizer tubes are disposable — single use reduces the risk of cross-contamination. And because homogenization occurs inside a closed tube, there is less chance of creating aerosols or other user exposure. This is especially important if samples contain toxic, infectious, or otherwise hazardous substances.

Use Fisherbrand Bead Mill Homogenizers for grinding, lysing, or homogenizing biological samples for molecular extractions. A unique tube carriage motion provides efficient bead movement, producing impact forces that match or exceed other bead mills in the market. And increased power reduces processing times to maximize your sample prep efficiency.

- Simultaneously homogenize up to 24 samples
- Front-loading tube holder for convenience, ease of use, and optimal bead/sample interactions
- Non-stop processing — no need to cool down between runs
- Rapid processing means less sample degradation
- Timer range: 1 to 99 seconds
- Performance range: 0.8 m/s to 6 m/s in increments of 1.5 m/s
- Cycles: 1 to 9
- Small footprint

Fisherbrand Bead Mill Homogenizers, coupled with sample-specific bead beating materials, can help you produce thorough homogenates regardless of the sample. Compared to conventional methods, bead mill homogenization increases the yields of nucleic acids, proteins, and other small molecules to improve the results of your downstream assays.



Fisherbrand Bead Mill 4 Homogenizer
Cat. No. 15-340-164



Fisherbrand Bead Mill 24 Homogenizer
Cat. No. 15-340-163



The Basics of Cleanroom Upkeep

By Mike Howie

Cleanroom upkeep is essential to protecting the quality of your finished products. In some cases, even one small contaminant can be destructive, potentially costing millions of dollars in wasted effort and material.

Every cleanroom is unique — built for a specific purpose and subject to specific guidelines. Your cleanroom upkeep routine should be just as unique, designed with your physical space and applications in mind. There are, however, general guidelines that can help you build the right routine for your needs. Broadly, there are four parts to cleanroom upkeep:

- Prevent as much contamination from entering the cleanroom as possible
- Clean as well as you can using the right tools for your environment
- Monitor surfaces and the air for contaminants
- Validate that your routine is keeping the environment adequately clean

Here are just a few tips for keeping your cleanroom clean.

Keep Contaminants Out

If you start with a clean cleanroom, upkeep is much simpler.

Anything entering the cleanroom, including people, is the largest source of potential contamination. Restricting items like cardboard and using sticky mats to remove particles from shoes can help prevent contaminants from entering the cleanroom in the first place. Take care when introducing any new people or items to the cleanroom. People should be wearing the right garments and donning them properly.

Use the Right Cleaning Products

Your choice of cleaning products is important because lower-quality products can introduce more contamination into your cleanroom. Keep in mind that there are always new mops, wipers, disinfectants, and other products entering the market and what you've been using may no longer be in line with best practices. String mops, for example, effectively spread disinfectant but don't reliably remove contaminants. A cleanroom mop should help disinfect surfaces and successfully remove contaminants.

While these cleaning products can be more expensive, ask yourself how much you're willing to pay for peace of mind. If a

0.5 µm particle could ruin a 1 million USD product, purchasing effective cleaning products is likely worth the investment.

Be Wary of Laundering Garments and Supplies

In some situations, laundering and reusing garments, mops, or other supplies can effectively reduce costs. But it also carries a risk: laundering may not remove all contaminants. If contaminants remain in a cleaning product after laundering, there's a risk that those contaminants could be reintroduced into the cleanroom. While using new cleaning products every time may be more expensive, it helps reduce the risk of introducing contaminants while cleaning.

Clean Before You Disinfect

Sterile is not the same as clean. While disinfecting may render contaminants inert, sterile contaminants can still adulterate your final product. By removing contaminants before disinfecting, you can be more confident that your cleanroom is truly clean.

Remain Vigilant

Continuously monitor surfaces and the air in your cleanroom for particles, viables, volatile organic compounds, and any other potential contaminants. By tracking conditions within your cleanroom, you'll have a better idea of the challenges you face and can evolve your cleaning procedures as needed.

After a while, you can analyze the collected data to identify any trends. Any increase in contamination is worthy of investigation — even if your cleanroom is still within defined limits.

Follow Applicable Guidelines

Multiple organizations publish guidelines for building, operating, and maintaining cleanrooms, including the Institute of Environmental Sciences and Technology (IEST), the Parenteral Drug Association (PDA), ASTM International, and the United States Pharmacopeia (USP). In some instances, these organizations even provide detailed direction on how to clean. Research which guidelines apply to your cleanroom and follow their guidance.

Cleanroom upkeep should be a constant priority and a daily consideration. By keeping your cleanroom as clean as possible at all times, you're protecting the quality of your finished products.

Mike Howie is a Thermo Fisher Scientific staff writer.

Curbside Recyclable Packaging



SONOCO ThermoSafe EOS Shippers


Make sure your temperature-controlled packaging waste won't end up in a landfill.

The outer and payload boxes of new Sonoco ThermoSafe EOS shippers are made from corrugate, using post-consumer recycled content. The EOS product offers high insulation, has cushioning to protect your precious samples and products, and is curbside recyclable. Just place the entire knock-down package into your recycling bin when you're finished.

Sonoco ThermoSafe is serious about sustainability and works hard to develop products that help minimize environmental impact while maintaining the quality and reliable transport conditions you've come to trust.



Introducing an elastomeric rubber reusable respirator WITHOUT an exhalation valve!



- More comfortable than wearing disposable style
- Two-way filtration when you breathe in and out
- Two sizes
- NIOSH approved N95 or P100 filters
- Protects the patient and HCP
- Airtight seal
- Safer, you can conduct a seal check each time
- 8-point adjustable harness, no earloops, no discomfort
- White allows you to write your name or a compassionate message
- Filters and masks can be disinfected between patient interaction
- Filters can be replaced, reducing waste versus disposable style respirators
- Immediate and long-term cost savings
- Will not fog glasses
- Fits under medical style faceshields
- Compact design reduces interference with other medical equipment during patient interaction



PROTECTING THOSE WHO PROTECT OTHERS




Setting science in motion to create a better world

From breakthrough discovery to agile delivery of innovative laboratory products and services, Avantor is your trusted partner, committed to providing exceptional convenience, collaboration, and customization to enable your success.

Visit the new Avantor homepage, featuring application-specific products, new literature, and product highlights, at fishersci.com/avantor or fishersci.ca/avantor.





Pharmaceutical Applications of Artificial Intelligence

By Iva Fedorka

The pharmaceutical industry has confidence in artificial intelligence (AI). They hope that it will help to reduce workload, shorten timelines, identify drug repurposing possibilities, improve industry productivity, and make clinical trials more successful.

Drug discovery and development for a single medication can take more than a decade and cost an average of 2.8 billion USD. Despite this time and resource commitment, nine out of ten still fail Phase 2 clinical trials and never get regulatory approval.

Increased data digitalization in the pharmaceutical sector lends itself to AI tools and networks that mimic human analytical thought processes and can be programmed to interpret, “learn,” and make “intelligent” decisions.

History

AI pioneers in the 1950s talked about machines that could sense, reason, and think like people. The rapid increase in computer processing power and large amounts of data and advanced algorithm development since then have significantly improved so-called machine learning (ML). This form of AI is task-focused, supports analysis, understands and generates text and speech (natural-language processing), and truly mimics the way we think.

More than 50 years later, on June 12, 2007, the robot Adam successfully identified the function of a yeast gene. Adam searched public databases and developed hypotheses about the genes that code for key enzymes in the yeast *Saccharomyces cerevisiae*. When these hypotheses were tested, researchers found nine genes that were new and accurate and only one that was wrong.

“Robot scientists using AI can test more compounds, and do so with improved accuracy and reproducibility, and exhaustive, searchable record-keeping,” said systems biologist Steve Oliver of the University of Cambridge, a member of the group that developed Adam. Later, the same team announced that Adam’s colleague Eve had discovered a possible new treatment for drug-resistant malarial parasites.

Success Stories

Researchers at Berg, a biotechnology company near Boston, Massachusetts, used AI to find potential treatments based on the cause of a disease. They tested more than 1,000 cancerous and healthy human cell samples, varying the growth conditions and then measuring cell production and output. The key differences between the diseased and healthy cells were demonstrated by the AI analysis, which incorporated patient biological and outcome data.

“We are turning the drug-discovery paradigm upside down by using patient-driven biology and data to derive more-predictive hypotheses, rather than the traditional trial-and-error approach,” said Niven Narain, Berg’s co-founder, president, and chief executive officer. His team used this approach to identify specific cancer metabolism molecules and determine how a new cancer drug would work. The drug (BPM31510) is now in a Phase 2 clinical trial for patients with advanced pancreatic cancer. The

company is using the same AI system to find drug targets and therapies for diabetes, Parkinson’s, and other diseases.

London-based Benevolent^{AI} has a cloud-based AI platform that analyzes data from research papers, patents, patient records, and clinical trials. The database includes over one billion documented or inferred relationships between genes, symptoms, diseases, proteins, tissues, species, and medications. Used like a search engine, it can produce graphs of medical conditions, associated genes, and effective compounds. “AI can put all this data in context and surface the most salient information for drug-discovery scientists,” said Jackie Hunter, Benevolent^{AI} board director.

When Benevolent^{AI} was asked to suggest treatments for amyotrophic lateral sclerosis (ALS), they identified about a hundred existing compounds with potential. Their scientists chose five to test at the Sheffield Institute of Translational Neuroscience, UK. At the December 2017 International Symposium on ALS/MND in Boston, Massachusetts, four of the compounds were reported to show promise and one even delayed neurological symptoms in mice.

Some think the potential for AI to pinpoint previously unknown causes of disease will accelerate the trend toward treatments designed for patients with specific biological profiles. “Personalized medicine has been talked about for a long time,” said Hunter. “AI is going to enable it.”

The Future

AI has already shown its ability to predict a drug’s pharmacokinetics, target receptors, physicochemical properties, solubility, binding affinity, bioactivity, toxicity, and other features that affect its effectiveness.

Industry leaders agree that the use of AI may create changes in the drug-discovery process. Some experts opine that future students will need an understanding of biology coupled with computer science, statistics, and machine learning. Some schools have already created undergraduate programs in biomedical computation, although demand for these degrees may change as new therapies emerge.

Others believe that predictions about AI’s ability to revolutionize drug discovery are too optimistic. Narain, who expects AI to drive major advances, agrees that claims may be exaggerated, but it won’t be long before they are proved or disproved. “The hype can’t last very long because over the next five years or so, the truth will come out in the data,” he said. “If by then we are creating better drugs, and doing it faster and cheaper, then AI will really take off.”

This content was inspired, in part, by “How artificial intelligence is changing drug discovery,” Nature, May 30, 2018; “Transforming Drug Discovery Through Artificial Intelligence,” Forbes, March 3, 2020; and “Artificial intelligence in drug discovery and development,” Drug Discovery Today, January 26, 2021.

Iva Fedorka is a Thermo Fisher Scientific staff writer.

Fact or Friction: Your Balance Static Questions Answered

Is your balance display slowly drifting while trying to take a weight? Are balance readings unstable?

These issues can be attributed to a well-known culprit: static electricity. Electrostatically charged samples can not only create frustration but also adversely affect weighing accuracy and waste material and lead to significant errors. Here are some common questions and answers about static.

Fact: Friction is by far the most common cause of electrostatic discharges.

Many laboratory tasks — like using a cloth to dry a glass beaker, picking up a measuring flask with disposable gloves, or filling a weigh pan with powder — can generate measurable electrostatic charges that interfere with weighing results.

Fact: The climate conditions of your lab affect your weighing results.

Relative humidity levels less than 40 to 50 percent increase electrostatic charge issues. The lower the humidity, the less chance that charged particles will find water molecules to “surf” on to get to a convenient place.

Depending on the relative humidity, a charge may take anywhere from a few seconds to several minutes to dissipate. In controlled and dry atmospheres with less than 20% relative humidity, charges on materials can cause weights to drift several hundred milligrams and may persist for many hours.

Fiction: There is no way to dissipate electrostatic charges.

METTLER TOLEDO anti-static kits use alternating current (AC) to generate ions

and immediately neutralize electrostatic charges from samples or containers. This method is effective, completely safe, and does not disturb air currents or increase stabilization time. The time savings, efficiency gains, and safety improvements make it a sound investment for your lab.

The XPR Solution

Advances in METTLER TOLEDO analytical weighing technology enable the latest generation of analytical balances to automatically detect, measure, and record electrostatic charges during a weighing operation. The integrated ionizer module can then eliminate these charges to avoid any influence on the weighing result. The static detection cycle is performed in just a few seconds while the balance settles, so there's no delay in obtaining the weighing result.



Description	Capacity	Readability	Mfr. No.	Cat. No
Balances				
XPR Balance XPR205, Single Range	220 g	0.01 mg	30355411	01-919-124
XPR Balance XPR225DR, Delta Range	121/220 g	0.01/0.1 mg	30594483	01-915-004
XPR Balance XPR225DU, Dual Range	121/220 g	0.01/0.1 mg	30594488	01-915-003
Ionizing Modules and Kits				
XPR Integrated Internal Ionizer Modules			30460823	13-940-060
Anti-Static Kits, with Large U-Electrode and Power Supply*			63052302	01-910-021

*Universal kit; works with all balances and weighing substances

Content provided by:

METTLER TOLEDO



Visit fishersci.com/mettlertoledo or fishersci.ca/mettlertoledo to learn more and invite excellence into your laboratory today.

HOMOGENIZING PERFECTED — for a healthier world

- ENCAPSULATION
- API WET-MILLING
- MICROSPHERES



Complete formulation unit for R&D process technology

APPLICATIONS INCLUDE:

- | | |
|---------------------------------|---|
| Eye drops | Ultrasound gel (contact gel) |
| Eye injection | Encapsulations |
| Wet-milling of API | Lyophilization |
| Extraction of API | Probiotics in pharma |
| Pharmaceutical aerosol | Analytic in quality control |
| Sterile liquid oral suspensions | Ointments & creams |
| Vaccines | Crystallizing of pharmaceutical ingredients |
| Medications | Liposomes |
| Ophthalmic suspensions | CMC solutions |
| Saline solutions | Inclusion bodies |
| Foamed patches | Fat particles |
| Phytopharmaceuticals | PLGA capsules |
| Micro-encapsulations | Cough syrup |
| Sterile suspensions | Clear dispersion |
| Infusion solutions | Suppository |
| Tablet coatings | |
| Foamed drug blisters | |



SWISS MADE DISPERSING SOLUTIONS FOR YOU. SINCE 1962.

Extend Disposable Glove Life

WELLS LAMONT
Industrial

Wells Lamont Glove Liners

Wells Lamont glove liners can help extend the life of disposable gloves, reduce usage, help with sizing issues, and reduce potential allergic reactions to glove materials.

- **Models M005 and M006:** Reusable (wash up to 30x), low lint, 100% continuous nylon filament, half- and full-finger options, color coded sizes, made in the USA
- **Model Y6701:** Disposable, 100% lightweight white cotton, full fingers, women's and men's sizes
- **Models M088 and M089:** Reusable, 100% white nylon knit, form fitting, flexible, and comfortable, half- and full-finger options, universal sizing, made in the USA




Model	Type	Material	Finger Style	Color	Size	Quantity	Cat. No.
M005	Reusable	100% Nylon Filament	Full	Green	Small	12/Pack	19-086-521
				Orange	Medium	12/Pack	19-086-520
				Blue	Large	12/Pack	19-086-522
				Red	X-Large	12/Pack	19-002-852
M006	Reusable	100% Nylon Filament	Half	Red	Women's	12/Pack	19-086-524
				Black	Men's	12/Pack	19-086-525
Y6701	Disposable	100% Cotton	Full	White	Women's	12/Pack	11-462-26A
				White	Men's	12/Pack	11-462-26B
M088	Reusable	100% Nylon Knit	Half	White	Universal	12/Pack	18-999-4591A
M089	Reusable	100% Nylon Knit	Half	White		12/Pack	18-999-4591B

SLi-6000

Fluorescence Microscope


A Giant Leap Forward in
Fluorescence Microscopy

SeBaLit Epi-Fluorescence Illumination (patent pending) developed by Laxco brings a simplified approach to fluorescence microscopy



Fluorescence

Enhanced Darkfield Contrast provides a perfect solution for locating the focal plane or area of interest when imaging in fluorescence



Enhanced Darkfield
Contrast



Maximizing Pipette Accuracy

One of the most frequently used tools in a typical laboratory is the air-displacement or air-cushion pipette. Convenient and easy to use, single-use pipette tips minimize the chance of cross contamination between samples.

Do you know how to select the correctly sized pipette for maximum accuracy? Let's first look at the mechanics.

Air displacement pipettes are based on a simple concept:

- A piston of consistent diameter is pushed through a seal to move a column of air; the air movement, in turn, moves liquids in and out of a disposable pipette tip
- An air cushion exists inside the pipette tip between the liquid and the end of the pipette; the liquid never touches the pipette itself
- After a pipette tip is affixed, the piston is moved from its resting position and through the piston seal to push air out of the pipette tip
- The tip is then immersed just below the surface of the liquid to be aspirated; the spring-action piston is returned to its resting position, creating a vacuum that fills the pipette tip with liquid

- Adjustable volume pipettes use a screw mechanism to change the stroke length to set the volume at 50% of the nominal capacity

For pipettes with nominal volumes of 100 μL or greater, the inaccuracy of an air displacement pipette is treated as a constant volume throughout the range of the pipette. ISO 8655 Part 2 lists this value at $\leq \pm 0.8\%$ of the nominal capacity (i.e. $\leq \pm 0.8 \mu\text{L}$ for a 100 μL pipette) and depends mostly on the mechanics of the pipette. For smaller volumes, the inaccuracy increases due to the compressibility of air and other factors.

Typically, pipette data sheets or catalogs present the accuracy as a percentage of the nominal volume, often with no additional information, so the user must determine the implications. When estimating the accuracy percentage for partial volumes, the numerator (the amount of inaccuracy) remains the same; the denominator (the set volume) changes.

Clearly, an air displacement pipette is most accurate at or close to its nominal capacity. Furthermore, accuracy decreases significantly as the set volume goes below 50% of the nominal capacity. Most pipettes therefore have a stated volume range with both upper and lower volume limits.

The amount of liquid to be pipetted often falls within the range of several different pipettes. For example, 20 μL could be measured with 2 to 20 μL , 5 to 50 μL , or 10 to 100 μL pipettes. The volume would be most accurately pipetted using the 2 to 20 μL pipette. If that size is not available, choose the pipette with the lowest capacity relative to the volume you need. In this case, the 5 to 50 μL pipette is a better choice than the one that measures from 10 to 100 μL .

If your lab has pipettes with only marginally overlapping ranges, the addition of one or two more pipettes with more overlap may increase the accuracy of your measurements.

The BRAND Transferpette S is a manual single-channel pipette available in ten different adjustable volume sizes (0.1 to 1 μL up to 1 to 10 mL); these sizes cover nearly the entire range from 0.5 μL to

Volume Setting	Percentage of Nominal Capacity	Volume of Inaccuracy	Pipetted Volume Accuracy
100 μL	100%	$\leq \pm 0.8 \mu\text{L}$	$\leq \pm 0.8\%$
75 μL	75%	$\leq \pm 0.8 \mu\text{L}$	$\leq \pm 1.1\%$
50 μL	50%	$\leq \pm 0.8 \mu\text{L}$	$\leq \pm 1.6\%$
20 μL	20%	$\leq \pm 0.8 \mu\text{L}$	$\leq \pm 4\%$
10 μL	10%	$\leq \pm 0.8 \mu\text{L}$	$\leq \pm 8\%$



Content provided by:



10 mL, overlapping the upper halves of the pipette ranges. Similarly, manual Transferpette 8- and 12-channel pipettes are available in five adjustable volume ranges (0.5 to 10 μ L to 30 to 300 μ L) for 1 to 300 μ L sample volumes.

The BRAND Transferpette is also available in an electronic version with similar options. Five different single-channel Transferpette electronic pipettes cover volume ranges from 0.5 to 10 μ L to 250 μ L to 5 mL; five multichannel models cover volumes from 0.5 to 10 μ L to 15 to 300 μ L.



NOTE: In most literature, "accuracy" is used to denote the amount of inaccuracy. To stay consistent with other literature, we use "accuracy" in this article when referring to percentages.



Boekel Scientific Programmable Heating Blocks

- Large color touchscreen with intuitive navigation
- Easy to use at single-temperature setpoints
- Create multistep temperature control programs

Choose between 1-, 2-, and 4- block models with many block options.

Visit fishersci.com or contact your Fisher Scientific sales representative to learn more.

Cat. Nos. 13-245-319, 13-245-321, 13-245-323



SMARTER SUSTAINABILITY

BAGS vs BOTTLES

- Less Storage Space
- Easy, Efficient Shipping
- User Friendly
- Less fossil fuel, emissions, and water used in manufacturing and transportation
- **Cost Effective**



Ditch your rigid plastic sample collection containers today!

Run Your Assays with Confidence

The BioTek Synergy Neo2 Multi-Mode Microplate Reader is designed to meet the needs of high-throughput laboratories, automating both cell-based and biochemical assays.

BioTek Synergy Neo2

The BioTek Synergy Neo2 Multi-Mode Microplate Reader features BioTek's patented Hybrid Technology with independent optical paths that help ensure uncompromised performance in all detection modes.

- CO₂/O₂ controls, incubation to 65°C, variable shaking for live cell assays, direct bottom illumination, and other advanced environmental features
- Powerful Gen5 Software for complete reader control, powerful data analysis, automation, and LIMS integration
- TRF and Alpha lasers for speed and improved sensitivity
- Patented Hybrid Technology with independent filter and monochromator-based optics
- Variable bandwidth selection to optimize fluorophore specificity and sensitivity
- Multiple PMT detectors for fast plate processing speeds



Description	Cat. No.
<ul style="list-style-type: none"> • Single Top PMT, Top Filter Fluorescence, Fluorescence Polarization, Time-Resolved Fluorescence, TR-FRET, Luminescence, Filter Luminescence, Monochromator UV-Vis Absorbance • For 1- to 1536-Well Microplates • Incubation to 65°C, Shaking 	BTNEO2S
<ul style="list-style-type: none"> • Dual Top PMT, Top Filter Fluorescence, Fluorescence Polarization, Time-Resolved Fluorescence, TR-FRET, Luminescence, Filter Luminescence, Monochromator UV-Vis Absorbance • For 1- to 1536-Well Microplates • Incubation to 65°C, Shaking 	BTNEO2
<ul style="list-style-type: none"> • Dual Top PMT, Top and Bottom Monochromator Fluorescence with Variable Bandwidth, Top and Bottom Filter Fluorescence, Fluorescence Polarization, Time-Resolved Fluorescence, TR-FRET, TRF Laser, Luminescence, Filter Luminescence, Monochromator UV-Vis Absorbance, Alpha Laser • For 1- to 1536-Well Microplates • Incubation to 65°C, Shaking 	NE2MALPHABT

Save 40% on a Microplate Stacker

When you purchase a BioTek Synergy Neo2 Multi-Mode Microplate Reader, you'll get 40% off a BioTek BioStack Neo Microplate Stacker.

Visit fishersci.com/synergy-neo2-lab-shot or fishersci.ca/synergy-neo2-lab-shot to learn more.

Restrictions apply. See page 48 for details.

Other configurations are available. Contact your Fisher Scientific sales representative for details.

Next-Generation Vacuum Pumps



KNF Neuberger LABOPORT Vacuum Pumps

New oil-free and chemically resistant LABOPORT diaphragm vacuum pumps improve everyday laboratory practices with integrated speed control, three-color status display, and exceptionally small footprints.

Compact LABOPORT vacuum pumps are easy to use, durable, and useful for rotary evaporation, degassing, filtration, SPE, fluid aspiration, gel drying, centrifugal concentration, vacuum ovens, and many other applications.

Visit fishersci.com/laboportvacuumpumps to shop.



Model	Application	Ultimate Vacuum	Flow Rate	ATEX Compliant	Integral Gas Ballast Valves	Cat. No.
N 96	Filtration, SPE, Aspiration	97.5 torr, 130 mbar	7 L/min.	No	No	13-880-904
N 820 G	Rotary Evaporation, Degassing, Fluid Aspiration, Centrifugal Concentration, Vacuum Oven, Gel Drying	4.5 torr, 6 mbar	20 L/min.	Yes	Yes	13-880-905
N 840 G	Rotary Evaporation, Filtration, Centrifugal Concentration, Vacuum Oven	4.5 torr, 6 mbar	34 L/min.	Yes	Yes	13-880-906




STAY COMPLIANT WITH EPA AND SPCC REGULATIONS

Youngstown Barrel & Drum Co. has a containment solution to fit your need. Hold up to four 55 gallon drums with Ultra-Spill Pallets, or get low-profile containment with Ultra-Spill Decks.



Description	Cat. No.
Ultra-Bung Access Funnels	19-130-4814
Ultra-Spill Deck Bladder Systems, Capacity: 77 gal.	19-058-005
Ultra-Drain Seal Mats, 36 x 36 in.	19-823-211
Ramp, Nonslip, Polyethylene, 55.5 x 28.5 x 8.75 in.	19-130-4953
Ultra-Spill Pallets Plus, 62 x 62 x 8.75 in.	19-130-4780
Ultra-Drum Rack P2 Racks, Holds Two (55 gal.) Drums	19-820-937
Ramp, 24 x 32 x 5.75 in.	19-058-033
Ultra-Spill Deck, Four Drum, 52 x 52 x 5.75 in.	19-058-074
Ultra-Spill Berm, Orange, 10 ft. x 4 in.	19-058-009
Ultra-Spill Trays, Capacity: 1.8 gal. (6.8 L), 17 x 12 x 4 in.	19-130-2734
Ultra-Spill Trays, Capacity: 2.9 gal. (11 L), 21 x 17 x 4 in.	19-130-2735



fisherbrand

FOR THE **LOVE** of **SCIENCE**

The Fisherbrand collection offers value-packed products that make your lab life easier — every day, when it matters most — guaranteed.

Visit fishersci.com/fisherbrand or fishersci.ca/fisherbrand to shop.

ELEVATE *your* EVAPORATOR

For years Heidolph has been redefining rotary evaporation. Collaborating with researchers, Heidolph designed a new generation of rotary evaporators.

The Hei-VAP Series includes superior ease of use, premium quality materials and uniquely designed features to overcome any chemical or pharmaceutical challenge.

 **heidolph**
research made easy




The Right Tools for Multicolor Flow Cytometry

Although the principles of good panel design match antigen density with fluorochrome brightness, about 20 percent of panels lack at least one required reagent. Improve resolution and generate high-quality data for better scientific insights with the BD reagents.

Visit fishersci.com/bdbiosciences to shop.



Mouse CD162 Comparison



Distributed by Fisher Scientific

Promotion Terms and Conditions

These promotions are not available to our current or potential healthcare customers. Healthcare customers include any healthcare practitioner, other healthcare provider or any individual or organization authorized to prescribe, dispense, purchase or influence the acquisition or use of medical devices or supplies for clinical use. Healthcare customers include, without limitation, any "Health Care Practitioner" under the Massachusetts Marketing Code of Conduct Law (defined at 105 C.M.R §907.004) or a "Health Care Provider" under the Vermont Prescribed Products Gift Ban and Disclosure Law (defined at Vt. Stat. Ann. tit. 18, §4631a(8)). By participating in any of these promotions, you warrant that you are not a current or potential healthcare customer. Thermo Fisher Scientific reserves the right to determine your status in the event of a dispute.

Customer acknowledges that these offers may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 U.S.C. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

The factual representations regarding the products accurately reflect the manufacturers' representations to Fisher Scientific. To the extent permissible, Fisher Scientific will make available to its customers the manufacturers' warranties. The sole and exclusive liability of Fisher Scientific shall be refund of the purchase price or replacement of the product. In no event shall Fisher Scientific be liable for incidental or consequential damages. **FISHER SCIENTIFIC EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS, IMPLIED, OR STATUTORY, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.** Prices may change without notice. Prices are U.S. list prices — pricing may vary in Canada. Content and printed prices subject to confirmation by your Fisher Scientific Sales Representative. Telephone calls may be monitored as part of Fisher Scientific's quality process.

© 2021 Thermo Fisher Scientific Inc.
All rights reserved.
Trademarks used are owned as indicated
at fishersci.com/trademarks.
Litho in U.S.A. 21-507-0049 KP 2/21
BN20211207