

## SOP for Periodic Repeatability Tests (Routine Tests)

<b>Title:</b>	SOP for Periodic Repeatability Tests
<b>Document No.:</b>	
<b>Author:</b>	METTLER TOLEDO Laboratory & Weighing Technologies
<b>Number of Pages:</b>	4
<b>Controls</b>	
<b>Replacement Document:</b>	N/A
<b>Reason for Revision:</b>	N/A
<b>Release Date:</b>	17. Oct. 2007
<b>Release for Routine Operation</b>	
<b>Reviewed by:</b>	
<b>Date:</b>	
<b>Signature:</b>	
<b>Approved by:</b>	
<b>Date:</b>	
<b>Signature:</b>	

## General

### Basic Rules for Handling Balances

- Before using a balance, make sure the balance was left on power for a sufficient period of time (mentioned in the balance operating instructions).
- Make sure the balance is leveled.
- Minimize environmental influences, e.g. open windows, direct sunlight or strong drafts.
- Do not enter the draft shield with hands. Use gloves or long tweezers.
- Place objects gently and in the center of the pan.

### Basic Rules for Handling Weights

#### Test Weights

- Only an external test weight with calibration certificate can make a balance a “traceable” piece of equipment.
- Test weights should always be placed gently on a clean weighing pan and put back immediately in their storage place after use.
- Test weights (since they are also part of measuring equipment) need to be re-calibrated at specified intervals (ISO 9001).
- Any incident, which might have affected the value of the test weight, should trigger an immediate re-calibration. METTLER TOLEDO’s calibration services will give advice on this.

#### How to Store Test Weights

- Test weights should be stored in their original box.
- Test weights should be stored in the same room as the balance they are used with, since temperature differences between test weights and their surrounding lead to measurement errors.
- Test weights that have not been stored at the same temperature need acclimatization, which can take several hours.

#### How to Move Weights

- Test weights should only be handled with appropriate tools such as tweezers, forks, handles or gloves (see METTLER TOLEDO’s accessories for weights).
- These tools should be exclusively used for transferring test weights, to avoid possible contamination.

# Repeatability Test

## Preparation

- Before performing the test, the test weight must be acclimatized to the ambient temperature of the balance.
- Prepare a sheet where you can note the readings of the repeatability test points. Prepare 3 columns titled "w/o Test Wght", "With Test Wght." and "Differences".

## Test Procedure

- Empty the pan
- If required, place tare load on the weighing pan
- Tare the balance (if required press zero)
- Read the stable value from the display and note it in the column "w/o Test Wght."
- Place the test weight in the center of the weighing pan
- Read the stable value from the display and note it in the column "With Test Wght."
- Remove the test weight
- Repeat the measurements from "tare the balance" to "remove test weight" until the customer defined number of "w/o Test Wght." and "With Test Wght." readings is performed

## Evaluation

- Subtract each "w/o Test Wght." from the corresponding "With Test Wght." to calculate the difference.
- Calculate the standard deviation of the differences.
- Evaluate whether standard deviation exceeds the defined "Control Limit" <sup>1)</sup>.

## Deviation

### Control Limit <sup>1)</sup>

- If the control limit is exceeded, report the problem to the laboratory supervisor or the person responsible of the balance.
- Mark the balance as "out of control limits".
- Contact METTLER TOLEDO service organization for advice.

<sup>1)</sup> – Values within the control limit: No action is necessary.

– Values beyond the control limit show that weighing process is no longer under control and immediate action is therefore required.

**Mettler-Toledo AG**

Laboratory & Weighing Technologies  
Im Langacher  
P.O. Box LabTec  
CH-8606 Greifensee  
Switzerland

Subject to technical changes  
© 06/2009 Mettler-Toledo AG  
Printed in Switzerland 11793057

[www.mt.com/GWP](http://www.mt.com/GWP)

---

For more information