

# LAB

R E P O R T E R

Science Innovations and Discoveries

NO. 1, 2023

# Ancestral Antiviral Proteins

## May Protect Against Infection

### Also Featured:

A New Commitment to **Renewable Energy**

How a Safety Specialist Can Help **Improve Safety in Your Lab**

Successful **Cell Therapy Manufacturing**


# CONTENTS



## 18

### COVER STORY

Ancestral Antiviral Proteins May Protect Against Infection



## 4

A New Commitment to Renewable Energy




## 10

How a Safety Specialist Can Help Improve Safety in Your Lab



## 26

First Plants Grown in Lunar Soil May Lead to Farming in Space



## 30

Successful Cell Therapy Manufacturing

**Note:** To check pricing, please sign in to your account.

**Supplier and Other Trademarks:** Agilent BioTek, Avantor, BAKERBOND, Brady, CELLMASTER, CiDehol, Cytation, Decon, DuPont, Erlab, FlaskScrubber, Gen5, Greiner Bio-One, J.T.Baker, KNF, Labconco, LabPartner, Latitude, Macron Fine Chemicals, METTLER TOLEDO, Mystaire, PROchieva, Saf-T-Zone, SevenDirect, SIMDOS, SmartCheck, SteamScrubber, Tychem, Whirl-Pak

**Trademarks of Thermo Fisher Scientific:** Fisher Scientific, Fisherbrand, Lab Reporter, MaximaDry, Thermo Scientific

For trademark ownership, see [fishersci.com/trademarks](https://fishersci.com/trademarks).

Visit [fishersci.com](https://fishersci.com) or [fishersci.ca](https://fishersci.ca) to shop.

# LAB REPORTER

**Lab Reporter** provides quick and easy access to today's cutting-edge products and trusted solutions for all of your scientific research and applications.

Science Innovations and Discoveries

## SUPPLIER ARTICLES



16

Cleanroom-Compatible  
CO<sub>2</sub> Incubators for  
Cell Therapy

Thermo Scientific



24

Sustainability in the  
Laboratory


Whirl-Pak



35

The Cost of Air and  
the Impact on  
the Laboratory

Erlab



32

Confocal Imaging  
and Analysis

Agilent BioTek

## SUPPLIER PRODUCT GUIDE


<b>Agilent BioTek</b> Cytation C10 Imaging Reader .....	15
<b>Avantor</b> J.T.Baker and Macron Fine Chemicals .....	7
<b>Brady</b> SmartCheck Safety Management System .....	6
<b>Decon</b> CiDehol 70 Isopropyl Alcohol .....	14
<b>DuPont</b> Tychem 2000 Garments and Accessories .....	29
<b>Fisherbrand</b> Bead Mills .....	17
<b>Fisherbrand</b> MaximaDry Pumps .....	34
<b>Fisherbrand</b> Multi-Pad pH Indicator Strips .....	22
<b>Fisher Scientific</b> Channel .....	22
<b>Fisher Scientific</b> Chromatography .....	23
<b>Fisher Scientific</b> Contamination Control .....	13


<b>Fisher Scientific</b> Edge Same-Day Shipping .....	17
<b>Fisher Scientific</b> Semiconductor Products .....	9
<b>Fisher Scientific</b> Spotlight Savings .....	34
<b>Fisher Scientific</b> Supplier Diversity .....	28
<b>Fisher Scientific</b> Sustainability .....	27
<b>Greiner Bio-One</b> CELLMASTER Cell Culture Roller Bottles .....	36
<b>KNF</b> SIMDOS Pumps .....	27
<b>Labconco</b> Glassware Washers .....	28
<b>METTLER TOLEDO</b> Instruments .....	14
<b>Mystaire</b> LabPartner, Latitude, and MY-PCR Enclosures .....	12



Visit [fishersci.com/labreporter-subscribe](http://fishersci.com/labreporter-subscribe)  
or [fishersci.ca/labreporter-subscribe](http://fishersci.ca/labreporter-subscribe) to  
subscribe for free.

 **fisher scientific**  
part of Thermo Fisher Scientific

 This publication was printed on recycled paper.



# A New Commitment to Renewable Energy

By Kylie Wolfe

States across the Midwest are leading the quest for renewable energy in America. After Texas and Iowa, Oklahoma is the third largest producer of wind energy in the United States. Its open landscapes are home to over a dozen wind farms, one of which is helping Thermo Fisher Scientific reach critical sustainability milestones.

In September 2022, Thermo Fisher furthered its commitment to renewable energy, purchasing a 90-megawatt portion of the Seven Cowboy wind project in western Oklahoma. This step supports the company's long-term plan to reduce emissions.

## Harnessing Wind in Western Oklahoma

Managed by Enel Green Power, the Seven Cowboy wind project is comprised of 107 wind turbines that are expected to generate roughly 1.3 terawatt hours of energy each year. This amount of power could keep over 120,000 U.S. households running and at least 758,000 tons of CO<sub>2</sub> emissions out of the atmosphere.

Through an eight-year virtual power purchasing agreement (VPPA), Thermo Fisher secured 90 megawatts of the overall system. This equates to 400,000 megawatt hours of renewable electricity each year, half of the company's electricity usage in the U.S.

"If we were to take every site and cover every roof, parking lot, and adjacent property with solar panels or wind turbines, we still wouldn't be able to match the total combined electricity used across the organization. These VPPAs, or off-site projects, give us that flexibility and the ability to fundamentally address our Scope 2 emissions at the scale that we're actually producing them," said Scott Self, director of procurement at Thermo Fisher Scientific.

## An Important Step Forward

The U.S. Environmental Protection Agency (EPA) has different classifications for emissions. Scope 1 emissions are defined as "direct greenhouse (GHG) emissions that occur from sources that are controlled or owned by an organization." Scope 2 emissions are those that result from the purchase of utilities like electricity. These are operational emissions associated with buildings and on-site manufacturing efforts. Scope 3 emissions are indirect and associated with an organization's supply chain, often the largest sources of emissions.

In 2021, Thermo Fisher joined the Business Ambition for 1.5°C campaign—led by the Science Based Targets initiative—and committed to net-zero emissions by 2050. This included reducing Scope 1 and 2 emissions by 30 percent compared to 2018 values. This goal was recently increased to 50 percent.

The Seven Cowboy wind project helps address Scope 2 emissions by generating renewable energy credits. It's those credits that are used to match emissions from sites across the country. By improving Scope 2 emissions, Thermo Fisher is not only supporting clean energy but also helping customers meet their Scope 3 reduction targets. That's because a company's Scope 1 and 2 emissions contribute to their customer's Scope 3 emissions.

"These VPPAs, or off-site projects, give us that flexibility and the ability fundamentally to address our Scope 2 emissions at the scale that we're actually producing them."

### Scott Self

*Director of Procurement  
Thermo Fisher Scientific*

"When you think about the relationships between us and our customers and even us and our suppliers, it's like links in a chain that all support each other," said Self.

The project is expected to be operational in the second half of 2023.

## Plans for a Sustainable Future

Today, Thermo Fisher has more than 60 sites around the world that exclusively use renewable electricity. All sites in Germany and Italy fall into this category, relying on sources like the Seven Cowboy project. Some sites also have their own solar panels and wind turbines, generating electricity physically on site. In Cork, Ireland, a wind turbine helps match a lot of the emissions at the facility.

Self says that the company is looking to expand on-site solar initiatives wherever possible. Over 20 projects are currently in development, including two physical installations at Massachusetts-based sites. His team is also looking to remove equipment that uses fossil fuels. By taking steps to add renewable systems on site and leverage long-term power purchasing agreements, Thermo Fisher is making strides toward a more sustainable future.

*Kylie Wolfe is a Thermo Fisher Scientific staff writer.*

# See the whole safety picture. Finally.

Safety software for the lab facilities manager

All your information. All in one place. Automate tedious tasks, prevent costly errors, and protect lab faculty and staff with the fully configurable SmartCheck™ Safety Software & Management System.



## Connect processes and people

Organize under a single online system to keep everyone you need in the know and on the go.



## Uncover and solve safety issues

Digitize and automate inspection and auditing to better identify, resolve, and recall hazards and incidents.



## Meet policies and regulations

Log and track all chemical SDSs and auto-populate compliant labels for printing, including GHS.



## Be a pro from day 1

A new system doesn't mean new headaches ... not with SmartCheck's included training and support.



## Own your safety program.

▶ Scan the code or call your Fisher Scientific representative to set up a free, no-obligation demo.





# Setting science in motion to create a better world



From breakthrough discovery to agile delivery of innovative laboratory products and services, Avantor® is your trusted partner, committed to providing exceptional convenience, collaboration and customization to enable your success.

Visit the new Avantor homepage, featuring application specific products, new literature and product highlights at [fishersci.com/avantor](https://fishersci.com/avantor) or [fishersci.ca/avantor](https://fishersci.ca/avantor).



**J.T. Baker® BAKERBOND® PROchieva™**

PROchieva™ will change your mind about protein performance.

Request a sample



*fisherbrand*

# Quick-Ship or Custom-Built

Choose from the Fisherbrand quick-ship furniture portfolio or work with us on a custom design.

The Fisherbrand line of lab furnishings provides options to meet the specific requirements of your lab environment, timetable, and budget.

Our specialists can work with you to design, select, and install furniture and equipment that suits your exact specifications. You can also order from our quick-ship portfolio to have items shipped within 14 days.\* Either way, you've found a fast, convenient source for getting the lab furnishings you need, when you need them.

Visit [fishersci.com/fisherbrand-furniture](https://fishersci.com/fisherbrand-furniture) or [fishersci.ca/fisherbrand-furniture](https://fishersci.ca/fisherbrand-furniture) to shop and learn more.

## Build Your Own Chair

Create chairs that are exactly right for your lab and your work. Just access our online tool and select the type, color, materials, base, height, and more. The tool will configure your chair and help you place your order.

Visit [fisherchair.com](https://fisherchair.com) to build your chair now.

*Available in U.S. only.*

\*Inventory of select items is subject to availability at the time of order. Large quantity orders may require longer lead times.

 **fisher scientific**  
part of Thermo Fisher Scientific



# Power the Future

## Products and Services to Fuel the Semiconductor Resurgence

Whether you're in R&D, fabrication and production, or quality control, inspection, and testing, the Fisher Scientific channel offers a comprehensive selection of supplies, equipment, and capabilities to drive your semiconductor work forward.

### Products to enable your routines:

- Apparel, wipers, respirators, and more to mitigate dust and particle contamination and support controlled and cleanroom environments
- Hand, foot, ear, and eye protection to help keep your people safe
- Materials storage, signs and tags, and electrical, industrial hygiene, and environmental monitoring solutions to promote safety for facilities, people, and products
- High-purity, electronic-grade solvents, etchants, coolants, acids, bases, analytical standards, and metals and salts for analysis and production
- Recirculating chillers, environmental chambers, ovens, electron microscopes, meters, as well as glass, plastic, and metalware to help meet exacting design, research, production, and inspection standards
- Furniture and storage for an array of safety and functional needs
- Antistatic bags and tapes to minimize electrostatic discharge

### Services and programs that scale with your operations:

- Distribution and warehousing capabilities
- Procurement and inventory services
- Third-party product sourcing
- Chemical stock and selection tools
- Safety services and expertise
- In-stock, same-day shipping programs
- Sustainable product and recycling resources



Visit [fishersci.com/semiconductor](https://fishersci.com/semiconductor) or [fishersci.ca/semiconductor](https://fishersci.ca/semiconductor) to learn more about offerings for semiconductor applications from the Fisher Scientific channel.



# How a Safety Specialist Can Help Improve Safety in Your Lab

By Christina P. Hooton

Underlying hazards exist in every facility, and the lab is no different. The key to maintaining a safe work environment for your employees, faculty, and students is identifying hazards and putting in place the personal protective equipment (PPE), training, and resources required to mitigate those risks.

Whether you're just starting your safety program or need to reinvigorate your existing safety processes, a safety specialist can serve as your trusted resource and offer the support you need to prioritize safety in your lab.

## Identifying Hazards

When designing a safety program, the lab manager or the environmental, health, and safety (EH&S) professional supporting the facility will perform a hazard assessment. During the assessment, they will identify all the hazards within a lab, prioritize each, and define the proper engineering controls, work practices, PPE, or other safety supplies to use in those situations.

"Safety specialists with a distributor or a PPE manufacturer can help support a hazard assessment, providing core information and then personalizing it to their particular organization and protocol," said Don Doyle, Fisher Scientific senior safety specialist, Thermo Fisher Scientific.

For example, hazardous waste may be present in the lab. In this situation, Charlie Fink, Fisher Scientific senior safety specialist, Thermo Fisher Scientific, recommends a site assessment. "It involves a walkthrough of your lab by a safety specialist who compiles a confidential report indicating where you may or may not be complying with applicable standards. I provide assessments to my customers who want to make sure they are properly storing hazardous materials and wastes and have the right products and supplies on hand to safely transport and dispose of wastes," said Fink.

## Selecting PPE

After PPE requirements are identified through the hazard assessment, a safety specialist can make product suggestions for your lab based on your specific needs. For example, hand protection is needed in every lab, but selecting the right option depends on the work being done, regulatory requirements, and your business performance goals.

"Most customers are dealing with sharps, which could be blades, broken glass, etc., so they need a good solution for cut protection. In these situations, I usually recommend an 18-gauge cut liner that can be worn underneath their exam gloves," said Doyle. "Cryo gloves are another common solution. Although they are generally large, puffy gloves, there are options that offer high dexterity, are suitable for -80°C cryo work, and are great for handling frozen cryo tubes."

## Optimizing Education and Training

Your lab personnel need to be aware of and comply with the standards and regulations of various government agencies, including the Occupational Safety and Health Administration (OSHA), the U.S. Department of Transportation (DOT), the U.S. Environmental Protection Agency (EPA), and the U.S. Food and Drug Administration (FDA), depending on their work. It's also important for them to be aware of your organization's internal standard operating procedures (SOP) and the requirements of state and local governments.

The initial assessment will help reveal which standards and regulations lab personnel need to adhere to and identify training and education opportunities. Safety specialists can make additional suggestions and connect labs with safety training and services as well as training in specific areas, such as contamination control.

## Focusing on Sustainability

In addition to maintaining safety standards, many labs need to meet sustainability goals that are intertwined with safety practices and products. Recycling and disposal can be challenging because of the many different waste streams—from hazardous and regulated medical waste to more traditional streams such as glass and plastic.

"From a lab and safety perspective, implementing sustainability measures could mean following green chemistry practices; recycling gloves, apparel, masks, and pipette tips; using biodegradable gloves and other products; and just properly disposing of and reducing waste of everyday safety supplies," said Josh Boyle, Fisher Scientific senior safety specialist, Thermo Fisher Scientific.

Relying on a safety specialist can make this process less daunting. They can provide greener product recommendations as well as suggest recycling programs that simplify the complexities.

No matter the hazard or area of focus, safety specialists can provide the extra support you need to improve safety in your lab every step of the way.

Visit [fishersci.com/contact-safety-specialist](https://fishersci.com/contact-safety-specialist) or [fishersci.ca/contact-safety-specialist](https://fishersci.ca/contact-safety-specialist) to find out who your Fisher Scientific safety specialist is or if you need assistance with any of your safety or PPE needs.

Visit [fishersci.com/safety-services-team](https://fishersci.com/safety-services-team) or [fishersci.ca/safety-services-team](https://fishersci.ca/safety-services-team) to find out more about the services we offer.

*Christina P. Hooton is a Thermo Fisher Scientific staff writer.*

Mystaire® manufactures fume hoods and laminar flow hoods to provide operator or process protection from toxic fumes, vapors and particulates.



## MY-PCR Workstations

MY-PCR workstations establish an ISO 5 clean work area with timed UV light. MY-PCR workstations create a “personal” clean work zone.

Standard Features:

- Combines vertical ISO 5 air with UV light sterilization
- Safety switch to eliminate UV light exposure
- Polycarbonate construction with polypropylene base

Description	Mfr. No.	Cat. No.
MY-PCR, W:24"	MY-PCR24	15-338-365
MY-PCR, W:32"	MY-PCR32	15-338-366
MY-PCR, W:48"	MY-PCR48	15-338-367

MY-PCR filters included in the purchase price

## LabPartner Workstation

LabPartner Workstations are compact, ductless solutions for capturing fumes and particulates during routine laboratory analysis. This lightweight model is made from materials that are durable in typical laboratory environments.

Standard Features:

- User-controlled variable air volume, 30-120 fpm
- 360° clear construction for optimal visibility
- Vapor proof fluorescent light standard on 24 and 48 inch models

Description	Mfr. No.	Cat. No.
LabPartner, W: 24"	MY-2400	15-338-184
LabPartner, W: 32"	MY-3200	15-338-185
LabPartner, W: 48"	MY-4800	15-338-188



## Latitude™ Series C Filtered Hoods

Latitude™ Series C filtered hoods offer superior performance for particulate and vapor containment while maximizing user accessibility and incorporating redundant HEPA filters for added safety and easy maintenance.

Standard Features:

- Saf-T-Zone™ technology allows safe HEPA filter change while hood is operating
- Exhaust air is HEPA-filtered to keep potential particulates from returning to the room
- Filter safety lockout to ensure correct filter replacement
- Deep blue base effectively contains accidental spills

Description	Mfr. No.	Cat. No.
Latitude, W: 48"	MY-LBE48	15-338-970
Latitude, W: 72"	MY-LBE72	15-338-971



# Minimize Risk for Your Products, Facilities, and People

Featured Suppliers:



As you work to meet the exact specifications and requirements of your critical environment, lean on our broad collection of contamination control products for efficiency and compliance. From apparel to cleaners and disinfectants, you'll find the products you need to minimize contamination.

**We have safety down to a science.**

Visit [fishersci.com/contamination-control](https://fishersci.com/contamination-control) or [fishersci.ca/contamination-control](https://fishersci.ca/contamination-control) to explore contamination control and controlled environments products.



Safety

| Protecting What Matters Most


 **fisher scientific**  
part of Thermo Fisher Scientific

# Earn Fisher Scientific Buck\$ and Experience Excellence

**METTLER TOLEDO**

As your needs continue to grow, so does METTLER TOLEDO's portfolio of laboratory solutions. From XPR balances to the new MyBrix handheld refractometer, we continue to develop innovative solutions that are easy to integrate and designed to meet the demands of the modern lab, providing you with an experience that can only be defined as extraordinary.

Find your perfect solution and earn Fisher Scientific Buck\$ today by visiting [fishersci.com/mettlertoledo](http://fishersci.com/mettlertoledo) or [fishersci.ca/mettlertoledo](http://fishersci.ca/mettlertoledo).



G20S Titrator

XPR Analytical Balance

SevenDirect pH Meter

MyBrix Refractometer

## 70% Isopropyl Alcohol



### Decon CiDehol 70

**Decon CiDehol 70 is a 70% (v/v) isopropyl alcohol (IPA) solution that meets USP specifications.**

Use CiDehol 70 to clean surfaces and hoods in microbiology and tissue culture labs or surfaces and process equipment in pharmaceutical, biotechnology, and medical device facilities. Since IPA solutions evaporate completely, no rinsing is required.


Filtered to 0.2 µm, CiDehol 70 is packaged in ready-to-use trigger spray bottles and other convenient sizes.

CiDehol 70 is also available in a sterile form as CiDehol 70 ST. It's made using WFI in ISO 5 cleanrooms, filtered to 0.2 µm, double-bagged, irradiated, and tested for endotoxin and sterility. Each case comes with lot-specific documentation of QC, irradiation, endotoxin, and sterility.

Visit [fishersci.com/decon](http://fishersci.com/decon) or [fishersci.ca/decon](http://fishersci.ca/decon) to find more Decon products.



Packaging	CiDehol 70			CiDehol 70 ST	
	Size	Quantity	Cat. No.	Quantity	Cat. No.
Trigger Spray Bottles	12 oz. (0.35 L)	—	—	16/Case	<b>04-355-500</b>
	16 oz. (0.47 L)	Each 12/Case	<b>04-355-63</b>	12/Case	<b>04-355-58</b>
	32 oz. (0.94 L)	Each 12/Case	<b>04-355-250</b>	12/Case	<b>04-355-85</b>
Plastic Bottles	1 gal. (3.8 L)	Each 4/Case	<b>04-355-71</b>	4/Case	<b>04-355-80</b>
Plastic Drum	55 gal. (209 L)	Each	<b>04-355-42</b>	—	—



## The Bench-Size Microplate Imaging and Analysis Workhorse

### Agilent BioTek Cytation C10 confocal imaging reader

From confocal, widefield, and live cell imaging to multimode reading, the Agilent BioTek Cytation C10 confocal imaging reader delivers crisp, high-resolution images and the ability to image and analyze 3D samples. Plus, without having to schedule time at a core lab, you'll have more control over your workflow. Small and powerful – a perfect fit in any lab.

Visit [fishersci.com/agilent](https://fishersci.com/agilent) or [fishersci.ca/agilent](https://fishersci.ca/agilent) to shop and learn more.



The new Agilent BioTek Cytation C10 confocal imaging reader

# Cleanroom-Compatible CO<sub>2</sub> Incubators for Cell Therapy

Cell therapy process development and manufacturing requires improved technologies for CO<sub>2</sub> incubators as well as other cleanroom equipment. These products must use materials optimized for stringent cleaning protocols and reduce the potential for microorganism and non-viable particulate contamination. CO<sub>2</sub> incubators must also provide optimal growth conditions for patient cell samples and adjustable shelving configurations to maintain high throughput per footprint.

## Important Considerations

### Cleaning and Disinfection Compatibility

CO<sub>2</sub> incubators with a brushed 304 stainless-steel exterior are recommended for cell therapy development and production. More resistant than a traditional painted steel exterior, this finish can tolerate chemical disinfectants and vaporized hydrogen peroxide (VHP) procedures.

Likewise, electropolished stainless steel components in the chambers have fewer microscopic structures that can become havens for microbial growth. Moisture-resistant housings and casings (rated IP54) and silicone-sealed touchscreen displays protect the electronics and can tolerate common cleanroom protocols.

### Low Particulate Emission and Microbial Contamination Prevention

Particulate management remains a primary concern in cell therapy manufacturing since cell-based products for human injection cannot be easily purified. Although not the only contaminant of concern, microbial contamination remains the number one reason for biological pharmaceutical recalls.

Choose a CO<sub>2</sub> incubator that offers active particle control for ISO Class 5 and Grade A/B cleanroom environments. Your CO<sub>2</sub> incubator should also offer contamination

control as well as control parameters to help cells express their optimal characteristics and potency.

- In-chamber circulating fans provide fast recovery, improve uniformity following door openings, and offer in-chamber HEPA filtration
- Carefully designed air circulation systems help parameters recover within 10 minutes after a thirty-second door opening
- Combined air circulation and HEPA filtration will help you capture particulates of all sizes and reach ISO Class 5 conditions in five minutes or less
- On-demand dry heat sterilization cycle tested by an independent, third-party laboratory will help to verify that it reaches a 12-log sterilization assurance level (SAL)
- Multiple-point temperature mapping demonstrates consistency throughout the chamber

### Certificates and Documentation

Certified cleanroom-compatible CO<sub>2</sub> incubators are tested against ISO 14644-1 requirements and adhere to air quality limits within grade A/B cleanroom standards. In addition, each unit should undergo complete end-of-line testing. The manufacturer should also provide certificates and documentation to facilitate and support the audit process for on-site qualification, including:

- Product certificates, specifications, and materials of composition

- Factory acceptance test documentation
- Clear guidance for cleaning, disinfecting, routine maintenance, accessory and part replacement, and expected performance

### Optimized Space to Help Achieve High Throughput

CO<sub>2</sub> incubators should accommodate appropriately sized culture vessels and produce an environment that promotes proper growth and the expression of desired receptors and optimum critical quality attributes (CQAs). Accessory incubator shelving systems can also help maximize limited incubator space, providing up to one and a half times more chamber capacity for higher output per footprint.

Visit [fishersci.com/co2-incubators](http://fishersci.com/co2-incubators) or [fishersci.ca/co2-incubators](http://fishersci.ca/co2-incubators) to learn more.

Content provided by:

**thermo** scientific





# Bead Mill 24 Homogenizer



- DNA/RNA extraction
- Tissue homogenization
- Enzyme isolation
- Protein purification

For best results, use with Fisherbrand pre-filled bead tubes.

Type	Particle Size	Material	Quantity	Cat. No.
Hard Tissue Grinding	2.4 mm	Metal	50/Pack	15-340-151
Tough Microorganism Lysing	0.5 mm	Glass	50/Pack	15-340-152
Soft Tissue Homogenizing	1.4 mm	Ceramic	50/Pack	15-340-153
Hard Tissue Homogenizing	2.8 mm	Ceramic	50/Pack	15-340-154



Same-Day Shipping  
on Even More Products



## Let Us Deliver for You

Look for the **edge** icon when shopping on our website to instantly identify items available for standard **same-day shipping** when ordered by 2 p.m.\*

Find even more ready-to-ship products with ongoing updates from real-time inventory.

Visit [fishersci.com/edge](https://fishersci.com/edge) or [fishersci.ca/edge](https://fishersci.ca/edge) to learn more and shop products.

\*Restrictions may apply.

FEATURED ARTICLE



# Ancestral Antiviral Proteins

May Protect Against Infection

By Gina Wynn

We can thank our ancestors for many of our inherited physical traits, including our eye color, height, and complexion. Now we may be able to add immunity to viral infections to that list, according to a recent Cornell University study.

By Gina Wynn

When examining cultured human cells in the laboratory, researchers found that proteins from ancient viral DNA were passed down from our ancestors. These human endogenous retroviruses (HERVs) may provide immunity by preventing modern viruses from entering cells that could become hosts. Previous studies have already documented this occurrence in mice, chickens, cats, sheep, and other primates.

“The results show that in the human genome, we have a reservoir of proteins that have the potential to block a broad range of viruses,” said research team member Cédric Feschotte, PhD, Cornell professor of molecular biology and genetics, as reported in the *Cornell Chronicle* article “Ancient viral DNA in human genome guards against infections.”

Feschotte’s former graduate student, John Frank, PhD, who is now a Yale University postdoctoral researcher, was the first author of the study. The group’s findings were published in the *Science* article “Evolution and Antiviral Activity of a Human Protein of Retroviral Origin.”

## A Pool of Protective Proteins

Further investigation could lead to the discovery of antiviral proteins that could help develop treatments that don’t cause autoimmune side effects. Because the proteins are already incorporated into human DNA, they aren’t seen as foreign to the body and don’t prompt an immune response.

Approximately eight percent of DNA in the human genome is made up of HERVs, which is at least quadruple the amount of DNA in the genes that code for proteins. Our genome could be harboring a vast defense system with significant implications for healthcare.

## Gaining Genomic Access

Viruses, in the form of proviruses, can introduce their RNA into a host cell, where it is converted to DNA and integrated into the host’s genome. The hijacked cell then follows the genetic

instructions from the virus to replicate it. If this occurs in germ or reproductive cells, the viral DNA becomes a permanent part of the genetic code and is passed down through generations.

“The results show that in the human genome, we have a reservoir of proteins that have the potential to block a broad range of viruses.”

**Cédric Feschotte**


*Professor of Molecular Biology and Genetics  
Cornell College of Agriculture and Life Sciences*

HERVs are evidence of viruses that became fixed in the human genome millions of years ago. They are remnants of the countless pandemics to which humans have been exposed throughout history. Scientists believe these viruses were widespread in ancient human populations since they have also been found in chimpanzee, gorilla, and other primate genomes, according to the Aidan Burn article in *Scientific American* “How the Ancient Viral DNA in Our Genome Affects Disease and Development.”

## The Key to Entry

To enter a cell, a virus first needs to be admitted by a viral envelope or spike protein. The envelope protein binds to a receptor on the cell surface, like a key in a lock, and opens the cell to admit the virus. For their study, Frank and team used computational genomics to scan the human genome to find viral envelope proteins that they could investigate.

After cataloging the potential retroviral envelope protein-coding sequences that potentially retained receptor binding activity, they identified which genes were actively expressing retroviral envelope gene products. The team found clear evidence of expression in the early embryo and germ cells. They also noted a subset of antiviral proteins that are expressed in immune cells upon infection.



## Ancestral Antiviral Proteins May Protect Against Infection

### Testing with Suppressyn

To test the susceptibility of different cell types to viral infection, the group focused on suppressyn (SUPYN), an antiviral envelope protein that is known to bind the alanine-serine-cysteine transporter 2 (ASCT2) receptor. SUPYN is prevalent in very early human embryonic development and placental cells, and ASCT2 is the cellular gateway for type D retroviruses.

Because placental cells are commonly targeted by viruses, Frank and his colleagues experimented with exposing human placental-like cells to RD114, a type D retrovirus that typically infects domestic cats and other felines. Although other human cell types not expressing SUPYN were easily infected, the placental and embryonic stem cells that did express SUPYN were not affected by the virus. When the team removed SUPYN from the cells, they became infected with RD114; when SUPYN was reintroduced, viral resistance returned.

Frank and his colleagues took their experimentation a step further and introduced SUPYN to embryonic kidney cells that don't typically express it. Normally the cells are susceptible to RD114, but with the added SUPYN, they were resistant to the virus.

### Infection Insight

These results demonstrate how a human retroviral protein can deny viruses entry into cells by using antiviral envelope proteins to block cell receptors. They provide insight into how ancient retroviruses in the human genome may protect developing embryos from infection by related viruses. Feschotte hopes to study other envelope-derived proteins in the human genome to see if they have the same antiviral effect on cells, according to the *Cornell Chronicle* article.

### Pushing the Envelope

In a similar study conducted in 2017, researchers recreated an envelope protein that helped the HERV-T virus infect human

cells by binding with the monocarboxylate transporter 1 (MCT1) receptor. Daniel Blanco-Melo, PhD, of The Rockefeller University in New York and colleagues published their findings in the *eLife* article "Co-option of an endogenous retrovirus envelope for host defense in hominid ancestors."

The HERV-T retrovirus spread among our primate ancestors for around 25 million years and became extinct roughly 11 million years ago. By analyzing the genetic remains of HERV-T in the genomes of humans and related primates, the team was able to recreate the ancHTenv envelope protein.

The group's analysis showed that the HERV-T provirus in hominid genomes includes an *env* gene (hsaHTenv) that has been uniquely preserved. Upon further investigation, they found that hsaHTenv caused the depletion of the MCT1 receptor. Because the MCT1 was not available to interact with the ancHTenv envelope protein, HERV-T was unable to enter the cells and they were not infected with the virus. These findings further illustrate the importance of the envelope/receptor relationship in viral transmission.

### Assistance from Our Ancestors

In both these studies, viral envelope protein and receptor pairs associated with HERVs in our ancestral DNA were shown to affect immunity to ancient viruses. The results of both investigations demonstrated that the lack of either the viral envelope or the receptor blocked the viruses from entering and reproducing in cells.

These findings may have major implications for researchers testing for and developing disease treatments, especially for HIV and other viruses that integrate into host cell DNA, some cancers, and possible future pandemics.

*Gina Wynn is a Thermo Fisher Scientific staff writer.*

**Introducing**

# Fisherbrand Multi-Pad pH Indicator Strips

**Dip and Read for Simple, Accurate pH Measurement**

Fisherbrand Multi-Pad pH Indicator Strips provide distinct color separations with clear, bright multicolor matches. Available in seven pH ranges, they're packaged in convenient flip-top vials that include color comparison charts.

Visit [fishersci.com/fisherbrand](http://fishersci.com/fisherbrand) or [fishersci.ca/fisherbrand](http://fishersci.ca/fisherbrand) to shop.



pH Range	Quantity	Cat. No.
0.0 to 14.0, 100/Vial	Each	13-640-521
0.0 to 6.0, 100/Vial	Each	13-640-522
2.0 to 9.0, 100/Vial	Each	13-640-527
4.0 to 7.0, 100/Vial	Each	13-640-523
4.5 to 10.0, 100/Vial	Each	13-640-524
6.0 to 10.0, 100/Vial	Each	13-640-525
7.0 to 14.0, 100/Vial	Each	13-640-526

## Unparalleled Choice and Convenience

The **Fisher Scientific** channel offers more than **1.7 million products** from over **9,000 suppliers globally** through a leading e-commerce platform, powerful distribution network, and trusted team of sales representatives, scientific specialists, and other experts. In addition to supplying products to advance science, our team provides equipment and instrument services, stockroom management, production solutions, and much more.

Whether you're engaged in scientific research, healthcare, safety, or education, the Fisher Scientific channel provides the products and support you need to succeed.

Visit [fishersci.com](http://fishersci.com) or [fishersci.ca](http://fishersci.ca) to shop and learn more.

---

# Complete Your Chromatography from End to End

No matter your chromatography method, find reagents, solvents, standards, columns, cartridges, and other products to support your workflow every step of the way.

## Liquid Chromatography

Shop sample preparation and separation products for small and large molecule analysis, HPLC/UHPLC and benchtop solutions, and reliable solvents to help you optimize performance.

## Life Sciences Chromatography

Determine protein characteristics with buffers, reagents, protein purification resins, columns from top brands, chemicals for separation, standards, and more.

## Gas Chromatography

Browse headspace vials, inserts, and closures, GC columns, and capping equipment.

## Chromatography Reagents

Shop solvents, buffers, standards, reagents, bulk chemicals, and so much more to support your work.

Visit [fishersci.com/chromatography](https://fishersci.com/chromatography) or [fishersci.ca/chromatography](https://fishersci.ca/chromatography) to shop these categories and learn more.



# Sustainability in the Laboratory

## Small Changes Are Adding Up to a Better Tomorrow

How will you effectively reduce waste, space, and costs while minimizing disruption in your laboratory?

You're under immense pressure to create a safer, healthier, and more productive process for your team, customers, and the environment—all while balancing costs and time. With Whirl-Pak, small changes can add up to help you meet the challenge.

Successfully reduce:

- Plastic waste by 5 to 10x
- Emissions by 5x
- Water used throughout life cycle by 10x
- Space by 95%
- Disposal costs

### It's Time for a Change

Sustainability is a word you hear often today, and companies are spending significant time and effort on this topic. Labs are notorious for using large quantities of plastic consumables and should be included in this discussion. Lab plastics are shatterproof, durable, lightweight, single-use products that can reduce cross contamination. They're easy to use and almost ubiquitous. Initiating changes in how they're used in a laboratory setting can be challenging. Standard operating procedures are commonplace, and scientists follow them day in and day out to maintain order and consistency.

### Reduce, Reuse, Recycle

Sustainability related to plastic consumables generally defaults to a discussion about recycling. Can products be safely recycled and is this activity convenient? Although important, recycling is not the only option to consider.

Another issue is the amount of plastic we use, which can often be reduced from the start. As a part of its overall waste management strategy, the U.S. EPA developed a Waste Management Hierarchy: source reduction is the most preferred tool, followed by materials reuse, and finally, recycling. Simply put: reduce, reuse, recycle. The order here is important because the greatest impact on sustainability is through reduction in overall use.

### Sample Collection Containers Are Leading the Way

One way of reducing overall plastic consumption in the lab is substituting flexible sampling containers for rigid plastic. Whirl-Pak sterilized sampling bags were developed over 60 years ago to transport liquid milk samples to the lab. Today, these flexible packages are preferred in many industries and across more than 75 countries. Not all industries, however, have embraced this sampling format.

One industry ripe for a sampling procedure change is the water industry. Of the hundreds of millions of water samples taken each year, the majority are collected in rigid, single-use plastic containers. On average, these containers use 5 to 10 times the amount of plastic in a comparably sized Whirl-Pak bag. While this more sustainable solution is readily available, initiating change has been difficult.

### Evaluating Environmental Impact

Whirl-Pak Filtration Group conducted an evaluation of the environmental impact of using 100 mL and 500 mL sterile sampling bags instead of the rigid plastic containers currently used. The study used Ecolmpact-COMPASS software, which is widely used for life cycle assessments in the packaging industry. A life cycle assessment characterizes the impact associated with sourcing, manufacturing, distributing, using, and disposing of a given product. It gives you an idea of which package may have a lower environmental footprint, allowing you to make better sustainability decisions.

The study compared Whirl-Pak bags made from linear low-density polyethylene (LLDPE) to rigid containers made from polypropylene (PP), polyethylene terephthalate (PET), or high-density polypropylene (HDPE). Environmental impact variables included



Content provided by:



fossil fuel use, greenhouse gas emissions, and water use. Both sizes of Whirl-Pak bags produced similar results; only results for the 100 mL size are discussed below.

The study assumed that one million containers were shipped a distance of 1,000 km. Data was segregated by material usage, manufacturing, transport, and end of life. In all product comparisons, the Whirl-Pak sterile sampling bags were significantly more sustainable than the three rigid containers commonly used in water collection, including less fossil fuel use, lower greenhouse gas emissions, less water usage, and an overall lower environmental impact throughout the product's life cycle.

## Cost Considerations

Even when a sustainable alternative is identified, cost considerations can affect the decision to implement changes. Budgetary constraints can limit even a highly motivated lab manager. One must consider overall costs and also hidden savings associated with shipping, storage, and waste disposal costs. Savings may be magnified if shipping long distances is required.

A focus on small, yet meaningful changes can make a measurable impact. Replacing rigid plastic sample collection vessels with flexible solutions like Whirl-Pak bags can be the beginning of your journey toward a more sustainable work environment.

Visit [fishersci.com/whirl-pak](https://fishersci.com/whirl-pak) or [fishersci.ca/whirl-pak](https://fishersci.ca/whirl-pak) to learn more.



# First Plants Grown in Lunar Soil May Lead to Farming in Space

By Mark Miller

Researchers at the University of Florida Institute of Food and Agricultural Sciences are the first to successfully grow plants in lunar soil brought back by the Apollo space missions. The lab-grown garden—part of research supporting the ambitious Artemis program at the National Aeronautics and Space Administration (NASA)—marks an important step toward growing plants and food on the moon and possibly worlds beyond.

## Tough Growing

Dirt from the moon is a challenging environment for plant growth. According to an article in *Science News*, lunar soil, known as regolith, is a fine powder full of metallic iron, which, unlike oxidized iron, isn't beneficial to plants. It also contains tiny shards of glass as a result of the moon's surface being continually struck by micrometeorites. What it doesn't contain much of is nitrogen, phosphorus, or other nutrients.

The research team chose thale cress or *Arabidopsis thaliana* to grow in this unwelcoming world. It's a relative of mustard greens and other cruciferous vegetables like broccoli and cauliflower. A small, hardy plant that grows easily, thale cress is also the model organism used in multiple studies of plant biology, which means scientists have a good understanding of its genetic makeup.

The team combined thale cress seeds and water with gram-sized samples of soil from the Apollo 11, 12, and 17 missions. These were placed in terrarium boxes under LED lights in a cleanroom. They also planted control samples in terrestrial volcanic material and added a nutrient to all of the samples daily.

"After two days, they started to sprout!" said Anna-Lisa Paul, PhD, a professor of horticultural sciences at the University of Florida and part of the research team, in a report from NASA. "Everything sprouted. I can't tell you how astonished we were! Every plant—whether in a lunar sample or in a control—looked the same up until about day six."

## Better Plots and Plants

As Paul's reaction indicates, both sample sets did not grow equally well. After about a week, the growth of the lunar samples slowed. They had stunted roots and leaves and an unhealthy, purplish pigmentation.

To better understand what had happened, the team conducted genetic analysis by studying the ribonucleic acid (RNA) from the plants grown in the moon soil. RNA transforms information from deoxyribonucleic acid (DNA) into proteins that carry out many

of a living organism's biological processes. Examination of the RNA sequence of the plants' cells revealed that the thale cress was reacting as it typically does when the environment contains too much salt or heavy metals.

The plants had the most trouble growing in the soil samples from Apollo 11, perhaps because these were drawn from the moon's Sea of Tranquility, according to an *NBC News* story. This region, according to the story, is an older surface of the moon and has experienced greater cosmic radiation and solar wind.

A correlation between the age of the lunar soil and the performance of the plants may help future moon gardeners know where—or where not—to collect their planting soil. It may also help to select plants better suited to moon dirt. "Maybe spinach plants, which are very salt-tolerant, would have no trouble growing in lunar regolith," wondered Paul in *Science News*.

## One Small Step

Although the lunar samples didn't produce growth as well as the volcanic soil, the research is a dramatic initial success that significantly advances the role of plant life in space exploration and habitation.

Plants can be used as model organisms to help study gravity, radiation, and other biological phenomena in space, according to a paper the University of Florida team published in *Nature*. It points out that plants help us understand how biological organisms adapt away from Earth and could provide food and oxygen to support long-term habitation of space and even extraterrestrial surfaces. Using moon regolith could prove much more practical than transporting enough water to support hydroponic agriculture.

## Returning to Space

The work done at the University of Florida is picking up where NASA left off. Over the course of their missions, Apollo crews returned to Earth with 842 pounds of moon rocks and soil. The plantings are making good use of these materials. "Here we are, 50 years later, completing experiments that were started back in the Apollo labs," said Robert Ferl, PhD, another member of the University of Florida team, in the NASA report.

The research plays a central role in the success of the Artemis program, which aims to return astronauts to the moon and establish a long-term presence there. Ultimately, its lunar missions will serve as a leaping-off point for human exploration of Mars.

*Mark Miller is a Thermo Fisher Scientific staff writer.*

# Accuracy, Simplicity, and Durability



## KNF SIMDOS Liquid Diaphragm Pumps

Accurately meter and dose with small, quiet, self-priming pumps that can operate dry without damage. Intuitive controls and splash-proof (IP65) housings combine straightforward operation with reliable durability. These pumps feature:

- Rotary knobs with displays to efficiently set dose volume, flow rate, and other parameters
- Simple calibration for accurate and repeatable liquid metering and dosing
- Chemically resistant construction to tolerate even highly aggressive chemicals




Model	Metering Rate (mL/min.)	Dose Volume (mL)	Head	Diaphragm	Valves	Cat. No.
FEM1.02KT.18S	0.03 to 20	0.03 to 999	PP	PTFE-Coated	FFKM	13-880-907
FEM1.02FT.18S	0.03 to 20	0.03 to 999	PTFE	FFKM	FFKM	13-880-909
UFEM1.10KT.18S2	1 to 100	1 to 999	PP	PTFE-Coated	FFKM	13-880-919
UFEM1.10FT.18S2	1 to 100	1 to 999	PTFE	PTFE-Coated	FFKM	13-880-921

## Reaching Net-Zero Goals Together



As part of Thermo Fisher Scientific, the Fisher Scientific channel is committed to achieving net-zero emissions by 2050 and helping you reach your goals, too. Our sustainability program offers an ever-expanding portfolio of products and resources to help you work more sustainably, including:

### The Greener Choice Program

Easily find lab products that align with your sustainability efforts. Simply look for the  leaf icon in search results to identify products that help protect and preserve the environment.

### Recycling Programs

Take part in recycling programs of all shapes and sizes that support your company's sustainability initiatives.

Visit [fishersci.com/sustainability](https://fishersci.com/sustainability) or [fishersci.ca/sustainability](https://fishersci.ca/sustainability) to learn more about sustainability-related products and resources available through the Fisher Scientific channel.

# 5,000 solutions. 1 washer.

Keeping up with cleaning all of your labware is a challenge. Labconco's glassware washers are up for it, with the ability to wash over 5,000 different labware configurations. Powerful software lets you store, customize, and protect programs and data. Labconco. Let us clean, so you don't have to.

Learn more at [fishersci.com](http://fishersci.com) or [fishersci.ca](http://fishersci.ca)  
Search: **FlaskScrubber** or **SteamScrubber**




**LABCONCO**



## Achieve Your Supplier Diversity Goals

The Fisher Scientific channel can help you manage and grow your spending with small and diverse businesses through our Supplier Diversity program. We'll work with you to build a custom plan based on your unique needs—and help you measure and track it across specific diversity categories based on federal designations.

You can also look for the Supplier Diversity program icon  on our website to quickly identify suppliers and products that meet federal requirements for small and diverse businesses.

Visit [fishersci.com/supplierdiversity](http://fishersci.com/supplierdiversity) to learn more and start building your customized plan.

*Available in U.S. only.*




# Trusted chemical protection

## Tychem® 2000 garments & accessories

DuPont™ Tychem® 2000 garments and accessories provide at least 30 minutes of protection against more than 40 chemical challenges. Select from various designs, including coveralls, aprons and sleeves. Protect lab workers with the trusted chemical protection of Tychem® 2000.





# Successful Cell Therapy Manufacturing

By Iva Fedorka

Cell- and gene-based therapies are transforming cancer treatment. In the United States and Europe, cell therapies debuted with the use of chimeric antigen receptor (CAR) T-cell therapies, which have become the standard of care for some types of blood cancers.

Creating these therapeutics involves collecting T-cells from patients or healthy donors, genetically modifying, and then administering them to patients to rid the body of cancer cells.

However, producing these products presents new challenges. The first cellular therapies were developed without considering the regular production of these advanced therapy medicinal products (ATMPs).

## The Manufacturing Process

Cell therapy manufacturing begins in an apheresis or clinical facility where the cells are collected and ends when the product is administered to the patient. Between those steps, a complex series of procedures takes place.

This process is very different from traditional pharmaceutical manufacturing. Most pharmaceutical products start with standardized materials and produce a fully characterized product. Cell therapy provides patients with personalized products based on their diseases, genetics, and medical histories, adding complexity and variability to their production.

## Manufacturing Challenges

Manufacturers must produce a consistently safe and effective living product from highly variable living source material. They must produce consistent quality and performance; comply with regulatory and documentation requirements; maintain safety and quality control; manage supply chains and logistics; and develop bioprocessing scale-up and scale-out capabilities.

The success of producing an ATMP depends on the raw materials, manufacturing process controls, and safety and quality control procedures.

## Using Appropriate Raw Materials

Perform a risk assessment of your raw materials to determine which ones may need to be altered or replaced. Choose raw materials that meet the standards and purity requirements for the end-product for all stages of development and manufacturing.

The responsibility for raw material quality and purity is shared by the supplier and manufacturer. Assess materials according to source and identity, purity, safety, and suitability as specified by USP <1043> Ancillary Materials (AM) for Cell, Gene, and Tissue-Engineered Products.

The Standard classifies materials as:

- Tier 1—low-risk, highly qualified, and regulated products
- Tier 2—low-risk, well-characterized, GMP, designed for use as AMs
- Tier 3—moderate risk, not intended for use as AMs
- Tier 4—high risk, biologically variable or contain harmful impurities and/or animal materials

As often as possible, choose Tier 1 and Tier 2 materials. Work with your supplier to identify alternative products for Tier 3 and Tier 4 materials, including substitutions or recombinant versions of animal components, additional material characterization testing, traceability documentation, and upgraded manufacturing methods like aseptic filling.

## Finding Manufacturing Solutions

Producing cell therapies is labor-intensive, subject to contamination at multiple points, complex, difficult to monitor, and has zero failure tolerance. Verify that any equipment and instruments for “research use only” (RUO) or for in vitro diagnostics (IVD) are suitable for producing injectable therapeutics.

Look for ways to introduce modular or multifunctional equipment and instruments, closed systems, and automation to help address these issues. Some manufacturers have developed products specifically designed for cell therapy manufacturing, including direct replacements for RUO or IVD versions.

## Characterization and Product Release

Cell therapy products have a short shelf life and cannot be fully characterized or purified pre-administration. However, some tests are essential for lot release and regulatory compliance for chemistry, manufacturing, and controls (CMC).

The U.S. Food and Drug Administration (FDA) Center for Biologics Evaluation and Research (CBER) Office of Tissues and Advanced Therapies (OTAT) oversees CMCs. The FDA requires information for Investigational New Drug Applications (INDs), including safety, identity, quality, purity, and strength/potency testing.

Chain of custody and identity must also be maintained throughout the manufacturing process. Failure to properly document and control sample provenance can have serious patient consequences.

## The Future of Cell and Gene Therapies

A new generation of CAR T-cell therapies has the potential to identify solid tumor cells in patients with prostate, breast, stomach, rectal, and other non-hematological malignancies. Stem cell therapies may one day be used to treat autoimmune, Alzheimer's, and Parkinson's diseases.

Visit [fishersci.com/production](https://fishersci.com/production) to learn more about Fisher Scientific production products and value-added services.

*Iva Fedorka is a Thermo Fisher Scientific staff writer.*

# Confocal Imaging and Analysis Using the Agilent BioTek Cytation C10

Since 2012, cells and tissues from many sources have been cultured in three dimensions (3D), enabling more complex biological models. The use of these 3D models has been increasing in medical research, precision medicine, disease modeling, and drug discovery efforts.

These 3D models are being used to simulate the native microenvironments found in organisms and are believed to provide a more accurate assay model in some instances. 2D models, consisting of a monolayer of cells, have been in use for decades and have provided extensive meaningful data.

However, 2D limitations became more apparent with the recent development of complex organoids, tissues, or tumoroids. These tissues, comprised of one or more cell types and often based on stem cells or patient-derived samples, have now become the focus of many studies.

## The Cytation C10

The ability to perform a variety of assays on a single, compact instrument is advantageous, especially given the limited bench space in many laboratories. The Agilent BioTek Cytation C10 confocal imaging reader can perform microplate reading, widefield imaging, and confocal imaging, offering the means to gather data for today's complex studies.

The Agilent application note “Confocal Imaging and Analysis of Spheroids for

Determination of Dose Response During Drug Treatment” by Peter J. Brescia describes the use of automated imaging to determine cell number in spheroids to evaluate drug dose responses.

## Study Materials and Methods

The study discussed in the application note used CT116-H2B-GFP cells; the H2B-GFP construct provided a constitutively expressing nuclear marker for quantification by cell count for comparison to a nuclear stain during live cell analysis.

The cells were seeded at a relatively low density of 500 cells/well in a 96-well, round bottom ULA imaging plate (Cat. No. 15-100-173) and grown for three to four days to produce spheroids of roughly 100  $\mu\text{m}$  in size.

Staurosporin, a potent, nonselective protein kinase inhibitor, was added to induce apoptosis. Propidium iodide (PI) was also added to monitor the kinetics to determine optimum timing (12 hours post-addition produced adequate apoptotic activity for analysis).

## Study Imaging Procedure and Analysis

To minimize the data set, spheroids were sized to be captured in a single field of view (20X objective) and a beacon was used to correct for positioning variability in the wells. A 60  $\mu\text{m}$  pinhole spinning disk provided a compromise between

acquisition time and signal intensity for confocal imaging. Widefield images were captured in both GFP and TRITC channels to compare to confocal images.

The instrument was set for “Scan,” followed by “Autofocus,” which identified the approximate center of the spheroid along the z-axis. Each spheroid was captured as a z-stack using 11 steps at 10  $\mu\text{m}$ , regardless of the z position in the well. Images were processed in several steps, including a z-projection for comparing single z-plane images at the approximate center of the spheroid and the entire z-stack for both imaging modes.

Automated cell counting was performed using confocal and widefield imaging mode and both single plane and z-stack projected images were analyzed.




Figure 1. Spheroid assay workflow incorporating automated image acquisition and analysis (Agilent BioTek Cytation C10 confocal imaging reader and Agilent BioTek Gen5 microplate reader and imager software).



Content provided by:




Figure 2. Spheroid captured as a z-stack using 11 steps at 10  $\mu\text{m}$  intervals. Analysis performed on either a z-projection of all images or a single image representing a plane through the center of the spheroid.




Figure 3. Examples of automated image acquisition and cell count analysis of spheroids by confocal and widefield fluorescent microscopy using z-projected images. Cell masking (yellow outlines) where GFP cell count indicates live cells and PI indicates dead cells.





Figure 4. EC<sub>50</sub> determinations for each imaging mode, single versus z-projected images.

Since spheroids are increasingly used to represent the complex microenvironments of biological systems, the ability to automate image capture using multiple imaging modes and perform automated image analysis can increase throughput and enable comparative analysis.

The study results confirmed that confocal optics enable more accurate segmentation of labeled cells within spheroid samples than widefield optics. Cell count

determination using a single z-plane image provided comparable results to methods relying on z-stack image acquisition and processing.

The study is just one example of the utility and flexibility of the Agilent BioTek Cytation C10 confocal imaging reader and Agilent BioTek Gen5 microplate reader and imager software to capture and analyze multiple parallel data sets.

Visit [fishersci.com/agilent](http://fishersci.com/agilent) or [fishersci.ca/agilent](http://fishersci.ca/agilent) to learn more.



Save All Year Long

# Spotlight

Get Up to **50% Off**

Spotlight is your destination for special offers from top brands available through the Fisher Scientific channel. Stop back often to find new savings on:

Lab essentials | Equipment & instruments | Chemicals & chromatography | Safety | Life sciences

Visit [fishersci.com/spotlight](http://fishersci.com/spotlight) or [fishersci.ca/spotlight](http://fishersci.ca/spotlight) to learn more and start saving.

## Oil Free, Low Maintenance

*fisherbrand*

### Fisherbrand MaximaDry Pumps

Eliminate downtime to add or change pump oil. Enjoy quiet, reliable operation without the mess of oil or the operating costs and environmental problems associated with water aspirators.

- **Corrosion-resistant materials**—The thick, solid PTFE resin heads will not peel or corrode when exposed to most aggressive solvents and the sturdy FFPM valve system handles occasional liquids without degrading performance
- **Quick and easy setup**—Be up and running in minutes; everything connects together without tools or fittings, just add 3/8 in. I.D. tubing (13-991-206)
- **Compact, stylish design**—For convenient access, your MaximaDry Oil-Free Pump fits easily where you need it—on your bench or inside your fume hood
- **Versatile**—Suitable for rotary evaporation, degassing, filtration, solid phase extraction, fluid aspiration, vacuum ovens, and other applications



13-880 Series

Flow Rate	Vacuum	Pressure	Dimensions (L x W x H)	Cat. No.
20 L/min.	75 torr	15 psig	13 x 6.5 x 8.5 in.	<b>13-880-14</b>
40 L/min.	65 torr	15 psig	13 x 6.5 x 8.5 in.	<b>13-880-16</b>
20 L/min.	6 torr	15 psig	15 x 6.5 x 8.5 in.	<b>13-880-18</b>
35 L/min.	6 torr	15 psig	16 x 7 x 9 in.	<b>13-880-20</b>
35 L/min.	1.5 torr	15 psig	16 x 7 x 9 in.	<b>13-880-22</b>

# The Cost of Air and the Impact on the Laboratory

Engineered controls, like fume hoods, are an integral part of most laboratory facility planning.

Fume hoods have a significant impact on building infrastructure, are energy-intensive, and basically, remove conditioned air from the building. The extracted air then needs to be brought back into the building for the HVAC to adequately balance temperature, humidity, and pressure. Ducted fume hoods impact the building's need for larger air handlers, boilers, chillers, cooling towers, and exhaust fans due to the amount of air that needs to be processed.

Studies have calculated that laboratories consume an average of 66 percent or two-thirds of a building's total energy load while occupying only about 25 percent of the facility's usable space.

How can we address this issue, what are some possible solutions, and how do we address these very real problems while achieving carbon neutrality?

Alternative engineered control systems (ECS), including ductless filtering fume hoods, are one answer.

## Cascade Effects

When a lab is designed with ductless filtering fume hoods (DFH) as the primary ECS, the impact includes carbon reduction. As we aim for carbon neutrality by 2030, the DFH plays a significant role in achieving this goal. Obvious savings can be found in reducing peak exhaust output, but that reduced exhaust air output allows buildings to achieve much more.

Let's look at the impact of an Erlab project at Bristol Community College in Massachusetts.

The Bristol high-performance laboratory was originally designed to include 22 ducted fume hoods, which had a combined intake and exhaust of 70,000 cubic feet per minute (CFM). It required three air-handling units (AHU), one with run-around loop heat recovery and one with enthalpy wheel energy recovery.

The proposed Erlab Zero Net Energy (ZNE) design replaced that equipment with:

- 13 filtered hoods
- 4 ducted hoods
- CFM: 24,000 intake and exhaust
- 2 AHUs
- Enthalpy wheel recovery

Along with the air reductions, the ZNE design included a combination of ground- and air-source heat pumps, enthalpy heat recovery wheels, fan coil units, centralized Indoor Air Quality (IAQ) monitoring, and natural ventilation. It reduced mechanical, electrical, and plumbing (MEP) to just 14 percent of the gross square footage (GSF). Additionally, the PV arrays, half of which would have supplied the 22 ducted hoods, would be used to recover energy to help achieve ZNE.

## A Holistic Approach to Safety

Fume hoods and containment devices should be designed with safety as the priority. When considering ductless filtering fume hoods, chemical handling needs must be assessed to determine whether they are suitable for filtration. If so, a filter life cycle must be proposed based on analysis using AFNOR NF X 15 211, ANSI Z9.5-2022, CSA Z316.5-2020, and NFPA 45-2023 edition standards.

The Bristol project considered many factors, including hood filtration performance. The hood's chemical listing booklet must provide the overall retention capacity of the filter prior to any release (under 1 percent) of the chemical threshold limit value (TLV). This reflects the life cycle of the filter and the maximum time between replacements.

Under normal operating conditions, a ductless filtering fume hood must also be able to guarantee user protection during two operating phases: detection and safety. One misconception about ductless filtering fume hoods is that they can be

used for only small amounts of chemicals or for odor mitigation. In reality, Erlab hoods are approved for 80 percent of pharmaceutical, organic chemistry, agriculture, flavor and fragrance, and other common applications.

Safety does not stop with the fume hood but continues throughout the chemical life cycle. Proper controls must be in place to protect the breathing zones of people working in the lab. These are achieved through a combination of filtering storage, whole room air filtration, integrated filtration packs for safety cabinets, and filtering fume hoods.

The various safety aspects should be monitored: face velocity, filtration efficiency, and ambient air pollution. Continuous monitoring provides critical safety metrics and helps establish more effective lab safety protocols.

## Conclusion

While we may be reluctant to change, we must consider what's best for the environment. Embracing new technologies can have a significant impact on the future. As we increase our understanding of the impact of our actions, we can make changes to improve safety and product performance.

Ductless filtering fume hoods can improve lab design efficiency and increase your flexibility as your lab changes and grows. Filtration protects you wherever chemicals exist and supports your facility's efforts to reduce its carbon footprint.

Erlab offers 54 years of air filtration technology experience, investment, and dedication to help protect your safety, and let you breathe easier.

Content provided by:



# Making a Difference in Your Workflows

## Greiner Bio-One CELLMASTER Cell Culture Roller Bottles

### Available Now

Choose from:

- Polystyrene or PET
- Growth areas from 850 cm<sup>2</sup> to 4,250 cm<sup>2</sup>
- Standard or filter screw caps

Visit [fishersci.com/greiner](https://fishersci.com/greiner) or [fishersci.ca/greiner](https://fishersci.ca/greiner) to learn more.



### Distributed by Fisher Scientific

#### Promotion Terms and Conditions

These promotions are not available to our current or potential healthcare customers. Healthcare customers include any healthcare practitioner, other healthcare provider or any individual or organization authorized to prescribe, dispense, purchase or influence the acquisition or use of medical devices or supplies for clinical use. Healthcare customers include, without limitation, any "Health Care Practitioner" under the Massachusetts Marketing Code of Conduct Law (defined at 105 C.M.R §907.004) or a "Health Care Provider" under the Vermont Prescribed Products Gift Ban and Disclosure Law (defined at Vt. Stat. Ann. tit. 18, §4631a(8)). By participating in any of these promotions, you warrant that you are not a current or potential healthcare customer. Thermo Fisher Scientific reserves the right to determine your status in the event of a dispute.

Customer acknowledges that these offers may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 U.S.C. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

The factual representations regarding the products accurately reflect the manufacturers' representations to Fisher Scientific. To the extent permissible, Fisher Scientific will make available to its customers the manufacturers' warranties. The sole and exclusive liability of Fisher Scientific shall be refund of the purchase price or replacement of the product. In no event shall Fisher Scientific be liable for incidental or consequential damages. **FISHER SCIENTIFIC EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS, IMPLIED, OR STATUTORY, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.** Prices may change without notice. Prices are U.S. list prices—pricing may vary in Canada. Content and printed prices subject to confirmation by your Fisher Scientific Sales Representative. Telephone calls may be monitored as part of Fisher Scientific's quality process.

© 2023 Thermo Fisher Scientific Inc.  
All rights reserved.

Trademarks used are owned as indicated at [fishersci.com/trademarks](https://fishersci.com/trademarks).

Litho in U.S.A. 22-758-2216 KP 12/22

BN224905594-NA