

## There are more reasons than ever to get your chemicals from Fisher Science Education.

- User-Friendly Reference Guide
- Easier-to-Read Labels
- Industry Standard Storage Plan
- Comprehensive Online Resources

### Classroom-Friendly Options

- Purity Grades for Every Application
  - **Lab (or Educational) Grade** — These economical chemicals are of reasonable purity for situations where no official standards for quality or impurity levels are required
  - **Reagent Grade** — The highest quality commercially available for this chemical; the American Chemical Society (ACS) has not officially set any specifications for this material
  - **ACS Reagent Grade** — Chemicals meet or exceed the latest ACS specifications; designed for analytical applications requiring tighter standards
- Quantities to Fit Your Requirements
  - Range of pack sizes, including smaller quantities to reduce cost and waste

Fisher Science Education is committed to providing you a comprehensive and reliable resource to use when selecting chemicals for your classroom. Each chemical listing includes the following information (where applicable):

Label Key

- 1 Chemical Name
- 2 Synonym(s); Chemical Formula
- 3 Formula Weight
- 4 CAS Number
- 5 Color/Appearance
- 6 Solubility
- 7 Shelf Life
- 8 Department of Transportation (DOT) Classification
- 9 ChemAlert™ Storage Code
- 10 Chemical Grade
- 11 Catalog Number
- 12 Available Sizes
- 13 Price

1	<b>Ethyl Alcohol</b>		
2	Ethanol; C <sub>2</sub> H <sub>5</sub> OH		
3	F.W. 46.07		
4	CAS: 64-17-5		9 <b>FLAMMABLE</b>
5	Color: Clear, colorless		
6	Solubility in Water: Slightly soluble		
7	Shelf Life: 2 years		
8	DOT Class 3, II: <i>Flammable Liquid</i>		
10	<b>70% DENATURED W/MIBK</b>		
11	S25306	500mL	10.15
	S25306A	4L	39.20
	<b>95% DENATURED W/IPA AND METHANOL</b>		
	S25309D	500mL	11.60
	S25309E	4L	52.65
	<b>95% DENATURED W/MIBK</b>		
	S25310	500mL	5.95
	S25310A	3.8L	38.20
	<b>ANHYDROUS, DENATURED W/MIBK</b>		
	S25307	500mL	7.00
	S25307A	1L	13.40
	S25307B	4L	31.50
	<b>ANHYDROUS, DENATURED, W/IPA AND METHANOL</b>		
	S25308	1L	21.65
	S25308A	20L	235.05
	<b>LAB GRADE - 95% DENATURED W/IPA AND METHANOL</b>		
	S25309	500mL	7.50
	S25309C	1L	9.75
	S25309A (SAFE-COTE)	3.8L	76.70
	S25309B	4L	33.05