


HEAVY DUTY BEAKERS

www.fisheredu.com/eisco

- Uniform, thick walls provide increased durability and stability
- Manufactured from high-quality borosilicate 3.3 glass
- Beaded rim and spout for easy mess-free pouring
- Easy to read, white printed graduations
- Superior heat tolerance and chemical resistance, autoclavable at 121°C
- Calibrated within a 5% tolerance



Eisco™ Heavy Duty Beakers		
Code	Capacity	Pack Size
S81276	50 ml	Each
S81277	50 ml	Pack/12
S81278	100 ml	Each
S81279	100 ml	Pack/12
S81280	250 ml	Each
S81281	250 ml	Pack/12
S81282	400 ml	Each
S81283	400 ml	Pack/6
S81284	600 ml	Each
S81285	600 ml	Pack/6
S81286	1000 ml	Each
S81287	1000 ml	Pack/6



Distributed by Fisher Science Education. Contact us today:

In the United States

Order online: fisheredu.com

Call customer service: 1-800-955-1177

Fisher Science Education

4500 Turnberry Drive

Hanover Park, IL 60133