


## TEST TUBE RACKS

[www.fisheredu.com/eisco](http://www.fisheredu.com/eisco)

### Epoxy Coated Test Tube Racks



Code	O.D. Tube (mm)	Format	Number of Places
S13445	18 - 20	4 x 5	20
S13446	15 - 16	4 x 6	24
S13447	18 - 20	4 x 10	40
S13448	20 - 22	4 x 10	40
S13449	10 - 13	6 x 6	36
S13450	15 - 16	6 x 12	72
S13451	18 - 20	8 x 10	80
S13452	22 - 25	8 x 10	80

# Wooden Test Tube Racks



**S07855**



**S07856**

Code	O.D. Tube (mm)	Format	Number of Places	Number of Drying Pegs
S02635	20-22	6x1	6	0
S02636	20-22	12x1	12	0
S07854	18-20	6x1	6	6
S07855	18-20, 22-25	7x1, 6x1	13	6
S07856	20-22	6x2	12	0
S07857	25-28	6x2	12	0

# Polypropylene Test Tube Racks



**S41727**

Code	O.D. Tube (mm)	Format	Number of Places
S41726	10-13	15x6	90
S41727	14-16	12x5	60
S41728	8-20	10x4	40
S27093	22-25	8x3	24
S27094	26-30	7x3	21

- Made of high quality, chemical resistant polypropylene
- 3 piece racks can be taken apart for thorough cleaning
- Snap together design, easy to assemble and disassemble for storage
- Autoclavable at 121°C

Distributed by Fisher Science Education. Contact us today:

**In the United States**

Order online: [fisheredu.com](http://fisheredu.com)

Call customer service: 1-800-955-1177

**Fisher Science Education**

4500 Turnberry Drive

Hanover Park, IL 60133

 **fisher science**  
education  
part of **Thermo Fisher Scientific**