



**Micronova**<sup>®</sup>  
MANUFACTURING INC.

CLINICAL &  
RESEARCH  
LABORATORY

Your strategic supplier for critical cleaning

*Excellence in cleaning since 1984*

Find it at [fishersci.com](http://fishersci.com) and [fishersci.ca](http://fishersci.ca)

 **fisher scientific**  
part of Thermo Fisher Scientific

# Critical Cleaning in Clinical and Research Labs

**C**leaning will vary for hoods, work stations, incubators or anterooms, yet in each case it is important to consider fabric options, hardware durability, and the right tools for the application. At Micronova, we understand the necessity of meeting strict laboratory requirements with the highest level of control.



## BSC Hood and Work Center Cleaning

*Safe and ergonomic cleaning*

Cleaning in BSC, powder containment hoods, RABS, Workstations and other hard to reach areas present a variety of challenges. Micronova's specialty cleaning tools and covers provide options to increase effectiveness, reduce time, and enhance ergonomics.



**SilSat™ Wipes** clean and lubricate gaskets and seals in isolators, incubators, refrigerators and safety cabinets, with a blend of silicone and IPA.



**PolyMitt™** Cleaning Mitt Wipe conforms to the hand for ergonomic wipe downs of Shelves, Equipment, Carts, Wire racks, Door knobs, Chairs.



**Hand Held Mop** is small enough for table top equipment. Use with mesh or scrub covers for stubborn residues or simple sanitizing.



**NovaSponges™** with braided polyester mesh meet the lab equipment and parts stubbing needs without shedding.

**Isolator Tool™** Hoods, RABS, and Workstations are safely and easily cleaned. Locking 360° Swivel joint allows for interior and sashing cleaning as well.



**Double CornerCleaner Mop** is ideal for sanitizing awkward spaces and deep corners of Pass-throughs, BSC hoods, chambers, and small spaces.



**MegaSwat™** flexible blade mop for reaching under, around and inside Freezers, Centrifuges, Balances, Incubators.



Durable small diameter **Brushes** are superior for cleaning tubing, glassware and small parts.



# Critical Lab Cleaning

Selecting a low particulating mop cover and cleaning tools that support proper technique is the basis for an efficient cleaning protocol. Multi surface flat head mops simplify cleaning and are an economical, versatile, choice for specialty room cleaning. A modular, compact, and lightweight cart and wringer system is important from an ergonomic standpoint for technicians, but is also vital if storage is limited.



## Modular Cleaning



### NovaSnap™

The snap on and release function for this polyester mop head produces a simple lightweight system.

Mop frame  
QPSL-14 19035241



### PocketMop™

This flat head mop cover envelopes the mop frame and is ideal for cleaning and sanitizing surfaces and equipment.

NovaPoly, Smooth polyester mop cover 14"  
PKM6-14 19045241

Mop frame  
QPSL-14 19035241



### SlideMop™

This quilted polyester pad mop cover effortlessly falls into place with a downward push on the collapsible mop frame.

Polyester pad mop cover 16"  
MZCR63-16 19116757

Stainless Steel mop frame  
QPSM-16 19116750

*Mop covers available gamma irradiated*

**MegaWringer™**  
Downpress wringer  
W-94  
1916750

**Bucket Liners™**  
BKL7-3119IR  
191304833



**BucketBinder™**  
Double bucket system  
BB-DBL-DBYE 19263630

## Disposal and Removal for Waste

*FDA compliant LDPE liners and bags*

### BinLiners™

Liners are available in various configurations and are used for:

- transfer raw material
- product storage
- garment collection
- material disposal



*Pictured to the right is a stainless steel binstand used with liners*

### Biohazard Bags

Gamma Irradiated Biohazard bags are ideal for disposal of:

- pipettes
- swabs
- rodac plate test samples
- medical and infectious waste
- hazardous drugs




# Hand Hygiene

## Touch Free Sanitation

The act of removing viruses/bacteria or other unwanted substance from your hands is simple with Micronova's soaps and sanitizers

**Touch Free MicroDispenser™**  
ICDISP-3 19035380



**AquaHol™**  
Sterile Alcohol Sanitizer  
IC-220  
19042094

**NovaClenz®**  
Ethanol Gel Hand Sanitizer  
IC-210  
19035378C

**HyGClenz™**  
Antimicrobial hand soap with persistent activity  
IC-420  
19035378D

**Dispenser and Floor Stand**  
ICSTAND-02D3

## Good Documentation

Markers with IPA resistance or pens with indelible ink are among the options for good documentation.



## Marking & Sealing

CSLB-2WHIR  
19306985



Low and medium tack tapes, irradiated tapes and autoclave indicator tapes on a plastic core, are designed for cleanroom applications.




**Find product numbers and more details in Micronova's Specialized Mini Catalogs (Click or tap to view)**

### Cleanroom Pens and Markers



### Cleanroom Tapes



### Mopping Systems



### Liners for Critical Containment



### Buckets, Carts, and Wringers



Contact us today:

#### In the United States

Order online: fishersci.com  
Fax an order: 1-800-926-1166  
Call customer service: 1-800-766-7000

#### In Canada

Order online: fishersci.ca  
Fax an order: 1-800-463-2996  
Call customer service: 1-800-234-7437

