



DWK LIFE SCIENCES  
DURAN WHEATON KIMBLE

# MULTI-EFFICIENT

DURAN®  
MULTI-NECK  
HPLC BOTTLES



Find it at [fishersci.com](http://fishersci.com) and [fishersci.ca](http://fishersci.ca)

 **fisher scientific**  
part of Thermo Fisher Scientific

# DURAN® MULTI-NECK HPLC BOTTLES

DURAN® bottles are ideal for the storage and delivery of high-performance liquid chromatography (HPLC) mobile phases and solvents. Increasing the number of GL 45 inlets on the bottles adds convenience.

One additional side neck allows extra solvent to be added without disturbing the main delivery cap. The five GL 45 necks on the purging bottle increase safety and reduce the risk of cross contamination and air intake.

## FEATURES AND BENEFITS

- **Convenient** – Side neck allows solvent filling without main delivery cap removal. Purge bottle accepts up to 5 HPLC caps
- **Safe** – Capped bottles reduce the risk of volatile organic compounds (VOC) escape and spillage
- **Compatible** – GL 45 threads fit standard HPLC connection caps
- Manufactured from the material of choice for HPLC solvent bottles: **borosilicate 3.3 glass** offering a very low inorganic ion extractables profile
- DURAN® is a borosilicate glass type 3.3 as specified by the international standard ISO 3585, and is classified as a "**Type I, Class A low-expansion borosilicate glass**" as defined by ASTM-E438-1992
- With Retrace Code, and **certificate for manufacturing lot traceability**
- Durable white enamel volume graduations

## APPLICATION AREAS

- Mobile phase reservoirs for HPLC
- HPLC mobile phase purging
- Reagent supply to liquid handling systems and analysers



### WHAT IS HPLC "PURGING"?

When setting up an HPLC system, the aim of the purge is simply to flush through all the lines so that any remaining solvent in them from a previous analysis or wash is replaced with the new mobile phase.

## ORDERING INFORMATION

Description	d (OD) (mm)	h (mm)	Pack Unit	Cat. No
DURAN® GL 45 Bottle 2,000 mL with GL 45 Side Neck	136	265.5	1	09841316
DURAN® GL 45 Hydra Square Bottle, Clear, 1,000 mL, with 4 GL 45 Side Necks	94	222.5	1	09841317

## RECOMMENDED ACCESSORIES

Description	Pack Unit	Cat. No
KIMBLE® ULTRA-WARE® Three Hole Delivery Cap, ¼"-28 threads top and bottom with 1/8" ID holes to accept 1/8" OD PTFE tubing Supplied complete with three ¼"-28 CTFE plugs and a TFE/propylene O-Ring	1	K9539130000



Distributed by Fisher Scientific. Contact us today:

### In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

### In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437

 **fisher scientific**  
part of **Thermo Fisher Scientific**