



Numerous configurations and sizes of mobile cabinets are available in powder coated steel, steel with wood fronts, wood and stainless steel. Heavy-duty locking casters make units easy to maneuver in tight lab spaces. Reconfigure the storage space as your tasks change.



Optima™ Mobile Workstation

Optima™ workstations offer a push button height adjustable work surface to meet the height demands of individual users and equipment. Add plumbing, shelving and accessories to meet your specific requirements. This workstation can be relocated and height adjusted with ease; simply move it...



Heavy Duty Mobile Cart

Carts can be easily integrated into existing work spaces or located as a self-supporting workstation. High and low cart configurations available in various sizes. Work surfaces and shelves can be adjusted in height to meet the specific requirements of individual users or equipment. Suspended lower cabinets can be added for your storage needs. Built to withstand 2600lb loads.



Custom Heavy Duty Workstations

Mott can design innovative custom heavy duty workstations to meet your specific load requirements. Push button adjustable height work surface and an extensive variety of other features can be added to meet your demanding applications.



QUALITY IN PRODUCT • PRECISION IN PROCESS • STRENGTH IN PEOPLE

MOBILE FURNITURE SOLUTIONS



Corporate Headquarters
452 Hardy Rd, Brantford ON, Canada N3T 5L8
T. 519.752.7825 • F. 519.752.2895

Wood Casework Division - 562 Industrial Park Rd,
P.O. Box 156, Maxwelton, WV, USA 24957
T. 304.497.2115

www.mott.ca

QUALITY BY DESIGN

Mott Manufacturing provides complete solutions for all your laboratory needs. We offer a comprehensive line of products featuring quality steel laboratory furniture, premium grade architectural wood casework, adaptable and mobile furniture systems, high performance and custom fume hoods, custom stainless steel work surfaces and sinks. We have been serving the industrial, pharmaceutical, education, health care and government markets since the early 1960's. Each laboratory's needs are unique, we can design and engineer products to meet your exact requirements.

Consistent Product Quality

Our ISO 9001 Registered Quality Management System, supported by quality materials and skilled employees, ensures only high quality product is shipped.

Quality Product On-Time

Competitive manufacturing lead times, along with our high on-time and complete shipment performance, provides our customers with reliable delivery of high quality products.

Adaptable Product Line

Extensive modular product line permits our customers to readily configure the most productive solution to meet their specific requirements. Our flexible manufacturing capabilities allow convenient modification of standard components to fit the unique requirements of each laboratory.

Competitive Value

Decades of experience allows us to provide productive furniture solutions to our customers. This experience, along with our high product quality, custom manufacturing capability, consistently strong service levels, and competitive pricing provide high value to our customers.

Extensive Dealer Network

Our dealers are the best. We build strong relationships with our dealers to provide, in partnership, industry-leading service levels. Our dealers' on-site leadership and years of experience combined with our design and manufacturing skills deliver the very best in laboratory furniture solutions.



JUNE 2015

MOBILE SOLUTIONS



ISO 9001:2008 REGISTERED

The Solution to Change...

Today's laboratories demand efficiency; mobile furniture solutions support ever changing work processes. Systems must be adaptable and flexible, and provide the best possible solution for every situation, to meet your changing requirements.

We're committed to understanding your needs and developing solutions that better support your laboratory. Our mobile solutions have been designed to be compatible with a broad spectrum of Mott solutions.

Our commitment to quality through our ISO Quality Management system assures the reliable, long life of our laboratory furniture.

Everybody's needs are unique. Mott provides solutions to meet your specific application demands; simply contact us with your needs.



Optima™ 2550 Series can be used as a single-sided or double-sided bench system. This system is 84" high and has 6" deep "D" shaped rear posts and round front legs. Work surfaces and shelving can be adjusted in 1" increments to accommodate the requirements of individual users. Integrated plumbing, electrical and data allows for quick and easy access to services.



Optima™ 2100 Series can be used as a single-sided system or benches can be placed back to back to create a double sided work bench. This system is 84" high and has 2" round rear posts and front legs. Work surfaces and shelving can be adjusted in 1" increments to accommodate the requirements of individual users. Integrated plumbing, electrical and data allows for quick and easy access to services.



Optima™ Mobile Solutions - Heavy-duty swivel casters have leveling capabilities for times when you need both mobility and absolute stability.



Mott offers a variety of cost-effective tables to fit the unique needs of every type of lab. Tables are available in powder coated steel, stainless steel or wood. We offer a wide selection of tables, fully welded, adjustable, mobile, hydraulic, vibration isolation and more...



Flexible Lab Spaces - Multi-functional areas need to be user-friendly and easily reconfigured for different users, equipment or changing process. Tables can be put to use to increase productivity and ultimately, saving money. They are perfect for sharing instruments between different work areas or labs. Tables can be equipped with casters or leveling feet. Tables can accommodate our full range of suspended cabinets, shelves and drawers.



Altus™ Table Systems are reliable, attractive and affordable. Flexible for changing laboratories; these tables easily integrate into your into your existing lab or use them as a basis of a freestanding easily modified new lab. The table frame is all-welded construction and 84" high. Work surfaces and shelving are adjustable in 1" increments to accommodate the requirements of individual users. Integrated plumbing, electrical and data allows for quick and easy access to services.



Altus™ Mobile Solutions - Heavy-duty swivel casters have leveling capabilities for times when you need both mobility and absolute stability.



Tailored Solutions - We're committed to understanding your needs and developing solutions that better support your laboratory. Mott has designed and manufactured thousands of tables to fit the unique requirements of every type of lab; simply contact us with your needs.



Sigma Carts™ are highly adaptable and can be conveniently integrated into existing work spaces or located as self-supporting work stations, making them perfect for sharing instrumentation, tools, or processes. Ergonomic work surfaces and shelves are height adjustable in 1" increments to easily adapt to new tasks, equipment and users.

