

QUALITY IN PRODUCT • PRECISION IN PROCESS • STRENGTH IN PEOPLE

STAINLESS STEEL

STAINLESS STEEL



QUALITY BY DESIGN

Mott Manufacturing provides complete solutions for all your laboratory needs. We offer a full line of products featuring quality steel laboratory furniture, premium grade architectural wood casework, flexible and mobile furniture systems, high performance and custom fume hoods, custom stainless steel work surfaces and sinks. We have been serving the industrial, pharmaceutical, education, health care and government markets since the early 1960's. Since each laboratory's needs are unique, we can design and engineer products to meet your exact requirements.

Customer Satisfaction is Mott's First Priority

Delivering industry-leading customer service is more than just words to the Mott Team. We know that satisfied customers drive our success. We therefore provide our customers with the following key values:

Consistent Product Quality

Our ISO 9001 Registered Quality Management System, supported by quality materials and skilled employees, ensures only high quality product is shipped.

Just-in-Time Delivery Capability

Our six week manufacturing lead time, along with our high on-time and complete shipment performance, provide exceptional reliability to our customers.

Adaptable Product Line

Our extensive modular product line permits our customers to readily configure the most productive solution to meet their specific requirements. Our flexible manufacturing capabilities allow convenient modification of standard components to fit the unique requirements of each laboratory.

Competitive Value

Our several decades of application experience allows us to provide productive furniture solutions to our customers. This experience, along with our high product quality, custom manufacturing capability, consistently strong service levels, and competitive pricing provide high value to our customers.

Extensive Dealer Network

Our dealers are the best. We build strong relationships with our dealers to provide, in partnership, industry-leading service levels. Our dealers' on site leadership and years of experience combined with our design and manufacturing skills deliver the very best in laboratory furniture solutions.



Corporate Headquarters
452 Hardy Rd, Brantford ON, Canada N3T 5L8
T. 519.752.7825 • F. 519.752.2895

Wood Casework Division - 562 Industrial Park Rd,
P.O. Box 156, Maxwelton, WV, USA 24957
T. 304.497.2115

www.mott.ca



AUGUST 2011

mott
MANUFACTURING

ISO 9001:2008 REGISTERED

Quality Assured...

Mott's extensive stainless steel product offering is designed and manufactured to deliver practical solutions to specialized laboratory environments. Mott's stainless steel furniture and fume hoods are highly durable, easy to clean, simple to maintain and have an extended product life providing an exceptional value.

Mott's entire standard product line is available in stainless steel or can be custom designed for specific applications. The manufacturing process begins with the selection of quality ASTM A 240 Type 304 or 316 stainless steel and is welded using the Tungsten Inert Gas (TIG) method. The result is truly attractive laboratory furniture that will justifiably complement any laboratory.

Numerous North American laboratories feature Mott's stainless steel furniture. Applications include clean rooms, operating rooms, animal facilities, wash up areas, pharmaceutical and food manufacturing environments that feature high use of water wash-down to eliminate contamination and maintain excellent cleanliness.

Mott's commitment to quality through our ISO Quality Management system assures the reliable, long life of our entire product line of stainless steel furniture and fume hoods.



- **Moisture & Corrosion Resistance** - Properly selected and maintained stainless steel provides excellent moisture and corrosion resistance in several applications.
- **Impact Resistance** - Stainless steel's high yield strength, even at elevated temperatures, provides high impact strength and therefore ensures durability and extended product life.
- **Easy to Sanitize** - Stainless steel work surfaces with smooth radius corners make it ideal for environments requiring high levels of cleanliness. This includes applications such as pharmaceutical manufacturing, food processing, health care, agricultural products and services, clean rooms, operating rooms, animal labs, autopsy areas, and radioactive processes.
- **Attractive & Aesthetic Appearance** - Polished stainless steel provides bright, easily maintainable surfaces resulting in contemporary furniture design.
- **Fume Hoods** - All Mott's standard and premium fume hood designs are readily available in stainless steel construction. Hoods can be custom designed and manufactured to meet specific application requirements. Features can include cove corner liners with integral work surfaces where high levels of cleanliness are required.
- **Friendly to the Environment** - Stainless steel products are 100% recyclable and have value even after a very long life in the laboratory. 60% - 85% of new stainless steel comes from old remelted stainless steel scrap; which contributes to LEED® credit MR4.



Mott's fume hoods meet UL 1805 requirements and are classified by Underwriters Laboratories Inc. as to fire, electrical, and mechanical hazards only. Control number 685N.



- **Full Range of Cabinetry & Components** - All furniture components are die formed, notched, welded, and tightly fitted to provide maximum strength and durability. Available in type 304 or 316 stainless steel manufactured to a #4 satin finish.



- **Specialty Cabinets** - Customizing products is regular practice at Mott, not an add-on process that can significantly increase your costs. Types 304 and 316 stainless steel can be used in custom cabinets, carts and tables designed to meet your exact needs.



- **Work Surfaces & Sinks** - Built to order, all surfaces are welded and finish polished with quality craftsmanship to produce one piece integral construction.



- **Canopy Hoods** - Designed to exhaust non-toxic steam, odors and heat. Available wall-mounted or suspended from the ceiling and several features can be added to meet specific requirements. Powder coated steel also available.



- **Scullery Sinks & Scrub Sinks** - These units are custom designed and manufactured to your exact requirements and can be configured with several features such as single or multiple compartments, drain boards and front panels.