

HUGHES



HUGHES-SAFETY-SHOWERS.COM

Justrite Safety Group



their safety is essential.

Products made to work hard—and built to last.

16-Gallon Gravity-Fed Eyewash Station

Don't risk safety, even in remote locations

This self-contained, standalone eyewash station is suitable wherever continuous water isn't available—no plumbing needed. The gravity-fed eyewash supplies a continuous, controlled flow that exceeds the ANSI Z358.1 standard of 1.5 LPM (0.4 GPM) for 15 minutes, providing a reliable safety solution no matter where you are located.

- Portable design with carrying handles and concealed mounting bracket make moving and relocating easy
- Can be wall mounted or used on a bench or tabletop
- High-visibility color is easy to identify in emergency situations
- Simple to use with pull-down activation handle
- Easy to maintain—wide fill opening enables quick inspection, cleaning and filling
- Meets ANSI Z358.1 and EN15154 requirements

Closed position allows easy access to pull-down activation handle.



The eyewash delivers an even spray to flush eyes completely.

Description	Capacity US Gallons / Liters	Operating Duration (min)	H x W x D	Water Flow per Minute US Gallons / Liters	Product Weight Full/Empty	Approv/Lstg Regulation*	Model
Gravity-fed eyewash station (one single unit)	16 / 60.5	15	25.12 x 21.26 x 11.3 in (638 x 540 x 287 mm)	0.55 / 2	141.7 lbs (64.3 kg) 19.8 lbs (9 kg)	AN, ENS	19 111 121
Gravity-fed eyewash station (one pallet of 12 units)			57 x 44 x 37 in (1449 x 1118 x 940 mm)		Empty weight: 283.52 lbs (128.6 kg)		19 111 122

*AN = Meets or exceed American National Standards Institute (ANSI) Z358.1; ENS = Complies with European standard for showers, EN15154

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437

 **fisher scientific**
part of Thermo Fisher Scientific

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

© 2020 Justrite Safety Group · Code HNA210r5-Fisher (0620)