



Safety

Fisher Science Education

2022 Product Catalog

Featured Suppliers



CORNING



thermo scientific



fisherbrand



Visit fisheredu.com/featuredsuppliers to learn more about these suppliers and their products.

Helpful Icons



New product



Must be shipped by truck for regulatory reasons



Meets Americans with Disabilities Act requirements



Protects against splashes from hazardous chemicals or potentially infectious materials



Applicable for remote learning

Guarantee

If you're not 100% satisfied with your purchase, contact our customer service team within 30 days of your invoice date and we'll either exchange, repair, or replace the product, or give you a credit for the full purchase price. Call us toll-free for a return authorization number.

Special order items, furniture, and closeouts cannot be exchanged or credited.

Phone: 1-800-955-1177 • 7 a.m. to 5:30 p.m.
Central Time, Monday through Friday

Fax: 1-800-955-0740 • 24 hours a day, 7 days a week

Email: fse.custserv@thermofisher.com

Website: fisheredu.com

Address: Fisher Science Education
4500 Turnberry Drive
Hanover Park, IL 60133

For international orders, see page 40.

Products featured in this catalog may have limited availability and extended lead times. Prices are subject to change.

Connect with Us on Social Media



fisheredu.com/facebook



twitter.com/fishersciedu



pinterest.com/fishersciedu

Safety

Preparing today's students to be the innovators of tomorrow isn't always easy, but finding the right teaching tools can be. From basic lab supplies to state-of-the-art classroom technology, the Fisher Science Education team has everything you need to create a 21st century STEM learning environment. Visit fisheredu.com to get started.

Want to customize aspects of your curriculum?

Explore custom kits to meet the unique demands of your classroom. Visit fisheredu.com/customkitting to learn more.

Concerned about shipping and handling fees?

We offer \$12 flat-rate shipping on most orders and our hazardous material fee is just \$15 — the lowest in the country.

Looking for Something Else?

Check out our other 2022 catalogs:

- Chemicals
- Equipment and furniture
- Lab supplies and consumables
- Life sciences
- Physical science
- STEM



Before planting, farmers plow fields to bring fresh nutrients to the surface and bury weeds and remains of old crops.

Visit fisheredu.com/headlinediscoveries for more articles and information for your classroom.

Fire, Emergency, and First Aid	4–9	Aprons and Lab Coats	22–24
Fire	4	Facilities	25–34
Emergency and First Aid	5–9	Social Distancing Aids	25–26
Personal Protective Equipment	13–24	Cleaners and Sanitizers	27–28
Face and Eye Protection	13–17	Sinks and Eyewash Stations	29
Waste Boxes	13, 18	Safety Cans	30–31
Hand Protection	19–21	Storage Cabinets	32–34

Visit fisheredu.com to view our full offering of products.



ADVANCING MEDICAL EDUCATION

The SAM Basic™ and Nikki the Nursing Manikin with **Auscultation** include an extensive sounds library to prepare students for a wide range of real-life healthcare situations.

For more information visit fisheredu.com



NEW
Featured
Products



Basic Lab Protection

Keep staff and students safe while in the lab with these goggles and aprons.

Description	Cat. No.	Price
United Scientific™ Safety Goggles	S111036	\$4.25
United Scientific Plastic Disposable Aprons	S111125	\$7.55



Recycle Used Gear

Easily recycle used gloves and other safety products with the TerraCycle™ Zero Waste Box™ system.

Description	Cat. No.	Price
TerraCycle Disposable Gloves Zero Waste Box System, 39.66 L	S88854124	\$135.00
TerraCycle Disposable Gloves Zero Waste Box System, 79.31 L	S88854125	\$255.00
TerraCycle Disposable Gloves Zero Waste Box System, 136.42 L	S88854126	\$410.00
TerraCycle Safety Equipment & Protective Gear Zero Waste Box System, 79.31 L	S475653	\$300.00
TerraCycle Safety Equipment & Protective Gear Zero Waste Box System, 136.42 L	S475654	\$492.00



Disinfect Tools and Other Objects

Choose between two sizes, both with removable shelves and 30-minute timers.

Description	Cat. No.	Price
Supertek Scientific UV Sanitizer Cabinet, L x W x H: 26 x 10 x 29 in.	S81791	\$801.00
Supertek Scientific UV Sanitizer Cabinet, L x W x H: 20.5 x 16 x 16 in.	S81792	\$420.00



Keep Your Classroom in Working Order

Find a variety of chemical-specific Eisco™ wash bottles for your classroom.

Description	Cat. No.	Price
Eisco Performance Plastic Chemical Wash Bottle for Soap, 500 mL, Each	S15096	\$5.65
Eisco Performance Plastic Chemical Wash Bottle for Bleach, 500 mL, Each	S27220	\$4.70
Eisco Performance Plastic Chemical Wash Bottle for Bleach, 500 mL, 2/Pk	S27221	\$9.35
Eisco Performance Plastic Chemical Wash Bottle for Iso-propanol, 500 mL, 2/Pk	S27188	\$13.00
Eisco Performance Plastic Chemical Wash Bottle for Iso-propanol, 500 mL, 6/Pk	S27189	\$38.25



Stay Informed

These double-sided, laminated charts provide coronavirus safety information in both English and Spanish.

Description	Cat. No.	Price
NewPath Learning CDC Coronavirus Safety Signs in English & Spanish	S81705	\$31.25

Fire, Emergency & First Aid

Fire

Fire Blankets

Provides protection against fire.



- Material: 80% wool
- Sold individually

Cat. No.	Description	Price
S47567A	Fire Blanket with Pouch	88.00/Ea
S47559C	Fire Blanket with Wall Case	131.00/Ea

Fiberglass Fire Blanket in Packet

Lightweight and stored in a highly visible packet.

- Smother fires, wrap a burn victim or shield against encroaching flames
- Stored in a highly visible, red vinyl packet that opens easily
- Woven glass fiber material is totally safe to handlers
- Stitched with flame-proof thread



Cat. No.	Length	Width	Price
S47560A	4 ft.	4 ft.	93.00/Ea
S47560B	6 ft.	4 ft.	100.00/Ea
S47560C	6 ft.	6 ft.	127.00/Ea

Bulk First Aid Kit

Designed for large facilities.

- Contains economical bulk supplies that can be readily accessed when needed
- Latex-free
- Weatherproof
- Supplies housed in waterproof metal or plastic case
- Sufficient for 25 Persons



ORDERING INFORMATION:
Rapid refill kits available separately. For complete lists of kit contents, visit www.fisheredu.com.

COMPLIANCE:
ANSI™ Z308.1-2003 compliant.

Cat. No.	Description	Price
S98414B	Bulk First Aid Kit	78.00/Ea

EVAC+CHAIR™ Evacuation Solutions

EVAC+CHAIR delivers universal evacuation solutions for emergencies such as fires, floods, earthquakes, power outages and more.

Evacuation Stair Chairs (S72764 and S72765)

- For transport up or down stairways.
- Foldable for easy storage
- Lightweight construction
- Heavy-duty wheels
- Lockable rear wheels
- Wall mounting bracket
- Adjustable patient head restraint
- Friction tracks
- Ergonomic adjustable handle
- Quick-release safety belts (Head, chest, thigh)
- Reflective decals
- Dual-position seating
- Vinyl storage cover
- Instructional DVD and user guide
- Photo-luminescent Indicator Sign
- Lifetime warranty
- **S72765** also features a fold away grip handle/footrest



Rescue Sheet (S72774)

Located permanently and unobtrusively under patient's mattress; does not interfere with bed making. Patient is rescued with their own mattress.

- Long pulling straps
- Strong, easy-slide material
- Rubber side panels
- Reinforced cross straps with press stud
- Installation label
- Plastic strap buckles with press studs

Cat. No.	Description	Weight Capacity	Price
S72764	Model 300H MK4 Evacuation Stair Chair	400 lb.	1495.00/Ea
S72765	Model 500H Two Person Operation Enhanced Evac+Chair	500 lb.	2195.00/Ea
S72774	Rescue Sheet	660 lb.	249.00/Ea

United Scientific™ Ultraviolet Light

Ideal for forensics or illuminating fluorescent and blacklight reactive items. United Scientific™ Ultraviolet Light is handheld.

- Includes built-in LED flashlight
- Requires 4 AA batteries (not included)
- 6 in. 4w blacklight



Cat. No.	Description	Price
S23993	Ultraviolet light	9.70/Ea

Emergency & First Aid

DR-Scan™ Thermometer

DR-Scan™ provides an easy and safe way to know your body's temperature



- DR-Scan™ FeverScan is an easy and safe way to know your body's temperature
- 95° to 104°F
- Hold the DR-Scan firmly by both ends and press against dry forehead for 15 seconds
- Read while it's on the forehead, green indicates correct temperature
- For best results, use in a room that is around 72°F
- Use it 30 minutes after eating, drinking, coming in from outdoors or exercising
- Thermometer should be cleaned frequently with a cloth
- Store in plastic sleeve out of heat or sunlight
- Always consult your physician when a fever is detected
- Pack of 10

Cat. No.	Description	Price
S81424	DR-Scan™ Thermometer	5.00/Pk

Disposable Alcohol Swabs

Easy to open and ready to use.

Be sure you always have a good supply of double-thick alcohol pads on hand for blood collecting. Impregnated with 13mL of isopropyl alcohol, these 1.2 x 1.6 in. (3 x 4cm) pads come individually sealed in foil envelopes.



- Pack of 100

Cat. No.	Description	Price
S17032	Disposable Alcohol Swabs	9.40/Pk

Adhesive Bandages

Easy handling and application.

- Easy to open, even with gloved hands
- Just pull firmly on the ends to release the outer wrap
- Apply the bandage directly without touching the pad or the wound
- Each bandage individually wrapped and sterile



Cat. No.	Description	Dimensions (L x W)	Price
S47527A	Sheer Bandage	0.75 x 3 in.	10.50/Pk
S47527B	Sheer Bandage	1 x 3 in.	11.75/Pk
S47527C	Flex Bandage	0.75 x 3 in.	16.50/Pk
S47527D	Flex Bandage	1 x 3 in.	15.75/Pk

fisherbrand Latex-Free Cohesive Flexible Bandages

Flexible Bandages free of natural rubber proteins and allergens



- Maintain consistent compression without slipping or restricting mobility
- High porosity, water-resistant
- Lightweight, porous material allows skin to breathe for ultimate comfort
- Adhere to themselves — no clips or fasteners needed
- Assorted colors
- Length: 5 yd. (4.6m)
- Case of 30

Cat. No.	Description	Price
S25053	Latex-Free Cohesive Flexible Bandages	65.00/Cs

Stainless Steel Bandage Scissors

Bandage Scissors are for sizing bandages and dressings, and cutting through bandages and medical gauze



Cat. No.	Description	Price
S72158	Rainbow Color	8.40/Ea
S72161	Black Oxide coating	8.40/Ea
S72160	One big Ring	8.25/Ea

Mini-Bandage Scissors

Mini-Bandage Scissors are 3.5in long with sharp blades



Cat. No.	Description	Price
S72159	Mini-Bandage Scissors	2.65/Ea

MobileAid™ Pro Easy-Roll™ High-Visibility Trauma First Aid Stations

Rolling high-visibility First Aid Bags with flag, stretcher, and 5 clear-view organizer boxes. Designed to allow EMT personnel to equip them with their own custom supplies.

- Both models 100 and 200 include an integrated gear and seating platform, an 8-foot extendable first aid flag, and flexible evac stretcher
- Model 100 has 5 clear-view organizer boxes utility pouches with vinyl lids for easy supply identification
- Model 200 has 5 Quick-Access utility pouches with customizable handle IDs and clear vinyl lids for easy supply identification
- Polyester with aluminum frame
- First aid consumables and other gear not included



Cat. No.	Description	Price
S29596	Model 100 Clear Organizer Boxes	336.00/Ea
S29597	Model 200 Organizer Pouches	403.00/Ea

MobileAid™ Modular Trauma First Aid Stations

Available in assorted configurations to meet your needs.

S29595: Hi-Visibility 260-piece mobile Modular Trauma First Aid Station

Includes:

- 5 Quick-Access first aid supply modules
- Crush-resistant frame
- Flexible stretcher
- First aid flag
- Flashing wheels
- High-visibility reflective safety stripe
- Integrated seat and gear platform.
- OSHA compliant: ANSI/ISEA Z308.1-2015 Class A, Type II first aid kit



S29686/S29687/S29688: First Aid Stations with Four Bleeding Control Kits

Extensive 260-piece mobile Modular Trauma First Aid Station with 4 additional Bleeding Control Kits to treat bleeding wounds from accidents, disasters, or active shooter situations.

- Compact dual-wheel mobile station Organized in quick-response first aid modules for rapid response and treatment
- Includes a flexible evacuation stretcher and First Aid Flag
- Integrated seat and gear platform allow patient to sit and be treated or can be used to haul additional equipment (e.g. AED)



S29685 and S29684: First Aid Stations with Bleeding Control Kits

Extensive 260-piece mobile Modular Trauma First Aid Station with an additional Bleeding Control Kit for greater ability to treat bleeding wounds from accidents, disasters, or active shooter situations. This compact dual-wheel mobile station is organized in quick-response first aid modules for rapid response and treatment.

Includes:

- MobileAid™ BleedStop Compact Bleeding Control Kit
- Tourniquet
- Mini-compress bandage
- QuickClot blood clotting gauze
- Petroleum gauze for chest wounds
- Other bleeding wound supplies
- Flexible evacuation stretcher
- First Aid Flag
- Integrated seat and gear platform



Cat. No.	Description	Price
S29595	Modular Trauma First Aid Station	519.00/Ea
S29684	Modular First Aid Station with BleedStop Compact 100 Bleeding Control Kit (SWAT-T Tourniquet)	578.00/Ea
S29685	Modular First Aid Station with BleedStop Compact 100 Bleeding Control Kit (SOFT-W Tourniquet)	614.00/Ea
S29686	Modular First Aid Station with 4 BleedStop 100 Bleeding Control Kits (SWAT-T Tourniquets)	927.00/Ea
S29687	Modular First Aid Station with 4 BleedStop 200 Bleeding Control Kits (SOFT-W Tourniquets)	1056.00/Ea
S29688	Modular First Aid Station with 4 BleedStop 300 Bleeding Control Kits (SOFT-W Tourniquets & Chest Seals)	1190.00/Ea

MobileAid™ Easy-Roll™ Emergency Incident Command Station

Easy-Roll Emergency Bag.

- 4-person team
- Includes supplies for a command shelter as well as for individual team members
- Includes safety vest, N95 mask, whistle, two-way radio, pen writing tablet, leather-palm gloves, flashlight, batteries, light stick, poncho, emergency blanket, and caution tape
- Material: Polyester, Aluminum



Cat. No.	Description	Price
S29608	MobileAid™ Easy-Roll™ Emergency Incident Command Station	843.00/Ea

MobileAid™ Reflex 200 Hi-Visibility Trauma First Aid Backpack Kit

A 218-piece trauma first aid kit with supplies packed in 3 OTS (Over-The-Shoulder) Medical Supplies Clear-View Pouches with Organizer Trays and 2 small clear-view pouches.



- Includes a high-visibility safety vest
- Carried in a lightweight MobileAid™ Reflex backpack
- Material: Rip-Stop, Polyester

Cat. No.	Description	Price
S29610	MobileAid™ Reflex 200 Hi-Visibility Trauma First Aid Backpack Kit	408.00/Ea

MobileAid™ BleedStop™ Bleeding Control and Gunshot Wound Trauma First Aid Kits

Immediately treat bleeding wounds from accidents, disasters, or shootings with tourniquets, trauma dressings, chest seals, and other essential bleeding control supplies.

BleedStop COMPACT 100 Bleeding Control and Gunshot Wound Trauma First Aid Kit (S29600)

- Supplied in a SecureLock clear plastic bag
- Includes SWAT-T tourniquet, mini-compress bandage, hemostatic blood clotting gauze, petroleum gauze for chest wounds, and other bleeding wound supplies
- Treats 1 bleeding wound



BleedStop IMMEDIATE RESPONSE 200 Bleeding Control and Gunshot Wound Trauma First Aid Kit (S29601)

- Single Bleeding Wound Trauma Kit supplied in an clear-view ballistic nylon pouch
- Includes SOFT-W tourniquet, mini-compress bandage, hemostatic blood clotting gauze, petroleum gauze for chest wounds, and other bleeding wound supplies
- Treats 1 bleeding wound



BleedStop DOUBLE 300 Bleeding Control and Gunshot Wound Trauma First Aid Kit (S29602)

- BleedStop DOUBLE bleeding control trauma first aid kit includes TWO SOFT-W tourniquets, BloodStopper compress bandages, chest seals, and petroleum gauge pads for torso wound.
- Treats 2 bleeding wounds
- Allows response to multiple casualties simultaneously



Cat. No.	Description	Price
S29600	BleedStop Compact 100	74.50/Ea
S29601	BleedStop Immediate Response 200	127.00/Ea
S29602	BleedStop Double 300	290.00/Ea

MobileAid™ XL™ Trauma First Aid Backpack Kits

Both kits include a flexible stretcher and ID placards.

S29612: High-capacity First Aid backpack with supplies (260 pieces) in clear-view organizer boxes.

- Includes a flexible stretcher and ID placards

S29611: High-capacity First Aid backpack with Trauma First Aid supplies (260 pieces) in pouches with organizer trays.



Cat. No.	Description	Price
S29612	XL 100 Trauma First Aid Backpack Kit	506.00/Ea
S29611	XL 200 Trauma First Aid Backpack Kit	624.00/Ea

Evacuation and Shelter-In-Place Emergency Kits

S29587: Emergency kit for up to 10 people come equipped with essential emergency supplies that can last for the first several hours of an emergency.

- Designed for small offices and work groups
- Packed in a durable duffel bag for easy storage and contains a Carry-Flat toilet that fits in the duffel
- Includes emergency instructions, search and rescue supplies and tools, dust masks for easier breathing, first aid, water (5 year shelf life), and items that provide warmth, shelter, communication, light and sanitation



S29583: Classroom lockdown emergency kit with supplies for 25 students and 1 teacher for first several hours after an emergency or disaster.

- Comes in a sealed bucket that converts to a toilet
- Meets the immediate needs of up to 25 students and 1 teacher for the first several hours after an emergency or disaster
- Includes supplies for breathing protection, hydration, first aid, warmth and shelter, communication and light, hygiene and sanitation, and search and rescue

S29588: Emergency kit with supplies for up to 25 people that is designed to last for the first several hours after an emergency or disaster

- Storage bucket converts to an emergency toilet
- Includes emergency instructions, search and rescue supplies and tools, dust masks for easier breathing, first aid, water, and items that provide warmth, shelter, communication, light, and sanitation

Cat. No.	Description	Price
S29587	10-Person Office Evacuation and Shelter-in-Place Emergency Kits	207.00/Ea
S29583	25-Student Classroom Lockdown and Emergency Response Kit	206.00/Ea
S29588	25-Person Shelter-in-Place Emergency Kit	206.00/Ea

Extended Emergency Support Kits

S29589: 10-person Compact survival kit with supplies for 24 hours of emergency support for up to 10 people.

- Use as a companion kit to the 10-Personal Office Emergency Kit
- Includes water, food, light sticks, and emergency thermal blankets
- Carry box has a handle and is stackable



S29590: Compact, waterproof emergency kit for up to 25 people comes equipped with supplies for a full 24 hours of emergency shelter-in-place support.

- Includes water, water purification tablets, food, light sticks, and emergency thermal blankets

S29584: Compact classroom emergency kit extends the supplies available for up to 25 students and 1 teacher in an emergency that lasts beyond a few hours.

- Includes supplies for a full 24 hours: food, water, water purification tablets, light sticks, and emergency blankets.

Cat. No.	Description	Price
S29589	Emergency Support Kit for 10 People, 24 hrs.	138.00/Ea
S29590	Support Shelter-in-Place Kit for 25 People, 24 Hrs	196.00/Ea
S29584	Extended Support Emergency Kit for 25 Students/1 Teacher, 24 Hrs.	196.00/Ea

Complete 25-Person Shelter-In-Place Combo Emergency Kits

Emergency kits with preparedness supplies for 25 people is equipped with enough supplies to last from 24 hours to three days after an emergency or disaster.



- Includes emergency instructions, search and rescue supplies and tools, dust masks for easier breathing, first aid, water, water purification tablets, food, and items that provide warmth, shelter, communication, light, and sanitation
- Storage bucket converts to emergency toilet

Cat. No.	Description	Price
S29591	24-hour Shelter-in-Place Kit	418.00/Ea
S29592	2-Day 25 person Shelter-in-Place Kit	598.00/Ea
S29593	3-Day Shelter-in-Place Emergency Kit	794.00/Ea

SchoolGuard™ Hi-Visibility Teacher and Staff School Emergency Kits

High-Visibility teacher and staff Mini-Backpack emergency response kit with the addition of a Bleeding Control Kit for treating bleeding wounds from accidents, disasters, or active shooter situations.

Compact storage on shelves or in drawers or cabinets.

Includes:

- MobileAid™ BleedStop Bleeding Control Kit
- Tourniquet
- Mini-compress bandage
- QuickClot blood clotting gauze
- Petroleum gauze for chest wounds
- Other bleeding wound supplies
- Other school emergency supplies:
- Identification
- Communication
- Breathing protection
- Light



Cat. No.	Description	Includes	Price
S29689	BleedStop Compact 100 Bleeding Control Kit	SWAT-T Tourniquet	145.00/Ea
S29690	BleedStop Compact 200 Bleeding Control Kit	SOFT-W Tourniquet	181.00/Ea

SchoolGuard™ Hi-Visibility Classroom Evacuation and Lockdown Kits

This high-visibility school emergency kit for up to 29 students and 1 teacher meets all 10 essential survival needs for lockdown or evacuation and includes a Bleeding Control Kit for treating bleeding wounds from accidents, disasters, or active shooter situations.

Portable fold-flat toilet stores in the survival backpack with other survival supplies: Water, food, breathing protection, first aid, shelter, warmth, communication, light, tools, and sanitation/hygiene.

Includes:

- MobileAid™ BleedStop Bleeding Control Kit
- Tourniquet
- Mini-compress bandage
- QuickClot blood clotting gauze
- Petroleum gauze for chest wounds
- Other bleeding wound supplies



Cat. No.	Description	Includes	Price
S29691	BleedStop Compact 100 Bleeding Control Kit	SWAT-T Tourniquet	346.00/Ea
S29692	BleedStop Compact 200 Bleeding Control Kit	SOFT-W Tourniquet	382.00/Ea

SchoolGuard™ Easy-Roll Classroom Evacuation and Lockdown Kit

This Easy-Roll high-visibility school emergency kit for up to 29 students and 1 teacher meets all 10 essential survival needs for lockdown or evacuation and includes a Bleeding Control Kit for treating bleeding wounds from accidents, disasters, or active shooter situations.

Portable fold-flat toilet stores in the survival backpack with other survival supplies: water, food, breathing protection, first aid, shelter, warmth, communication, light, tools, and sanitation/hygiene

Includes:

- MobileAid™ BleedStop Compact Bleeding Control Kit
- Tourniquet
- Mini-compress bandage
- QuickClot blood clotting gauze
- Petroleum gauze for chest wounds
- Bleeding wound supplies



Cat. No.	Description	Includes	Price
S29693	BleedStop Compact 100 Bleeding Control Kit	SWAT-T Tourniquet	366.00/Ea
S29694	BleedStop Compact 200 Bleeding Control Kit	SOFT-W Tourniquet	402.00/Ea

How to Choose Respiratory Protection

Protecting What Matters Most

Keeping employees and customers safe is a top priority for all of us, and protecting them from airborne particles and pathogens has never been more important. Use this guide to help choose respiratory protection for you and your colleagues. Find related definitions and distinctions between types of masks and respirators and PPE options based on the operational task.



Cloth Masks

Cloth masks are reusable (washable) and nonsterile. Commonly worn to reduce aerosols that can spread viruses in social situations, they are soft, breathable, and comfortable for everyday wear. Most have front panel pleats, soft elastic ear loops, and an adjustable nose piece for better fit. Cotton masks offer low levels of filtration efficiency and are not recommended for surgical use. They do not require fit testing and are not airtight.

Eisco™ Cloth Face Masks, Adult or Child Size

Cat. Nos. S27420, S27421

Available in U.S. only



Disposable Non-Woven Masks

Disposable face masks are designed for single or limited use and are similar to surgical masks in filtration and fluid resistance. Universally sized, they are latex free with soft, elastic ear loops. Not intended for medical use, disposable face masks have 2- or 3-ply construction, do not require fit testing, and are not airtight.

Fisherbrand™ Disposable Face Masks with Ear Loops

Cat. No. 12-888-001A

Available in U.S. and Canada




Procedure and Surgical Masks

Surgical masks are tested to meet levels of barrier protection as specified in ASTM F2100-19. Procedure masks are FDA-approved medical devices for use in healthcare procedures. They do not require fit testing and are not airtight.

Medicom™ SafeMask™ FreeFlow Earloop Mask, Level 3

Cat. No. 19-910-717

Available in U.S. only



N95 Respirators

An N95 mask or respirator filters at least 95% of airborne particles as defined by the U.S. National Institute for Occupational Safety and Health (NIOSH, 42 CFR Part 84). They are functionally equivalent to FFP2 respirators (European Union) and KN95 respirators (China), with slightly different certification criteria.

In the United States, the Occupational Safety and Health Administration (OSHA) requires that healthcare workers use N95 respirators if they are working with patients suspected or confirmed to be infected with COVID-19, Ebola, Mycobacterium tuberculosis, pandemic or avian influenza, or severe acute respiratory syndrome (SARS).

Fit testing is required by OSHA for N95 respirator wearers to both identify the right model, style, and size, and to verify fit thereafter. These respirators can make breathing more difficult, so the wearer must also be medically evaluated for heart conditions, lung disease, and claustrophobia.

Moldex™ N95 Disposable Particulate Respirators

Cat. No. 18-999-1449

Available in U.S. and Canada

Surgical N95 Respirators

Surgical respirators (also called Surgical N95, medical, or healthcare respirators) are N95 masks used specifically in healthcare; they meet both FDA and NIOSH requirements. In response to shortages during the COVID-19 pandemic, the Families First Coronavirus Response Act changed liability and certification requirements to allow the use of industrial N95 respirators.

Fit testing is required by OSHA for N95 respirator wearers to both identify the right model, style, and size, and to verify fit thereafter. These respirators can make breathing more difficult, so the wearer must also be medically evaluated for heart conditions, lung disease, and claustrophobia.

Moldex™ N95 Healthcare Particulate Respirators Surgical Masks

Cat. No. 19-150-0429

Available in U.S. only

Available Resources

- <https://www.cdc.gov/coronavirus/2019-ncov/prevent-getting-sick/about-face-coverings.html>
- https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/default.html
- <https://www.cdc.gov/niosh/npptl/pdfs/UnderstandDifferenceInfographic-508.pdf>
- https://www.osha.gov/pls/oshaweb/owadisp.show_document?p_id=12716&p_table=standards
- <https://www.nelsonlabs.com/en14683-harmonizes-bacterial-filtration-efficiency-and-differential-pressure-with-astm-f2100/>
- <https://buyandsell.gc.ca/specifications-for-COVID-19-products#200>
- <https://www.canada.ca/en/health-canada/services/drugs-health-products/covid19-industry/medical-devices/personal-protective-equipment/medical-masks-respirators/face-covering-classifications-notice.html>

Useful Links

- <https://www.cdc.gov/niosh/npptl/hospresptoolkit/default.html>
- <https://www.cdc.gov/coronavirus/2019-ncov/communication/print-resources.html>

ASTM Levels

- **Level 1:** Low barrier protection; for low-risk, non-surgical procedures and exams that do not involve aerosols, sprays, or fluids
- **Level 2:** Moderate barrier protection; for low-to-moderate levels of aerosols, sprays, and fluids
- **Level 3:** Maximum barrier protection; for situations with potential for heavy aerosols, sprays, and fluid exposure

Due to recent unprecedented global demand, some respiratory protection products have limited availability and extended lead times.

Visit fishersci.com/respiratory-protection or fishersci.ca/respiratory-protection to find other respiratory PPE products.

Visit fishersci.com/safety-specialist-contact or fishersci.ca/safety-specialist-contact to contact your safety specialist.



Eye Protection

The Fisher Science Education brand offers a full line of goggles and spectacles that meet OSHA and ANSI standards for eye protection. The information provided below will help you choose the eye protection that is best suited to your application. For additional information and assistance in selecting eye protection, contact Fisher Science Education customer service at **1-800-955-1177**.

OSHA requires that eye and face protective devices meet the ANSI Z87.1 standard for Occupational and Educational Personal Eye and Face Protection Devices.

The ANSI Z87.1 standard sets forth requirements for the testing, marking, selection, care and use of eye protection devices, including spectacles, goggles, and face shields. Rigorous testing including high-mass and high-velocity impact tests are required to verify impact and penetration resistance. ANSI Z87.1 also emphasizes that protective eye devices do not provide complete protection and should not be used as substitutes for protective measures such as mechanical hoods and operational controls.

To comply with OSHA regulations, eye protection equipment must:

- Provide adequate protection against the hazards for which it is designed
- Be marked to identify the manufacturer and that it meets ANSI Z87.1
- Be reasonably comfortable under the conditions of use
- Be durable
- Fit securely without interfering with vision or movement
- Be kept clean and in good repair

Protective Spectacles

Use these glasses for impact protection against flying debris and fragments. Protective spectacles feature plastic, nylon, graphite, or metal frames and removable, vented, or solid sideshields for added protection from particles, such as sand and dirt. All protective spectacles have strong polycarbonate lenses that provide impact resistance and durability, as well as protection from ultraviolet (UV) radiation. Because pits or scratches in lenses may affect impact resistance, lenses should be inspected regularly and replaced if lenses are damaged.

Impact Goggles

These goggles provide more peripheral protection from flying particles than protective spectacles. Perforations in goggle frames provide airflow for comfort and help prevent fogging. Some lenses also feature an antifog coating. Many goggles are designed to fit over prescription eyewear and some will accommodate extra-large eyeglass frames.

Chemical Splash Goggles

Using indirect vents, these goggles protect against splashes of hazardous chemicals or potentially infectious materials. A labyrinth ventilation system rejects contaminants and provides excellent splash protection. Strong polycarbonate lenses provide durability, clarity of vision, and UV protection.

Face Shields

Face shields feature acetate, propionate, or polycarbonate windows for full facial protection against sparks, splashes, and splatter. Polycarbonate windows also provide UV protection, while some specially treated or coated windows offer heat reflection or protection from infrared radiation (IR). Faceshields are secondary protection and must be worn with protective spectacles or goggles.

Personal Protective Equipment

Face and Eye Protection

Eisco™ Disposable Face Mask

NEW

Disposable face masks for everyday use. Secured with elastic ear loops.

- Disposable paper face masks
- Elastic hooks around ears
- Pack of 50



Cat. No.	Description	Price
S27422	Adult: 170 x 95 mm	45.00/Pk
S27423	Child: 144 x 95 mm	45.75/Pk

3M™ N95 Particulate Respirator

Lightweight and comfortable with snug, secure seal.

At least 95% NIOSH tested filter efficiency for use where solids or non-oil-based liquid aerosols are present.

- Synthetic elastic headband
- Adjustable nose clip helps provide a custom fit and reduces the potential for eyewear fogging
- Soft nose foam
- Size: Universal
- Pack of 20



These safety products must be used in accordance with OSHA regulations and the user instructions, warnings and limitations accompanying each product.

COMPLIANCE:
Meets NIOSH 42CFR84 N95 requirements; NIOSH approval TC-84A-0007.

Cat. No.	Description	Price
S47062A	N95 Particulate Respirator	42.50/Pk

Eisco Safety Face Shield

NEW

Fits securely to deep molded brow guard with adjustable crown strap for protection against chemical splashes and impact from flying objects

- For use in any home, school, or professional laboratory
- Fits securely to deep molded brow guard with adjustable crown strap
- Protects against chemical splashes and impact from flying objects
- Steel bound plastic visor
- Non-flammable
- Dimensions (L x W): 7 x 8 in.



Cat. No.	Description	Price
S27083	Safety Face Shield	11.00/Ea

fisherbrand Disposable Face Shield

High-quality, optically clear, distortion-free 7mil film.

- Polyester lens is treated with antifog and antistatic coating
- Headpiece is held in place by a stretch strap and includes a latex-free foam band for comfortable extended wear
- Pack of 25

COMPLIANCE: Conform to OSHA 29CFR Part 1910.1030 for bloodborne pathogens and the CDC's recommendations for protective eyewear.



Cat. No.	Description	Price
S98127	Disposable Face Shield	50.00/Pk

TerraCycle™ Disposable Masks Zero Waste Box™

Recycle disposable masks. Each box includes signage, a built-in liner, shipping to and from your location, and processing of the materials.

- Select a Zero Waste Box for a specific waste stream and place your order.
- Collect the accepted items to fill the system.
- Send it back to TerraCycle with the prepaid return label on your system.
- Order your next Zero Waste Box system to continue to recycle.
- Includes: One Zero Waste Box system, One prepaid return label affixed to your box, Outbound shipping to your location, Processing and recycling of the materials, Customer support from TerraCycle's team of experts
- **Types of masks that can be included:** disposable face masks, including N95, KN95, 3-ply nonwoven, and dust masks; cannot include paper masks or reusable fabric masks.

⚠️ ALERTS: Not suitable for medical waste. This includes any 'red bag waste' (i.e. material contaminated with blood or bodily fluids that originate from health care facilities, hospitals, physicians' offices, dental offices, blood banks, veterinary clinics, research laboratories, etc.) Does not accept household items that may be considered hazardous waste. Do not ship dripping packages. This product is not available in Hawaii, Alaska, or Puerto Rico.



Cat. No.	Dimensions (L x W x H)	Capacity	Price
S88854115	11 x 11 x 20 in.	39.66 L	86.00/Ea
S88854116	11 x 11 x 40 in.	79.31 L	125.00/Ea
S88854117	15 x 15 x 37 in.	136.42 L	219.00/Ea

Astrospec 3000™ Eyewear

Quality, comfort and styling make these a popular choice.

- Adjustable temple length and lens inclination system
- Duoflex™ comfort cushions relieve pressure behind the ear
- Optically correct lens has Uvextreme™ AF anti-fog and anti-scratch coating



COMPLIANCE:

Meet ANSI™ Z87.1-2003 standards.

Cat. No.	Frame Color	Lens Type	Price
S476251	Patriot (Red, White, Blue)	Anti-fog	13.25/Ea
S476252	Black	Scratch-resistant, Anti-fog	13.00/Ea

fisherbrand 200 Series Safety Glasses

Economical, lightweight protection.

- One-piece lens provides a full and unobstructed field of vision
- Absorb >99.9% of harmful UV light
- Molded-in sideshields and browguards
- Adjustable spatula temples
- Pack of 12



COMPLIANCE: Meet ANSI™ Z87.1-2003 standards.

CERTIFICATIONS: CSA-approved.

Cat. No.	Description	Price
S01821	Blue Frame	16.50/Pk
S01822	Black Frame	16.50/Pk

fisherbrand Racer Series Qualifier Eyewear

Wraparound design offers superior side coverage.

- Scratch-resistant, clear polycarbonate lenses block >99.9% of harmful UV radiation
- Wraparound design offers superior side coverage
- Lightweight frames are durable and comfortable
- Flared single-lens design provides sleek protection
- A soft nose bridge and soft, bayonet-style side arms grip the face and temples lightly but firmly for a comfortable fit
- Choose from black or red nylon frames
- Pack of 12



Cat. No.	Frame Color	Lens Tint	Price
S98423	Red	Clear	83.00/Pk
S98420	Black	Clear	83.00/Pk

fisherbrand Visitorspec Safety Glasses

Ideal for protecting visitors in industrial or medical settings.

Economical safety spectacles available with coated or uncoated lens.

- Pack of 12



Uncoated Lens

- 100% impact-resistant polycarbonate wraparound unilens design with blue temples
- Clear lens absorbs >99.9% of harmful UV light

Coated Lens

- 100% impact-resistant lens absorbs >99.9% of harmful UV light
- Polycarbonate wraparound style with molded-in browguard and single lens
- Coated to resist scratching, fogging and static
- Adjustable-length spatula temples ensure a secure, comfortable fit

COMPLIANCE:

Meet ANSI™ Z87.1-2003 standards.

Cat. No.	Description	Price
S01819	Coated Lens	44.00/Pk
S01820	Uncoated Lens	28.00/Pk

fisherbrand Impact Goggles

Lens coated to minimize scratching and fogging.

- Built for comfortable protection with a soft, flexible PVC body that conforms to facial contours
- Large frames fit over most prescription eyewear
- Coated, clear polycarbonate lens resists impact, scratching and fogging
- Headstrap easily adjusts for a custom fit
- Protect eyes from dust, mist and other airborne particles
- perforations on goggle top and sides provide direct ventilation



COMPLIANCE:

Meet ANSI™ Z87.1-2003 and CSA standards.

Cat. No.	Description	Price
S01823	Impact Goggles	22.75/Ea

Magnifying Safety Glasses with Clip-on, Flip-up Lens System



Carson™ Magnifying Safety Glasses with Clip-on, Flip-Up Lens System include a pair of 1.5X power (+2.5 diopter) clip-on, flip-up magnifying lenses to help protect the eyes from all sorts of hobbies, crafts, tasks and other low vision needs.



- Distortion-free clip-on, flip up magnifying lenses contain a 1.5X power magnification which is equivalent to +2.5 diopters.
- Scratch-resistant polycarbonate lenses and cushioned rubber temple pads for maximum user comfort.
- Clear lens design provides UV A/B/C protection along with protective side shields.
- Rated for eye protection against low to moderate impact hazards but are not unbreakable and are not recommended for any high impact hazards such as but not limited to shooting, racket sports, paintball or any other similar activities.
- Do not provide unlimited protection.
- Must be worn for the whole duration of exposure to hazards. (Leave the work area if any irritation occurs or if the glasses become damaged or cracked.)
- Meet the ANSI/ISEA Z87.1-2015 personal eye protection standard for impact resistance.
- Not rated to protect from hazards including, but not limited to chemical, heat, welding, dust, light (IR/UV/Laser/etc.) and splash hazards.
- Additional protective eyewear such as face shields and helmets should be worn to provide additional protection depending on the application and the manufacturer's specifications.
- Come with a protective PVC storage pouch with a magnetic enclosure system.
- Ideal for multiple hobbies, crafts, tasks and other low vision tasks.
- Limited lifetime warranty

Cat. No.	Description	Price
S81167	Magnifying Safety Glasses with Clip-on, Flip-up Lens System	13.75/Ea

Eisco™ Clear Vented Safety Glasses

Clear Safety Glasses fit over most prescription lenses and have vents to prevent fogging

- Made of light weight polycarbonate and impact resistant lenses with UV protection
- Made of durable polycarbonate
- CE marked, conforms to ANSI Z87.1 2010 sideshield protection.



WARRANTY: 1 year limited manufacturer warranty

Cat. No.	Description	Price
S41140	Clear Vented Safety Glasses	5.00/Ea

Eisco™ Adjustable Clear Vented Safety Goggles

Adjustable Clear Safety Goggles are designed to protect eyes from chemical splashes and dust.

- Fit over most prescription lens
- Includes side vents to prevent lenses from fogging up
- Features black adjustable strap to fit most



WARRANTY: 1 year limited manufacturer warranty

Cat. No.	Description	Price
S41141	Adjustable Clear Vented Safety Goggles	8.05/Ea

fisherbrand Economy Chemical-Resistant Goggles

Affordable chemical-resistant goggles.

- Clear polycarbonate lens
- Soft, flexible PVC body molds to facial contours for a comfortable fit
- Adjustable headstrap ensures custom fit
- Resist a broad range of chemicals
- Circular vents allow air to flow freely through goggle body



COMPLIANCE:

Meets ANSI/ISEA Z87.1-2015 American National Standard for Occupational and Educational Personal Eye and Face Protection Devices, Section 8.1 Droplet and Splash Hazard (D3)

Cat. No.	Description	Price
S68660A	Chemical Splash Goggles	4.50/Ea

fisherbrand Chemical Splash Goggles

Lens coated to minimize scratching and fogging.

- Built for comfortable protection with a soft, flexible PVC body that conforms to facial contours
- Large frames fit over most prescription eyewear
- Coated, clear polycarbonate lens resists impact, scratching and fogging
- Headstrap easily adjusts for a custom fit
- Resist a wide range of chemicals and feature indirect ventilation to direct airflow over the lens



COMPLIANCE: Meet ANSI™ Z87.1-2003

Cat. No.	Description	Price
S01817	Chemical Splash Goggles	6.65/Ea

fisherbrand V-Maxx™ Chemical Splash and Impact Goggles

Comfortable fit ensures protection — even for contact wearers.

Offer protection against impact, dust, mist and chemical splash.



- FogBan™ lens coating
- Wraparound design provides maximum coverage and 180° field of vision
- Soft rubber construction ensures secure, comfortable fit
- Large enough to be worn over prescription glasses

COMPLIANCE:

Meet ANSI™ Z87.1-2003 and CSA™ Z94.3 standards.

Cat. No.	Description	Price
S68660	Chemical Splash and Impact Goggles	13.25/Ea

Eisco™ Green Safety Goggles

General-purpose goggles featuring side vents to prevent fogging.

- With adjustable elastic strap
- Fit over prescription lens
- Protect from chemical splash and dust



Cat. No.	Description	Price
S13975	Green Safety Goggles	1.95/Ea

Visorgogs™ Goggles

Designed to protect eyes from flying particles or liquids.

- Offer complete protection from glaring overhead light
- Lenses are one-piece, impact-resistant plastic with a down-angle vision

COMPLIANCE:

Meet ANSI™ Z87.1-2003 (Z87+ High Impact).



Cat. No.	Description	Price
S47619S	Full Vision Goggle	7.50/Ea

Eisco™ Adjustable Green Safety Goggles, Indirect Vent

Designed to protect eyes from chemical splash and dust.

- Fit over most prescription lenses
- Indirect vents prevent fogging
- Borosilicate glass with clear lens
- Includes adjustable elastic strap
- General purpose



WARRANTY AND SERVICES:

1 year limited manufacturer warranty

Cat. No.	Description	Price
S88323	Adjustable Green Safety Goggles	1.95/Ea

Standard Goggles

Protective goggles with perforated styling for maximum ventilation.

- Small
- Latex-free elastic strap
- Made of PVC for comfortable fit over prescription eyewear
- Polycarbonate lens filters 99.9% of harmful UV rays
- Duramass™ scratch-resistant coating



Cat. No.	Description	Price
S14246	Standard Goggles	1.85/Ea

Standard Size Chemical Goggles

Soft vinyl goggle frame conforms to facial contours for a secure fit.

- Fog-free coated lens and four indirect vents minimize fogging
- Fit comfortably over most prescription eyewear
- Fit most adults and high school students

COMPLIANCE: Meet ANSI™ Z87.1 standards.



Cat. No.	Description	Price
S63095	Standard Size Chemical Goggles	10.25/Ea



Prepare Your Classroom

Find glass and plastic labware in our Lab Supplies & Consumables catalog.

Uvex™ Futura™ Goggles

180° panoramic view for better peripheral vision

- Firm, comfortable seal protects against chemicals, dust, and grindings
- Impact-resistant, chemical-splash polycarbonate lens
- Permanently bonded uvextreme™ AF antifog coating
- Wide, neoprene strap is slotted to prevent slippage



COMPLIANCE: Meet ANSI Z87.1-2010 and CSA Z94.3 standards

Cat. No.	Description	Price
S47619C	Safety Goggles	14.25/Ea
S47619CS	Replacement Strap	3.05/Ea

Colored Safety Goggles

Can be used any time students need reliable eye protection.

- Sized just right for students and adjustable for adults
- Set of six pairs of goggles in three bright, translucent colors (red, blue and green)
- Pack of 6



ORDERING INFORMATION:
Grades K+.

COMPLIANCE:
Meet ANSI™ Z87.1 safety standards.

Cat. No.	Description	Price
S90821	Colored Safety Goggles	30.00/Pk

Safety Goggles

With durable plastic offering good protection to the eyes.

- Flexible body increases wearing comfort
- Great for kids and adults

ORDERING INFORMATION:
Grades K+.

COMPLIANCE:
Meets ANSI™ Z87.1 standards



Cat. No.	Description	Price
S98484	Safety Goggles	5.15/Ea

Allegro™ Eyewear Cleaning Station

Keep eyewear clean and clear at all times.

- Small Disposable Cleaning Station
- Includes an 8 oz. anti-static, anti-fog, non-silicone lens cleaner that removes dirt and oil from glasses and lenses and 600 count easy-dispensing non-abrasive, low-lint 3x 8 in. lens tissues
- Heavy-duty, recyclable cardboard dispenser
- Can be mounted to the wall or used as a tabletop dispenser
- Dimensions:
7W x 8 1/4H x 5 1/8 in.D



Cat. No.	Description	Price
S47601	Disposable Station	9.65/Ea
S47631	100 Pre-moistened Cleaning Wipes	8.50/Pk



Launch Your School's Biotech Program

Find the curriculum, products, and support you need for every step of your journey.

Visit fishersci.com/biotechresources to learn more.

 **fisher science**
education
part of Thermo Fisher Scientific

Waste Boxes

TerraCycle™ Zero Waste Boxes™

Recycle Safety Equipment, Protective Gear, and Consumables. Each box includes signage, a built-in liner, shipping to and from your location, and processing of the materials.

How It Works

- **Order:** Select the appropriate Zero Waste Box system for your discards and place your order.
- **Collect:** Collect the accepted items to fill the system.
- **Ship:** Send it back to TerraCycle with the prepaid return label on your system.*
- **Reorder:** Order your next Zero Waste Box system to continue to recycle everything.

What's Included?

- One Zero Waste Box
- One prepaid return label, affixed to your box
- Processing and recycling of the materials
- Customer support from our team of experts

What You Can Recycle:

Disposable Gloves Box

- Disposable vinyl, nitrile, and latex gloves
- Do not send in other types of protective gear like hair nets, earplugs, and lab coats or garments. Do not include biodegradable gloves.

Hair Nets and Beard Nets Box:

- Any brand and size of used hair and beard nets
- Do not send in gloves, lab coats, or garments

Protective Eyewear Box:

- Swimming, tanning, and safety goggles
- Do not send in other types of protective gear, industrial waste, sunscreen bottles, ear plugs, gloves, lab coats, garments, or glass lenses.

Safety Equipment & Protective Gear Box:

All brands and sizes of disposable safety equipment and protective gear including:

- Ear plugs
- Beard nets
- Hair nets
- Gloves
- Safety glasses
- Disposable Garments
- Masks

Important Shipping Restrictions

For the Zero Waste Box system to work effectively and safely, please follow these simple shipping guidelines on how and what you can ship back to TerraCycle:

- Only send in items this system accepts. Sending in items that fall outside the selected category causes contamination and the need for added separation, which increases costs.
- Do not send in household items that may be considered hazardous waste. These include items such as lithium-ion batteries, pressurized canisters, pesticides, oil-based paint or medical sharps, organics, broken glass, biomedical waste, soiled diapers, or expired medication. Do not send in anything sharp, flammable, reactive, corrosive, ignitable, radioactive, toxic, infectious, or pathogenic.
- Do not send in medical waste. This includes material contaminated with blood or bodily fluids that originate from health care facilities, hospitals, physicians offices, dental offices, blood banks, veterinary clinics, research laboratories, etc.
- Do not ship dripping packages. Please make sure all excess product has been removed before putting it in your Zero Waste Box, Pouch, or Pallet. If you choose to rinse your items, please note that the system must be completely dry prior to shipping.

Cat. No.	Description	Price
<i>Disposable Gloves</i>		
S88854124	Small	135.00/Ea
S88854125	Medium	255.00/Ea
S88854126	Large	410.00/Ea
<i>Hair Nets and Beard Nets</i>		
S88854138	Small	101.00/Ea
S88854139	Medium	130.00/Ea
S88854140	Large	163.00/Ea

Cat. No.	Description	Price
<i>Protective Eyewear</i>		
S88854144	Small	108.00/Ea
S88854145	Medium	143.00/Ea
S88854146	Large	215.00/Ea
<i>Safety Equipment & Protective Gear</i>		
S475652	Small	157.00/Ea
S475653	Medium	300.00/Ea
S475654	Large	492.00/Ea

⚠ Alert: These boxes are not suitable for medical waste, including 'red bag waste', glass lenses, household items that may be considered hazardous waste, or dripping packages. Not available in Hawaii, Alaska, and Puerto Rico.



Hand Protection

fisherbrand Powder-Free Nitrile Gloves

Choose a latex-free glove with exceptional durability and puncture resistance. Fisherbrand™ Powder-Free Nitrile Exam Gloves provide a superior fit, comfort, tactile sensitivity and dexterity.

- Provide superior fit, comfort, tactile sensitivity and dexterity
- Contains no latex proteins that can cause allergic reactions
- Textured surface for solid grip in wet or dry conditions
- 3.5mil thick at palm, 4.7mil at fingertip
- Pack of 100



Cat. No.	Size	Price
S79585XS	X-Small	75.00/Pk
S79585S	Small	75.00/Pk
S79585M	Medium	75.00/Pk
S79585L	Large	75.00/Pk
S79585XL	X-Large	75.00/Pk
S79585XXL	2X-Large	75.00/Pk

fisherbrand Comfort Nitrile Exam Gloves

Choose a glove with enhanced comfort and tactile sensitivity. Fisherbrand™ Comfort Nitrile Exam Gloves are powder-free and latex-free with a beaded fit and textured fingers.



- Powder-free, latex-free
- Exceptional fit, feel and durability
- Ambidextrous with textured fingers and beaded cuff
- Nonsterile
- Length: 9.5 in. (24cm)
- Thickness: 3mil (palm)
- Color: Periwinkle blue
- Pack of 200

Cat. No.	Size	Price
S05009A	Small	93.00/Pk
S05009B	Medium	93.00/Pk
S05009C	Large	93.00/Pk

fisherbrand Powder-Free Nitrile Exam Gloves with Aloe

Help moisturize and soothe dry, chapped hands. Fisherbrand™ Powder-Free Nitrile Exam Gloves with Aloe are the perfect choice for situations that require gloves for an extended period of time. Textured surface for a superior grip.

Offer superior puncture and chemical resistance.

- Unique, 100% latex-free, odor-free synthetic formulation reduces hand fatigue
- Textured surface for superior grip
- Nonsterile, ambidextrous
- Length: 9.5 in. (24.1cm)
- Thickness: 5.6mil
- Green
- Pack of 100



CERTIFICATIONS: 510K approved for medical use

Cat. No.	Size	Price
S66382	Small	85.00/Pk
S66383	Medium	85.00/Pk
S66384	Large	85.00/Pk
S66385	X-Large	85.00/Pk

Kimtech™ Purple Nitrile™ Gloves

Purple nitrile glove offers maximum protection for higher-risk procedures. Unique color differentiates from latex and other synthetic gloves.

- Textured fingertips offers dexterity and tactile sensitivity
- Latex-free
- Ambidextrous
- 9.5 in. length
- Beaded cuff prevent roll-down and allow for easy donning
- Powder-free



Cat. No.	Size	Quantity	Price
S78942	X-Small	100	47.50/Pk
S78943	Small	100	47.50/Pk
S78944	Medium	100	47.50/Pk
S78945	Large	100	47.50/Pk
S78946	X-Large	90	47.50/Pk

MidKnight™ Black Nitrile Exam Gloves

Hide stains and provide a contrasting backdrop for identifying lighter colored material with a distinct black color. MidKnight™ Black Nitrile Exam Gloves offer a high tensile strength compared to both ASTM D6319 and EN455 standards, allowing the hands to move freely and comfortably.

- Fully textured design and non-foaming formula enable wearers to have a firm wet grip when handling objects
- Contain no natural rubber latex and are powder-free to help protect wearers from Type I skin allergies, skin irritation, and dryness
- Recommended for use in environments that demand firm grip and reliability, and where precise contact requires enhanced tactile sensitivity
- Single use only
- Ambidextrous
- Beaded cuff, standard examination length
- Pack of 100



Cat. No.	Size	Price
S28986	X-Small	36.75/Pk
S28987	Small	38.00/Pk
S28988	Medium	37.75/Pk
S28989	Large	36.50/Pk
S28990	X-Large	37.75/Pk
S28991	2X-Large	38.00/Pk

Disposable Polyethylene Food Service Gloves

Resistant to most chemicals, reagents and acids.

For non-precision food preparation and service applications

- Economical, lightweight
- Pack of 100
- Made of polyethylene-blended film
- Extra-long cuffs
- Medium weight (1mil thickness)
- Tough and durable, premium-weight, soft embossed glove
- Easy-on and off for jobs requiring frequent changes
- Textured for enhanced grip and hand comfort



Cat. No.	Size	Price
S47330	Medium	5.00/Pk
S47332	Large	5.00/Pk

SP Bel-Art™ Autoclave Gloves

Excellent for handling hot objects out of the autoclave or oven.

- Heat-resistant up to 232°C (450°F); replace asbestos gloves
- Terry cloth knit gloves
- Provide superior protection in the school laboratory
- Machine wash, air dry
- Bright orange color keeps gloves visible



CERTIFICATIONS:

CE marked.

⚠️ ALERTS:

Product is not waterproof and should not be used to handle dry ice.

Cat. No.	Length	Price
S47334	13 in.	36.50/Pair
S473341	18.5 in.	49.75/Pair

SP Bel-Art™ Bellows Type Glove Box Gloves



Made from natural latex with a removable hand, the sleeves fit 20 cm (8 in.) circular ports and 18 x 25 cm (7 x 10 in.) oval ports.




- Bellows allows cooling air circulation
- Forms a leak-proof unit with O-ring construction integral with the box
- When the glove and sleeve are separated, the sleeve remains attached to the box
- Glove sizes are interchangeable, and almost any standard laboratory, surgical or disposable glove can be used
- Length of the extended sleeve and glove is approximately 76 cm (30 in.)

Cat. No.	Description	Price
S81344	Bellows Type Glove Box Gloves: Size SM	508.00/Ea

United Scientific™ Silicone Rubber Hand Protector



Constructed of silicone rubber with a specially textured grip for secure handling of hot objects such as beakers, flasks, and bottles.



- Silicone rubber is capable of withstanding high temperatures
- Specially textured grip prevents containers from slipping
- Ideal for gripping hot objects such as beakers, flasks, and bottles
- Red
- One size
- Latex-free

Cat. No.	Description	Price
S27814	Silicone Rubber Hand Protector	11.25/Ea

SP Bel-Art™ Hot-Hand Protector Glove/Mitt



Mitt protects hands when handling extremely hot or cold objects, including dry ice



- Use for handling hot and cold objects from -57°C (-70°F) to 260°C (500°F)
- Soft and compliant silicone rubber naturally conforms to hand and object. Gripping surfaces studded with small concave nubs for positive grip
- Easier and faster to use than glove or tongs for many applications
- Tab with hole for hanging hand protector near work station
- One size fits all hands, left or right 10 × 19cm (3.875 × 7.5 in.)

Cat. No.	Color	Price
S47371D	Red	28.25/Ea
S42988	Sky Blue	31.25/Ea
S42989	Lime Green	31.25/Ea

Silicone Hand Protector



Silicone Hand Protector for safe handling of hot and cold objects between -70 and 500°F.

- Pockets for thumb and fingers
- Fits every hand size
- Gripping surface studded for positive grip
- Tabs for hanging
- Size: 4 × 7.5 in.
- Case of 20
- Color: Red



Cat. No.	Description	Price
S81046	Hand Protector	142.00/Cs

Stay Connected with the Fisher Science Education Team

Follow us on



fisher science
education
part of Thermo Fisher Scientific

Aprons & Lab Coats

fisherbrand Disposable Apron

Lightweight plastic apron is sturdy enough to be reused.

- Thickness: 1.0mil
- White
- Sold in handy dispenser pack
- Pack of 50



Cat. No.	Description	Price
S98418	Disposable Apron	30.00/Pk

Eisco™ Clear Vinyl Bib Apron

Clear vinyl bib apron protects clothes during laboratory exercises. Ideal for use in labs, kitchens, or shops.

- Basic transparent vinyl apron
- Perfect for biology, chemistry, and general science laboratories
- Protects clothes from liquid spills and splashes
- Acid resistant



Cat. No.	Size	Price
S24922	Medium	5.50/Ea
S24921	Large	5.50/Ea

Eisco™ Clear Vinyl Apron

Simple clear vinyl apron that is perfect to protect clothes from splashes during science lab activities, washing dishes, crafts or other similar activities. Easily slips over head and ties around the back. Extra small size is suitable for children.

NEW

- Simple clear vinyl apron for labs, dishwashing, or crafts
- Protects clothes from water and chemical splashes
- Fits around neck, ties in back
- Size: extra small (15 in. width, 22 in. length)
- Suitable for children (X-Small)



Cat. No.	Description	Price
S27425	X-Small Clear Vinyl Apron	4.10/Ea

United Scientific™ Vinyl Apron

Chemical-resistant vinyl fabric that stays flexible

- Chemical-resistant
- Flexible, vinyl fabric
- Secure cloth ties at the neck and waist
- Reinforced at points of strain

NEW

Cat. No.	Size	Price
S27819	Small	5.60/Ea
S27820	Medium	6.70/Ea
S27821	Large	7.55/Ea
S27822	X-Large	8.20/Ea

Fisher Scientific™ PVC Aprons

Reusable blue PVC apron with cotton ties at the waist and neck.

Keep a supply of these durable, attractive aprons on hand to protect your body and clothing against any accidental spills in the lab or classroom.



Cat. No.	Size	Price
S47381A	Small	8.65/Ea
S47382A	Medium	9.75/Ea

Fisher Science Education™ PVC Apron of the Elements

Displays current periodic table with properties of 114 elements.

- Reusable
- Blue
- Made of PVC (polyvinyl chloride)
- Cotton ties at the waist and neck



Cat. No.	Size	Price
S32134A	Small	8.70/Ea
S32135A	Medium	9.65/Ea

Fisher Scientific™ Transparent PVC Aprons

Reusable transparent PVC apron with cotton ties at the waist and neck.

Keep a supply of these durable, attractive aprons on hand to protect your body and clothing against any accidental spills in the lab or classroom.



Cat. No.	Dimensions (L x W)	Price
S475651A	36 x 27 in.	8.25/Ea
S475652A	42 x 27 in.	7.10/Ea

Eisco™ Black Rubber Apron

Simple black rubber apron that is perfect to protect clothes from splashes during science lab activities, washing dishes, crafts or other similar activities. Easily slips over head and ties around the back. Extra small size is suitable for children.



- Simple black rubber apron for labs, dishwashing, or crafts
- Protects clothes from water and chemical splashes
- Fits around neck, ties in back
- Size: extra small (15 in. width, 22 in. length)
- Suitable for children (X-Small)
- Contains Latex



Cat. No.	Description	Price
S27424	X-Small Black Rubber Apron	5.90/Ea

Rubberized Chemical Resistant Aprons

Flexible and lightweight.

- Made of chemically resistant rubberized material
- Reinforced at points of strain
- Ties at waist and neck
- Latex-free
- Color: black



Cat. No.	Length	Width	Price
S475611D	30 in.	24 in.	7.50/Ea
S475611A	36 in.	27 in.	7.95/Ea
S475611B	42 in.	27 in.	10.50/Ea
S475611E	42 in.	36 in.	10.50/Ea
S475611C	46 in.	36 in.	11.00/Ea

Eisco™ Rubber Apron Bib

For general lab work, labs, kitchens, or shops.

- Comfortable and flexible
- Waterproof
- Resistant to most chemicals



Cat. No.	Size	Dimensions (L x W)	Price
S24915	Small	30 x 24 in.	5.30/Ea
S24916	Medium	36 x 27 in.	6.50/Ea
S24914	Large	42 x 27 in.	9.60/Ea

Rubberized Cloth Apron

Provides basic protection from chemicals in laboratory or classroom settings. Made of rubberized cloth and available in several sizes.

- Great for laboratory and classroom safety
- Material is rubberized and resistant to most common chemicals
- Cloth tie strings provide steadiness and flexibility
- Available in several sizes

⚠ ALERTS:
Intended for ages 13 and up. May contain latex.



Cat. No.	Dimensions (L x W)	Price
S13706	24 x 30 in.	3.25/Ea
S13707	27 x 36 in.	3.85/Ea
S13708	27 x 42 in.	4.10/Ea
S13709	36 x 42 in.	5.40/Ea
S13710	36 x 46 in.	5.80/Ea

DR Uniforms™ Kids Lab Coat

Fine quality kid lab coats for classroom and/or home use. These lab coats are sturdy, moisture absorbent, and ideal for any class or project.

- Knee length
- Two front pockets, one breast pocket
- Side slits for pant pocket access
- Soft and comfortable twill fabric
- Lead and Phthalate free buttons
- Soft and smooth feel, highly sweat-absorbent



Cat. No.	Size	Price
S30777	7	15.75/Ea
S30778	8/10	15.75/Ea
S30779	12/14	15.75/Ea
S30780	16/18	14.25/Ea

Eisco™ White Lab Coats

White Lab Coats with long sleeves feature concealed press stud fastenings, backvents, and three large pockets

- White, modern style lab coat
- Polyester/cotton twill
- Variety of sizes available

WARRANTY:

1 year limited manufacturer warranty



Cat. No.	Size	Price
S41146	Small	18.50/Ea
S41142	Medium	18.50/Ea
S41145	Large	18.50/Ea
S41147	X-Large	18.50/Ea

DR Uniforms™ Unisex Lab Coat

DR Uniforms™ Unisex Lab Coat, 60/40 CottonPoly Blend features three patch pockets, two side slits and four buttons. Cool, comfortable and breathable throughout the day.

- Unisex
- 60% Cotton, 40% Polyester blend - the higher cotton content makes the coat breathable, moisture absorbent, soft and comfortable
- Stays cool throughout the day
- Sizes XS - 5XL
- Full-length (42 inches)
- 3 patch pockets
- 2 side access slits



Cat. No.	Size	Price
S72134	2X-Small	22.75/Ea
S72135	X-Small	22.75/Ea
S72136	Small	22.75/Ea
S72137	Medium	22.75/Ea
S72138	Large	22.75/Ea
S72139	X-Large	22.75/Ea
S72140	2X-Large	24.75/Ea
S72141	3X-Large	24.75/Ea
S72142	4X-Large	28.50/Ea
S72143	5X-Large	28.25/Ea

United Scientific™ Lab Coats

Well fitted and comfortable.

- Made of easy-care white cotton blend fabric (65% polyester, 35% cotton)
- Knee-length
- One upper left breast pocket and two roomy lower pockets
- With side slit openings and adjustable belt
- Button Closures



Cat. No.	Size	Price
S04775	X-Small	29.50/Ea
S04776	Small	29.50/Ea
S04777	Medium	29.50/Ea
S04778	Large	29.50/Ea
S04779	X-Large	30.25/Ea
S04780	2X-Large	36.75/Ea

FISHER SCIENCE EDUCATION NURSING AND CTE EDUCATION PROGRAMS

Career and Technical Education (CTE) gives students a head start on a future career in a STEM field. The Fisher Science Education brand offers the tools educators need to get them ready.

Visit fisheredu.com/nursingcte for a full list of products available to support your CTE program for nursing, medical technology, and more.



fisher science
education
part of Thermo Fisher Scientific



YOUR STEM PROGRAM
MAKE IT AMAZING!

fisheredu.com/STEM: Your source for products, professional development, resources and everything STEM.

Facilities

Social Distancing Aids

Personal Protection Kit



Contains face masks and gloves for individual protection against dangerous germs in addition to antimicrobial wipes and hand sanitizer for disinfection of hands or surfaces.

- Protection Kits contain Face Masks, Disposable Gloves, Hand Sanitizer Packets or Liquid Hand Sanitizer, and Antimicrobial Wipes
- Available for 1 person, or 10 people
- 8 x 10 in. clear bag included




Cat. No.	Description	Includes	Price
12888601	1 Person Kit	1.75 oz Liquid Sanitizer	18.00/Ea
12888600	1 Person Kit	2 Hand Sanitizer Packets	9.05/Ea
12888602	10 Person Kit	1.75 oz Liquid Hand Sanitizer	73.00/Ea

BenchPro™ FeatherLite™ Panels



Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.



Clean with soap and water. Antibacterial sprays are ok to use as well.

Cat. No.	Height	Price
Acrylic		
S81594	35 in.	257.00/Ea
S81593	52.5 in.	322.00/Ea
S81592	70 in.	441.00/Ea
Vinyl		
S81591	35 in.	148.00/Ea
S81590	52.5 in.	165.00/Ea
S81589	70 in.	178.00/Ea

BenchPro™ Edge Mounted Sneeze-Guard



Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.



Cat. No.	Height	Width	Price
S81604	18 in.	24 in.	43.75/Ea
S81606	18 in.	30 in.	60.25/Ea
S81605	24 in.	24 in.	54.00/Ea
S81607	24 in.	30 in.	74.50/Ea

Dynalon Countertop Shield



Dynalon Countertop Shield are as clear as glass to keep workers and customers safe.

- Standard protection shield
- Bottom opening for exchange of materials
- Protect important office personnel from splash contamination
- Bottom opening for exchange of materials



- Available in two sizes: 36 and 47 in. width
- Removable feet for flat shipment and storage
- Manufactured of 1/4 in. clear acrylic, adding strength and stability along with the mobility to move where necessary.
- Ideal for office check-in and receptionist staff needing splash protection from the public.

Cat. No.	Dimensions (L x W)	Price
S27930	36 x 24 in.	258.00/Ea
S27931	48 x 30 in.	399.00/Ea

Log in to your account on fisheredu.com to see your discounts.

Eisco™ Polycarbonate Safety Shield

Polycarbonate Safety Shield includes adjustable outer panels to accommodate a variety of lab set ups.

- The width of the outer panels is about 8.9 in. (22.5 cm) and the center panel is 11.8 in. (30 cm) wide.
- All panels are 4mm thick and provide excellent visual clarity
- Adjustable
- For classrooms and labs
- Provides clear, easy viewing during experiments
- Mounting clamps are not required for use, but will add stability (not included)



WARRANTY: 1 year limited manufacturer warranty

Cat. No.	Description	Price
S41719	Polycarbonate Safety Shield	85.00/Ea

BenchPro™ Upright Mounted Sneeze-Guard

Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.



Cat. No.	Height	Width	Price
S81595	14 in.	48 in.	61.50/Ea
S81598	14 in.	60 in.	116.00/Ea
S81601	14 in.	72 in.	168.00/Ea
S81596	26 in.	48 in.	77.00/Ea
S81599	26 in.	60 in.	144.00/Ea
S81602	26 in.	72 in.	210.00/Ea
S81597	38 in.	48 in.	94.00/Ea
S81600	38 in.	60 in.	171.00/Ea
S81603	38 in.	72 in.	251.00/Ea

BenchPro™ Free Standing Sneeze-Guard

Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.



Cat. No.	Height	Length	Price
S81575	24 in.	18 in.	40.25/Ea
S81576	24 in.	24 in.	57.50/Ea
S81578	24 in.	30 in.	67.75/Ea
S81580	24 in.	48 in.	141.00/Ea
S81579	38 in.	30 in.	139.00/Ea
S81577	40 in.	24 in.	101.00/Ea

BenchPro™ Free Standing Sneeze-Guard with Pass Through



Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.



Cat. No.	Height	Length	Price
S81581	24 in.	18 in.	46.25/Ea
S81582	24 in.	24 in.	57.75/Ea
S81584	24 in.	30 in.	66.75/Ea
S81586	24 in.	48 in.	141.00/Ea
S81585	38 in.	30 in.	138.00/Ea
S81583	40 in.	24 in.	99.00/Ea

BenchPro™ Ceiling Sneeze-Guard



Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.



Cat. No.	Height	Price
S81587	24 in.	57.00/Ea
S81588	48 in.	142.00/Ea

BenchPro™ Sliding Sneeze-Guard



Get protected while you are working with this BenchPro™ Sneeze Guard Panel. It is easy to setup and easy to clean.

Clean with soap and water. Antibacterial sprays are ok to use as well.

- White ABS



Cat. No.	Height	Length	Price
S81608	24 in.	18 in.	6.70/Ea
S81609	24 in.	24 in.	8.95/Ea
S81611	24 in.	30 in.	18.00/Ea
S81613	24 in.	48 in.	15.00/Ea
S81612	38 in.	30 in.	17.75/Ea
S81610	40 in.	24 in.	11.25/Ea

Cleaners and Sanitizers

United Scientific™ Lab Clean 70% IPA Wipes



General purpose disposable 70% Isopropyl Alcohol (IPA) wipes are ideal for a variety of cleaning applications, including electronics, tools, workspaces, and other surfaces that require quick removal of contaminants. Suitable for Class 6-7 Cleanrooms.



- General purpose disposable 70% Isopropyl Alcohol (IPA) wipes
- Ideal for a variety of cleaning applications, including electronics, tools, workspaces, and other surfaces that require quick removal of contaminants
- Suitable for Class 6-7 Cleanrooms
- 18 cm x 14 cm (7 x 5.5 in.) perforated wipes are constructed of durable, lint-free, tear-resistant polyester
- Suitable for use in Class 6-7 Cleanrooms.
- Wipes are presaturated in a solution of 70% Isopropyl Alcohol (CAS: 63-67-0) and 30% Deionized Water
- Foil-sealed canister contains 100 wipes
- Canister includes a screw-top and a convenient flip lid with pull-through dispenser

Cat. No.	Description	Price
S83788	Lab Clean 70% IPA Wipes	15.25/Ea

Antimicrobial Hand Sanitizer



Antimicrobial hand sanitizer is intended for use as a no-rinse hand sanitizer and meets or exceeds the requirements of the USDA to be used as an E3 antimicrobial hand sanitizer.



- Complies with the W.H.O. specifications and requirements for an effective hand sanitizer
- Consistency similar to water
- This high-quality hand sanitizer is a pH balanced formula that is a mild blend of skin emollients and alcohols and contains glycerine
- Leaves hands feeling soft and moisturized
- The product is mild on skin, yet tough enough to exceed FDA/USDA efficacy guidelines for E3 hand sanitizer and antiseptic hand treatments
- Viscosity: 5 to 15 cSt

Cat. No.	Description	Price
S27937	Size: 100 mL Poly Bottle with Pump Spray	7.10/Ea
S27938	Size: 500 mL Poly Bottle with Flip Cap	10.75/Ea

Feel Welcome Hand Sanitizer Gel and Pump

Gel formula hand sanitizer (70% ethyl alcohol) with aloe and vitamin E



- Kills germs and protects your hands
- Made in the USA according to FDA guidelines
- Independently tested and verified by a third party lab
- 12-888-610 is a dispenser pump

Cat. No.	Size	Price
12888608	8 oz.	8.00/Ea
12888609	1 gal.	62.50/Ea
12888610	Hand Sanitizer Pump; 1 Gal.	6.00/Ea

Diversified Spaces™ Hand Sanitizer Unit



A complete hand sanitizer unit that includes the stand and dispenser. The stand is suitable for indoor and outdoor use, with a durable powder-coated base and mounting plate to withstand the elements. The dispenser mounts on the stand with an universal mounting plate that fits most dispensers.



- Dispenser Includes a Refillable Plastic Bottle (1000 mL Capacity)
- Easy to Assemble
- Heavyweight, Free-standing Design
- Indoor/Outdoor Stand & Dispenser
- Powercoated Base
- Requires 4 x C Batteries (Not Included)
- Stands 49.5 in. Tall When Assembled
- Touchless, Motion-Activated Gel Dispenser
- Universal Mounting Plate with Theft Resistant Screws or Double-sided Tape (not included)

Cat. No.	Description	Price
S81425	Hand Sanitizer Unit (Stand and Dispenser)	183.00/Ea
S81426	Hand Sanitizer Dispenser Only	110.00/Ea

Diversified Spaces™ Hand Sanitizer Dispenser



A motion-activated, touchless gel dispenser that provides clean, accessible hand sanitization. It comes with a refillable, 1000 mL capacity plastic bottle (hand sanitizer is not included).



- 1000 mL Refillable Plastic Bottle Included
- Hand Sanitizer Gel is Not Included
- Hands-free, Touchless, Motion Activated Design
- Requires 4 x C Batteries
- Universal Gel Dispenser
- Works with Standard Hand Sanitizer Gels

Cat. No.	Description	Price
S81427	Hand Sanitizer Dispenser	86.00/Ea

Eisco™ Goggle Sanitizer Cabinet

UV light kills the germs on goggles after laboratory use.



- Wall-mounted or bench top unit runs on 120VAC
- Heavy gauge CRC steel treated with white enamel finish
- Unit has doors and includes 15 minute timer with UV bulb
- UV light operates only when doors are closed
- Light leak avoided by tight fitting piano style hinges
- Sterilization process is indicated by indicator light
- Use to sanitize 36 pairs
- 120VAC

Germicidal Replacement Bulb

- For use with Goggle sanitizer cabinet (S07817).
- Philips T8 TUV bulb
- Fluorescent

Cat. No.	Description	Price
S07817	Goggle Sanitizer Cabinet	711.00/Ea
S24853	T8 TUV Fluorescent Replacement Bulb	65.50/Ea

Economy Ultraviolet Goggle Sanitizers

Sanitize and store 30 to 40 goggles or safety glasses in the same cabinet.



- High-quality, 24-gauge, white baked enamel sheet metal
- 15-minute mechanical timer gives controlled treatment
- Two safety door switches, key lock and pilot light
- Long-life ultraviolet bulbs
- S47604B has a reflective aluminum interior
- Replacement bulb S47604BLB for Economy Ultraviolet Goggle Sanitizers
- For use with Goggle Sanitizers (S47604A and S47604B)

ORDERING INFORMATION:

Shipped assembled and ready for use. UL-approved. Electrical requirements: 110V 50/60Hz.

Cat. No.	Description	Price
S47604A	Goggle Sanitizer for 30 Goggles	1154.00/Ea
S47604B	Goggle Sanitizer for 40 Goggles	1211.00/Ea
S47604BLB	Replacement Bulb; Length: 18 in.	48.00/Ea

Premoistened Lens Cleaning Towelettes

Fast and easy way to keep lenses clean.

- For use on glass, plastic or polycarbonate safety goggles, shields and computer screens
- Help reduce eyestrain caused by dirty eyewear
- Excellent cleaning without streaking
- Silicone-free, antistatic, antifog formula
- Pack of 100

ORDERING INFORMATION:

Individually foil wrapped. Packed in ready-to-use dispenser pack with built-in eye hook for easy wall hanging.



Cat. No.	Description	Price
S01814A	Premoistened Lens Cleaning Towelette	12.10/Pk

Sinks and Eyewash Stations

Diversified Spaces™ Stainless Steel Mobile Sinks



Mobile Stainless-Steel Sink provides a way for people to wash their hands when there is no plumbing around. Requires electricity to operate.

- Easy to Refill Fresh Water – No Disassembly Needed
- 115 Volt 3.3 GPM, 45 PSI Water Diaphragm Pressure Pump (CSA Listed)
- 15 ft. GFCI Power Cord with Standard Wall Plug (UL Listed)
- 18 Gauge Stainless Steel Body
- 4 Gauge Stainless Steel Top
- 5 Gallon Removeable Fresh Water Tank
- 7 Gallon Removeable Waste Water Tank
- Drain Hose for Easy Draining
- Four Locking 5 in. Casters for Easy Mobility (2 Swivel, 2 Locking)
- Heavy Duty Stainless Steel Construction
- Hose Connections
- Lockable Access Door
- Sink Size: 10 x 14 x 5 in.



- Soap & Towel Dispenser
- Soap Dispenser (Top-Loading, Refillable)

Cat. No.	Description	Price
S81500	Cold Water Sink with Sensor Faucet	5499.00/Ea
S81501	Hot Water Sink with Sensor Faucet	6324.00/Ea
S81498	Metered Cold Water Sink	5053.00/Ea
S81499	Metered Hot Water Sink	5888.00/Ea

Facewash/Eyewash

Access quick, effective and dependable water coverage for the face and eyes.

- Provides a flow of 8 gpm to wash eyes and face with the extra-large, ABS plastic dual-stream head with convenient float-off covers
- Highly visible, yellow ABS plastic bowl
- Chrome-plated brass valve stays open until manually closed
- Self-adjusting regulator assures even flow
- Includes wall brackets and safety sign



COMPLIANCE: Meets ANSI Z358.1-2009 standards.

Cat. No.	Description	Price
S477014	Face/Eye Wash	523.00/Ea

Diversified Spaces™ Shower/ Eyewash Station

Provides convenient storage for lab safety supplies.

- Emergency drench shower head
- High impact ABS plastic
- Chemical resistant earth friendly UV finish
- Stay open chrome-plated brass valve
- Eyewash unit: Stainless steel bowl (Dia.: 10 in.)
- Lever action ball valve on plastic laminate top
- Two adjustable shelves

COMPLIANCE:

SEFA

CERTIFICATIONS:

MAS Certified Green



Cat. No.	Description	Price
S38038	Shower Station	2353.00/Ea
S38039	Shower Station with Locks (for security or access for a wheelchair)	2584.00/Ea

Faucet-Mount Eyewash

Eyewashes attach easily to any standard laboratory faucet; no special plumbing needed.

- Made of chrome-plated brass with float-off dust covers
- Back-mount unit provides 1.5 gpm
- Angled unit has flow control

INCLUDES:

Common adapters that fit most standard and gooseneck faucets

COMPLIANCE:

Meet ANSI Z358.1-1998 standards

⚠ ALERTS: To avoid scalding, cold water or tempered water connection is recommended



Cat. No.	Description	Price
S47621	Faucet-Mount Eyewash	163.00/Ea



FEATURED SUPPLIERS MORE FOR YOU

Visit fisheredu.com/featuredsuppliers to learn more about enhanced products and promotional information.

Safety Cans

SP Bel-Art™ Benchtop Biohazard Disposal Can with Motion Sensor Lid

Compact benchtop biohazard disposal can contains odors and is perfect for small items such as pipette tips and microcentrifuge tubes.

- Automatic opening lid provides hands free use; closes automatically and securely for a more sanitary environment
- Holds a 25 × 25cm (10 × 10 in.) autoclavable biohazard disposal bag (S43506)
- Includes a polypropylene liner with handle for easy emptying of contents
- Aesthetic and durable; stainless steel body and chrome lid are resistant to bleach and easy to clean
- Capacity 1.5L



REQUIRES: (3) AAA batteries (not included)

COMPLIANCE: CS & RoHS

Cat. No.	Description	Price
S72309	Benchtop Biohazard Disposal Can with Motion Sensor Lid	52.50/Ea

Justrite™ Nonmetallic Oval Type I Safety Cans For Flammables



30% more storage capacity than round shape models.

- Constructed of high-density polyethylene
- Trigger mechanism enables pouring of small quantities of liquid
- Leak-tight, self-closing cap prevents spillage and automatically vents to guard against explosion
- Dual-density flame arrester protects against external heat sources
- Current carrying carbon insert imbedded into can for proper grounding
- Nonmetallic, stainless-steel hardware
- Red



Cat. No.	Description	Price
S04768MF	1 Gallon Safety Can	98.00/Ea
S04769MF	1/2 Gallon Safety Can	112.00/Ea

Eisco™ Laboratory Safety Chart

Laboratory safety chart with colored illustrations can be used in any school laboratory.

- Laboratory safety chart with colored illustrations
- Printed on 29 × 39 in. polyart plastic sheet fitted with plastic roller
- Communicates vital safety information
- Ideal for school laboratories
- Dimensions (L × W): 29 × 39 in.



Cat. No.	Description	Price
S26685	Laboratory Safety Chart	26.50/Ea

Universal Mini Chemical Spill Kit

Be prepared for spills in your lab. Universal sorbent eliminates guesswork — safe for use on all chemicals except hydrofluoric acid. Kit contains enough material to absorb 1L of solvent.



Includes:

- 1, 1kg Universal sorbent
- 1 pair Nitrile gloves
- 1 pair splash resistant goggles
- 1 dustpan and brush set
- 2 Poly bags 9 x 16
- 2 twist ties
- 2 blank tags
- Instruction sheets
- SDS

DISCLAIMERS:

WARNING: This product can expose you to chemicals including respirable crystalline silica, which is known to the State of California to cause cancer. Not for use on hydrofluoric acid.

Cat. No.	Description	Price
S28619*	Mini Chemical Spill Kit	27.75/Ea



* California - WARNING: Cancer www.P65Warnings.ca.gov

** California - WARNING: Reproductive Harm www.P65Warnings.ca.gov

*** California - WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Acid, Caustic and Solvent Combination Spill Kit

Kit contains enough material to clean up a liter of four solvents in your lab



Kit Includes:

- 1100 g Diatomaceous Earth (Flammable Absorbent Material)
- 1.0 kg Acid Spill Cleanup Kit (Acid Neutralizing Mixture)
- 1.2 kg Spill Kit, Caustic Cleanup (Caustic Neutralizing Mixture)
- 200 g Vermiculite (Absorbent)
- Splash Resistant Goggles
- 24 Nitrile Gloves
- 1 Dustpan and Brush Set
- 1 Spill Kit Storage Container

DISCLAIMERS:

WARNING: This product can expose you to chemicals including Respirable Crystalline Silica, which is known to the State of California to cause cancer.

Custom-Built Kits for Your Classroom



The Fisher Science Education team is proud to bring you the convenience of custom educational kits for elementary- to college-level instruction.

Visit fisheredu.com/customkitting to learn more.



Cat. No.	Description	Price
S20636*	Combination Spill Kit; Acid Caustic and Solvent	148.00/Ea

Engage Your Students Wherever They Are

Keep your students safe and engaged with a selection of remote learning tools and materials.

You'll find:

- Products for real-time instruction
- Activities for independent learning
- Physical distancing aids
- And more

Visit fisheredu.com/remotelearning to learn more.



Storage Cabinets

Justrite™ Sure-Grip™ EX™ Safety Cabinets with Handles for Flammables



Nonslip gripper allows fingertip opening.

- Sure-Grip EX flush-mounted door handle with two keys for security
- Handle can be padlocked (padlock not included), no interference with traffic
- Exclusive self-latching feature ensures safe closing of cabinet doors without manually engaging door handle
- Double-wall construction provides 1.5 in. (3.8cm) of insulating air space to help prevent spread of flames
- High-visibility Haz-Alert™ reflective warning labels offer added safety during power outages or fires
- Dual vents with flame arresters at top and opposite bottom allow hazardous vapors to be piped away
- Easy-closing, self-latching doors facilitate fire protection of the cabinet
- Self-latch engages critical three-point latching system for protection under fire conditions
- Two-door models hinged left and right with self-closing mechanism ensures proper closing sequence
- Fusible links hold self-closing doors wide open, flush with front of cabinet for easy access and resupply; automatically release and slam doors shut at 165°F (74°C)
- Concealed self-close mechanism provides optimum storage space and obstruction-free access to contents
- SpillSlope™ safety shelves direct incidental spills to the back and bottom of leakproof sump
- Adjustable leveling feet ensure solid positioning and alignment for tight door closure
- Exterior side panel has built-in grounding connector for easy grounding
- Finished inside and out with durable lead-free powder paint
- Meet OSHA and NFPA requirements
- Color: Yellow



WARRANTY: 10 years

COMPLIANCE:

Meet NFPA Code 30 and OSHA regulations; models with self-close doors also comply with the International Fire Code and NFPA1, FIRE CODE™.

CERTIFICATIONS: FM-approved

Cat. No.	Capacity	Door Style	No. of Doors	No. of Shelves	Price
S04760MF	30 gal.	Manual-close	2	1	1414.00/Ea
S04761MF	45 gal.	Manual-close	2	2	1738.00/Ea
S04763MF	45 gal.	Self-closing	2	2	1915.00/Ea
S04762MF	45 gal.	Self-closing	1	2	2890.00/Ea

Justrite™ Sure-Grip™ EX™ Countertop Safety Cabinets for Flammables



For storing small quantities of flammable liquids at work station or adding to existing storage capacities.

- Countertop style cabinet
- Supplied with one adjustable shelf
- Allows increased access to often used liquids, without compromising safety
- Portable size keeps work areas flexible
- Color: Yellow
- 13.75D x 13.75W x 18.5 in.H
- 4 gal. capacity



INCLUDES:

- Optional wall hanger assembly for wall mounting

Cat. No.	Weight	Certifications/Compliance	Price
S04772MF	62 lb.	NFPA, OSHA	943.00/Ea
S04773MF	64 lb.	FM, NFPA, OSHA, U (NFPA)	996.00/Ea

Justrite™ Sure-Grip™ EX™ Blue Steel Safety Cabinets for Corrosives



Epoxy baked-on powder coat finish provides increased chemical resistance and resists solvents used in labs.

- Double-wall construction and dual vents
- Grounding wire connections and leakproof sills
- Three-point self-latching doors and leveling feet
- Available in different capacities and varying numbers of adjustable shelves
- Polyethylene trays attached to galvanized steel shelves and separate polyethylene liner for bottom sump
- Extra polyethylene top work tray in 30 gal. (15L) two-door cabinet



Cat. No.	Load Bearing Capacity	No. of Doors	No. of Shelves	Price
S04767MF	45 gal. (170 L)	2	2	2081.00/Ea
S04766MF***	30 gal. (113.5 L)	1	1	2779.00/Ea

Acid Cabinets

Polypropylene lock assembly, padlock with two keys.



- Constructed of 1 in. (2.5cm) thick plywood with acid-resistant epoxy paint
- Wooden dowel hinge assembly eliminates the problem of hinge deterioration
- Single fixed shelf is positioned for storage of large 4L bottles on top shelf and 2.5L bottles on bottom shelf
- Recessed top with polypropylene tray, holds 12L of spilled liquid
- Floor includes 2 in. (5.1cm) liquid-tight trough
- Storage capacity: thirty 2.5L or ninety 500mL bottles
- Partially lined model has polypropylene on the floor, shelf and top tray
- Fully lined model has polypropylene on the floor, interior walls, interior top and top tray

ORDERING INFORMATION:

- Exterior: 20D × 31W × 36.5 in.H (50.8 × 78.7 × 92.7cm)
- Interior: 17.13D × 29W × 28.5 in.H (43.5 × 73.7 × 72.4cm)
- Additional polypropylene slide-in isolation compartment available as an accessory (Cat. No. S51074)

Cat. No.	Description	Price
S51072MF	Partially-Lined Acid Cabinet	1490.00/Ea
S51073MF	Fully-Lined Acid Cabinet	1980.00/Ea



FISHER SCIENCE EDUCATION
MAKING SCIENCE MATTER

Visit us online at fisheredu.com to learn more about our programs, services and resources.



* California - WARNING: Cancer www.P65Warnings.ca.gov
 ** California - WARNING: Reproductive Harm www.P65Warnings.ca.gov
 *** California - WARNING: Cancer and Reproductive Harm -www.P65Warnings.ca.gov



Justrite™ Polyethylene Corrosives/Acid Cabinet

Sturdy, completely metal-free poly cabinet for storing acids or bases without the worry of corrosion



- No metal parts to corrode
- Holds up to 36, 2.5L (90L) bottles
- Offers chemical resistance against spills or potential damaging vapors when storing acids or bases
- Shelf adjusts on 3-inch (76mm) increments, and has a load capacity of 125 lb. (57kg)
- Removable yet sturdy sump cover, which also serves as a bottom shelf, conceals the fully integrated, leakproof sump with an 8.5 gal (32L) capacity
- Recessed channels in shelf and base help direct spills to the back of the cabinet and down into the sump
- Three application-specific labels identify contents to safely store and segregate incompatible chemicals — ACID, BASE, and CORROSIVE
- Doors open nearly 180 degrees for easy access, engage firmly with the body when closed, and accept a padlock (not included) for added security
- Cabinet can be used for stand-alone storage
- Grooved top work surface
- Recessed toe-kick allows the user to stand closer to the counter for improved ergonomics
- Two cabinets can be stacked — top cabinet offers convenient eye-level access to contents
- Four pre-formed vent locations at rear provide a guide for an aftermarket vent kit (not included) in the event the application requires a pressurized vapor removal system
- Holds up to 36, 2.5L (90L) bottles
- Exterior: 36W x 35H x 25D
- Interior: 31.25W x 28H x 18D
- Adjustable shelf

⚠ ALERTS:

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices

Cat. No.	Description	Price
S27960MF	Polyethylene Corrosives/Acid Cabinet	2001.00/Ea

Flammables Cabinets

Four adjustable leveling feet.

- Top and bottom vents are fitted with 2 in. (5cm) pipe thread on sidewalls. Double walls of 18-gauge steel provide 1.5 in. (3.8cm) air space
- 2 in. (5cm) Raised sill contains spills
- Adjustable galvanized steel shelves supported by four brackets hold over 500 lb. (227kg) each
- Self-closing doors have fusible link that melts at 71°C (160°F)
- Three-point key lock
- Labeled “FLAMMABLE — KEEP FIRE AWAY”



Justrite™ Wood Laminate Storage Cabinets for Corrosives and Acids

Suitable for storing highly corrosive liquids such as phenol, nitric acid and sulfuric acid.

- Long-lasting wood construction provides exceptional strength
- Stainless-steel handles and hinges resist corrosion from harsh acids
- Shelves adjust on 32mm centers to accommodate wide variety of containers
- Sturdy clip system adjusts shelves if required
- Two doors
- Color: blue



INCLUDES:

Adjustable leg levelers, patent pending Haz-Alert™ reflective labels and double key set for safety and security

COMPLIANCE: FM approved, except S51021MF; meet OSHA and NFPA Code 30 requirements

Cat. No.	Load Bearing Capacity	Door Style	No. of Doors	No. of Shelves	Price
S51015MF	30 gal.	Manual-close	2	1	1653.00/Ea
S51020MF	45 gal.	Manual-close	2	2	2049.00/Ea
S51021MF	45 gal.	Self-closing	2	2	2224.00/Ea

Countertop Model

- Ideal for smaller quantities used at workbench level
- No shelves

Cat. No.	Description	Dimensions (D x W x H) Interior	Load Bearing Capacity	Price
S04764MF	Countertop Storage Cabinet	15.28 x 22.48 x 19.44 in. (38.8 x 57.1 x 49.4 cm)	6 x 2.5 L Bottles	1067.00/Ea

HEADLINE Discoveries

Get Real-World Science
in **Your Inbox**

Get the latest science news, classroom activities, and cutting-edge teaching tools delivered straight to your inbox every month.

Visit fisheredu.com/HDSUBscribe to sign up for this free email.

 **fisher science**
education
part of Thermo Fisher Scientific

Numeric	F	U
3M™	Fisher Science Education™	United Scientific™
N95 Respirators 13	Aprons 22	Alcohol Wipes 27
D	Fisher Scientific	Aprons 3, 22
Diversified Spaces	Aprons 22-23	Goggles 3
Eyewash Stations 29	Fisherbrand™	Lab Coats 24
Sanitizer Dispensers 27-28	Aprons 22	Ultraviolet Light 4
Sinks, Mobile 29	Bandages 5	
E	Face Shields 13	
Eisco	Gloves 19	
Aprons 22-23	Goggles 14-16	
Facemasks 10, 13	Safety Glasses 14	
Cabinets, Sanitizer 28	S	
Goggles 15-16	SP Bel-Art	
Lab Coats 24	Biohazard Disposal Cans 30	
Safety Glasses 15	Gloves 20	
Safety Shields 26	Hand Protectors 21	
Safety Wall Chart 30	Mitts 21	
Wash Bottles 3	Supertek	
	Sanitizer Cabinets 3	

A		
Adhesive Bandages		5
Alcohol		
Swabs		5
Wipes		27
Aprons		3, 22-23
Autoclave Gloves		21
B		
Bandages		5
Bibs		22-23
Biohazard Disposal Cans		30
C		
Cabinets		
Acid Storage		32-34
Flammable Storage		32, 34
Goggle Sanitizer		28
Cans		
Biohazard Disposal		30
Flammable Storage		30
Coronavirus Safety Signs		3
D		
Disposal Boxes		
Eyewear		18
Facemasks		13
Gloves		3, 18
Hair Nets & Beards		18
PPE		3, 18
E		
Evacuation Chairs		4
Eyewash Stations		29
Eyewear		
Cleaning Station		17
Disposal Box		18
Eyeglasses		14-15
F		
Face Shields		13
Facemasks		10, 13
Fire Blankets		4
First Aid		
Alcohol Swabs		5
Backpack Kits		7-9
Bandages		5
Bleeding Control Kits		7
Evacuation Kits		8
Lockdown Kits		9
Mobile Incident Command Stations		7
Mobile Trauma Stations		6-7
Scissors		5
Shelter-in-Place Kits		8-9
Thermometer		5
Trauma Kits		7-8
Trauma Stations		7
Kits		4
G		
Gloves		
Aloe		19
Autoclave		20
Blue		19
Disposal Box		3, 18
Exam		19, 20
Food Service		20
Glove Box Gloves		20
Nitrile		19-20
Polyethylene		20
Powder-free		19
Purple		19
Goggles		3, 14-17
Sanitizer Cabinet		28
H		
Hair Nets & Beards		
Disposal Box		18
Hand Protectors		
Hand Sanitizers		
Dispensers		27-28
Solutions and Gels		27
L		
Lab Coats		
Adult		24
Childrens		23
Lens Cleaning Towelettes		28
M		
Magnifying Eyewear		15
Mask Disposal Box		13
Mitts		21
N		
N95 Respirators		9
P		
PPE		
Disposal Boxes		3, 18
R		
Respirators		11, 13
S		
Safety Glasses		13-15
Safety Wall Chart		30
Sanitizer Cabinets		3
Scissors		5
Sinks, Mobile		29
Social Distancing Aids		
Countertop Shields		25
Free-standing Shields		26
Mounted Shields		26
Panels		25
Pass-through Shields		26
Personal Protection Kits		25
Shields		26
Sneeze Guards		25-26
Suspended Shields		26
Spill Kits		
Swabs		
T		
Thermometer		5
U		
Ultraviolet Light		4
W		
Wash Bottles		3
Waste Boxes		13, 18
Wipes		27

TERMS AND CONDITIONS OF SALE

Unless otherwise expressly agreed in writing, all sales are subject to the following terms and conditions:

GENERAL: Fisher Scientific Company L.L.C. ("Seller") hereby offers for sale to the buyer named in the order ("Buyer") the products listed on the face hereof (the "Products") on the express condition that Buyer agrees to accept and be bound by the terms and conditions set forth herein. Any provisions contained in any document issued by Buyer are expressly rejected and if the terms and conditions in this Agreement differ from the terms of Buyer's offer, this document shall be construed as a counter offer and shall not be effective as an acceptance of Buyer's document. Buyer's receipt of Products or Seller's commencement of the services provided hereunder will constitute Buyer's acceptance of this Agreement. This is the complete and exclusive statement of the contract between Seller and Buyer with respect to Buyer's purchase of the Products. No waiver, consent, modification, amendment or change of the terms contained herein shall be binding on Seller unless in writing and signed by Seller and Buyer. Notwithstanding the foregoing, Seller reserves the right to make changes to these Terms and Conditions of Sale at any time by posting such changes here, and such changes shall apply to orders of Products made by Buyer thereafter. Seller's failure to object to terms contained in any subsequent communication from Buyer will not be a waiver or modification of the terms set forth herein. All orders are subject to acceptance in writing by an authorized representative of Seller. Seller reserves the right to reject any order for any reason. By purchasing products from Seller, Buyer agrees that it is purchasing the Products only for its own use and not for purposes of resale. Products purchased pursuant to these Terms and Conditions may not be sold or transferred to any person or entity for purposes of resale. Seller reserves the right to take any action it deems necessary to address Buyer's violation of these terms, including, but not limited to, the suspension of Buyer's ability to purchase Products.

PRICE: All prices published by Seller or quoted by Seller's representatives may be changed at any time without notice. All prices quoted by Seller or Seller's representatives are valid for thirty (30) days, unless otherwise stated in writing. All prices for the Products will be as specified by Seller or, if no price has been specified or quoted, will be Seller's price in effect at the time of shipment. All prices are subject to adjustment on account of specifications, quantities, raw materials, cost of production, shipment arrangements or other terms or conditions which are not part of Seller's original price quotation. List prices may vary from the price displayed on the website for certain products.

TAXES AND OTHER CHARGES: Prices for the Products exclude all sales, value added and other taxes and duties imposed with respect to the sale, delivery, or use of any Products covered hereby, all of which taxes and duties must be paid by Buyer. If Buyer claims any exemption, Buyer must provide a valid, signed certificate or letter of exemption for each respective jurisdiction.

TERMS OF PAYMENT: Seller may invoice Buyer upon shipment for the price and all other charges payable by Buyer in accordance with the terms on the face hereof. If no payment terms are stated on the face hereof, payment shall be net thirty (30) days from the date of invoice. If Buyer fails to pay any amounts when due, Buyer shall pay Seller interest thereon at a periodic rate of one and one-half percent (1.5%) per month (or, if lower, the highest rate permitted by law), together with all costs and expenses (including without limitation reasonable attorneys' fees and disbursements and court costs) incurred by Seller in collecting such overdue amounts or otherwise enforcing Seller's rights hereunder. Seller reserves the right to require from Buyer full or partial payment in advance, or other security that is satisfactory to Seller, at any time that Seller believes in good faith that Buyer's financial condition does not justify the terms of payment specified. All payments shall be made in U.S. Dollars.

DELIVERY; CANCELLATION OR CHANGES BY BUYER: The Products will be shipped to the destination specified by Buyer, F.O.B. Seller's shipping point. Seller will have the right, at its election, to make partial shipments of the Products and to invoice each shipment separately. Seller reserves the right to stop delivery of Products in transit and to withhold shipments in whole or in part if Buyer fails to make any payment to Seller when due or otherwise fails to perform its obligations hereunder. All shipping dates are approximate only, and Seller will not be liable for any loss or damage resulting from any delay in delivery or failure to deliver which is due to any cause beyond Seller's

reasonable control. In the event of a delay due to any cause beyond Seller's reasonable control, Seller reserves the right to terminate the order or to reschedule the shipment within a reasonable period of time, and Buyer will not be entitled to refuse delivery or otherwise be relieved of any obligations as the result of such delay. Products as to which delivery is delayed due to any cause within Buyer's control may be placed in storage by Seller at Buyer's risk and expense and for Buyer's account. Orders in process may be canceled by Buyer only with Seller's written consent and upon payment of Seller's cancellation charges. Orders in process may not be changed by Buyer except with Seller's written consent and upon agreement by the parties as to an appropriate adjustment in the purchase price therefore. Credit will not be allowed for Products returned without the prior written consent of Seller.

TITLE AND RISK OF LOSS: Notwithstanding the trade terms indicated above and subject to Seller's right to stop delivery of Products in transit, title to and risk of loss of the Products will pass to Buyer upon delivery of possession of the Products by Seller to the carrier; provided, however, that title to any software incorporated within or forming a part of the Products shall at all times remain with Seller or the licensor(s) thereof, as the case may be.

WARRANTY: Seller warrants that the Products will operate or perform substantially in conformance with Seller's published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended usage by properly trained personnel, for the period of time set forth in the product documentation, published specifications or package inserts. If a period of time is not specified in Seller's product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of shipment to Buyer for equipment and ninety (90) days for all other products (the "Warranty Period"). Seller agrees during the Warranty Period, to repair or replace, at Seller's option, defective Products so as to cause the same to operate in substantial conformance with said published specifications; provided that Buyer shall (a) promptly notify Seller in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after Seller's review, Seller will provide Buyer with service data and/or a Return Material Authorization ("RMA"), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to Seller with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of Seller. All replaced parts shall become the property of Seller. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Seller's Terms and Conditions of Sale.

Consumables are expressly excluded from this warranty. If Seller elects to repair defective medical device instruments, Seller may, in its sole discretion, provide a replacement loaner instrument to Buyer as necessary for use while the instruments are being repaired.

Notwithstanding the foregoing, Products supplied by Seller that are obtained by Seller from an original manufacturer or third party supplier are not warranted by Seller, but Seller agrees to assign to Buyer any warranty rights in such Product that Seller may have from the original manufacturer or third party supplier, to the extent such assignment is allowed by such original manufacturer or third party supplier.

In no event shall Seller have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Products in a manner for which they were not designed, (v) causes external to the Products such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Products or (vii) use of the Products in combination with equipment or software not supplied by Seller. If Seller determines that Products for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse Seller for all costs of investigating and responding to such request at Seller's then prevailing time and materials rates. If Seller provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay Seller therefor at Seller's then prevailing time and materials rates. ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, RELOCATION OR ALTERATION TO OR

TERMS AND CONDITIONS OF SALE

OF, OR OTHER TAMPERING WITH, THE PRODUCTS PERFORMED BY ANY PERSON OR ENTITY OTHER THAN SELLER WITHOUT SELLER'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY SELLER, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE AFFECTED PRODUCTS.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, SELLER DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCTS, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. SELLER DOES NOT WARRANT THAT THE PRODUCTS ARE ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

INDEMNIFICATION BY SELLER: Seller agrees to indemnify, defend and save Buyer, its officers, directors, and employees from and against any and all damages, liabilities, actions, causes of action, suits, claims, demands, losses, costs and expenses (including without limitation reasonable attorney's fees) ("Indemnified Items") for (i) injury to or death of persons or damage to property to the extent caused by the negligence or willful misconduct of Seller, its employees, agents or representatives or contractors in connection with the performance of services at Buyer's premises under this Agreement and (ii) claims that a Product infringes any valid United States patent, copyright or trade secret; provided, however, Seller shall have no liability under this Section to the extent any such Indemnified Items are caused by either (i) the negligence or willful misconduct of Buyer, its employees, agents or representatives or contractors, (ii) by any third party, (iii) use of a Product in combination with equipment or software not supplied by Seller where the Product would not itself be infringing, (iv) compliance with Buyer's designs, specifications or instructions, (v) use of the Product in an application or environment for which it was not designed or (vi) modifications of the Product by anyone other than Seller without Seller's prior written approval. Buyer shall provide Seller prompt written notice of any third party claim covered by Seller's indemnification obligations hereunder. Seller shall have the right to assume exclusive control of the defense of such claim or, at the option of the Seller, to settle the same. Buyer agrees to cooperate reasonably with the Seller in connection with the performance by Seller of its obligations in this Section.

Notwithstanding the above, Seller's infringement related indemnification obligations shall be extinguished and relieved if Seller, at its discretion and at its own expense (a) procures for Buyer the right, at no additional expense to Buyer, to continue using the Product; (b) replaces or modifies the Product so that it becomes non-infringing, provided the modification or replacement does not adversely affect the specifications of the Product; or (c) in the event (a) and (b) are not practical, refund to Buyer the amortized amounts paid by Buyer with respect thereto, based on a five (5) year amortization schedule. The foregoing indemnification provision states Seller's entire liability to Buyer for the claims described herein.

INDEMNIFICATION BY BUYER: Seller agrees to indemnify, defend and save Buyer, its officer, directors, and employees from and against any and all damages, liabilities, actions, causes of action, suits, claims, demands, losses, costs and expenses (including without limitation reasonable attorney's fees) ("Indemnified Items") for (i) injury to or death of persons or damage to property to the extent caused by the negligence or willful misconduct of Seller, its employees, agents or representatives or contractors in connection with the performance of services at Buyer's premises under this Agreement and (ii) claims that a Product infringes any valid United States patent, copyright or trade secret; provided, however, Seller shall have no liability under this Section to the extent any such Indemnified Items are caused by either (i) the negligence or willful misconduct of Buyer, its employees, agents or representatives or contractors, (ii) by any third party, (iii) use of a Product in combination with equipment or software not supplied by Seller where the Product would not itself be infringing, (iv) compliance with Buyer's designs, specifications or instructions, (v) use of the Product in an application or environment for which it was not designed or (vi) modifications of the Product by anyone other than Seller without Seller's prior written approval. Buyer shall provide Seller prompt written notice of any third party claim covered by Seller's indemnification obligations hereunder. Seller shall have the right to assume exclusive control of the defense of such claim or, at the

option of the Seller, to settle the same. Buyer agrees to cooperate reasonably with Seller in connection with the performance by Seller of its obligations in this Section.

Notwithstanding the above, Seller's infringement related indemnification obligations shall be extinguished and relieved if Seller, at its discretion and at its own expense (a) procures for Buyer the right, at no additional expense to Buyer, to continue using the Product; (b) replaces or modifies the Product so that it becomes non-infringing, provided the modification or replacement does not adversely affect the specifications of the Product; or (c) in the event (a) and (b) are not practical, refund to Buyer the amortized amounts paid by Buyer with respect thereto, based on a five (5) year amortization schedule. THE FOREGOING INDEMNIFICATION PROVISION STATES SELLER'S ENTIRE LIABILITY TO BUYER FOR THE CLAIMS DESCRIBED HEREIN.

Buyer shall indemnify, defend with competent and experienced counsel and hold harmless Seller, its parent, subsidiaries, affiliates and divisions, and their respective officers, directors, shareholders and employees, from and against any and all damages, liabilities, actions, causes of action, suits, claims, demands, losses, costs and expenses (including without limitation reasonable attorneys' fees and disbursements and court costs) to the extent arising from or in connection with (i) the negligence or willful misconduct of Buyer, its agents, employees, representatives or contractors; (ii) use of a Product in combination with equipment or software not supplied by Seller where the Product itself would not be infringing; (iii) Seller's compliance with designs, specifications or instructions supplied to Seller by Buyer; (iv) use of a Product in an application or environment for which it was not designed; or (v) modifications of a Product by anyone other than Seller without Seller's prior written approval.

SOFTWARE: With respect to any software products incorporated in or forming a part of the Products hereunder, Seller and Buyer intend and agree that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "Buyer" or similar or derivative words are understood and agreed to mean "licensee". Notwithstanding anything to the contrary contained herein, Seller or its licensor, as the case may be, retains all rights and interest in software products provided hereunder.

Seller hereby grants to Buyer a royalty-free, non-exclusive, nontransferable license, without power to sublicense, to use software provided hereunder solely for Buyer's own internal business purposes on the hardware products provided hereunder and to use the related documentation solely for Buyer's own internal business purposes. This license terminates when Buyer's lawful possession of the hardware products provided hereunder ceases, unless earlier terminated as provided herein. Buyer agrees to hold in confidence and not to sell, transfer, license, loan or otherwise make available in any form to third parties the software products and related documentation provided hereunder. Buyer may not disassemble, decompile or reverse engineer, copy, modify, enhance or otherwise change or supplement the software products provided hereunder without Seller's prior written consent. Seller will be entitled to terminate this license if Buyer fails to comply with any term or condition herein. Buyer agrees, upon termination of this license, immediately to return to Seller all software products and related documentation provided hereunder and all copies and portions thereof.

Certain of the software products provided by Seller may be owned by one or more third parties and licensed to Seller. Accordingly, Seller and Buyer agree that such third parties retain ownership of and title to such software products. The warranty and indemnification provisions set forth herein shall not apply to software products owned by third parties and provided hereunder.

LIMITATION OF LIABILITY: NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF SELLER UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF SELLER FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER SECTION 7 ABOVE)) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THERETOFORE PAID BY BUYER TO SELLER WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE MILLION DOLLARS (\$1,000,000). NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN,

TERMS AND CONDITIONS OF SALE

IN NO EVENT SHALL SELLER BE LIABLE FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF USE OF FACILITIES OR EQUIPMENT, LOSS OF REVENUE, LOSS OF DATA, LOSS OF PROFITS OR LOSS OF GOODWILL), REGARDLESS OF WHETHER SELLER (a) HAS BEEN INFORMED OF THE POSSIBILITY OF SUCH DAMAGES OR (b) IS NEGLIGENT.

EXPORT RESTRICTIONS: Buyer acknowledges that each Product and any related software and technology, including technical information supplied by Seller or contained in documents (collectively "Items"), is subject to export controls of the U.S. government. The export controls may include, but are not limited to, those of the Export Administration Regulations of the U.S. Department of Commerce (the "EAR"), which may restrict or require licenses for the export of Items from the United States and their re-export from other countries. Buyer shall comply with the EAR and all other applicable laws, regulations, laws, treaties, and agreements relating to the export, re-export, and import of any Item. Buyer shall not, without first obtaining the required license to do so from the appropriate U.S. government agency; (i) export or re-export any Item, or (ii) export, re-export, distribute or supply any Item to any restricted or embargoed country or to a person or entity whose privilege to participate in exports has been denied or restricted by the U.S. government. Buyer shall cooperate fully with Seller in any official or unofficial audit or inspection related to applicable export or import control laws or regulations, and shall indemnify and hold Seller harmless from, or in connection with, any violation of this Section by Buyer or its employees, consultants, agents, or customers.

PROMOTION CODE TERMS AND CONDITIONS: Discounts are valid on the purchase of list-price items only through fishersci.com orders. One-time use only. All promotions and discounts are on Fisher Scientific cart items only. Valid in the U.S. only. To redeem, enter the promo code in the box marked promo code at checkout. Free shipping offers apply to standard shipping only. Not valid with any other promotion. Terms are subject to change. Offer void where prohibited, licensed or restricted by federal, state, provincial or local laws or regulation or agency/institutional policy. Other restrictions may apply.

ACCEPTABLE PAYMENT METHODS: The Seller prefers to receive payment via ACH or other electronic interface methods that directly exchange funds between the Buyer's and Seller's bank accounts. The Seller also accepts checks mailed to one of its lockbox remittance locations. Although the Seller does accept credit card payments at the time of purchase, it does not accept credit card payments after the point of sale.

PROMOTIONAL TERMS FOR HEALTHCARE PROVIDERS AND FEDERAL, STATE AND LOCAL GOVERNMENT CUSTOMERS: Many promotions are not available to our current or potential healthcare customers or to federal, state and local government customers.

Healthcare customers include any healthcare practitioner, other healthcare provider or any individual or organization authorized to prescribe, dispense, purchase or influence the acquisition or use of medical devices or supplies for clinical use. Healthcare customers include, without limitation, any "Health Care Practitioner" under the Massachusetts Marketing Code of Conduct Law (defined at 105 C.M.R. §907.004) or a "Health Care Provider" under the Vermont Prescribed Products Gift Ban and Disclosure Law (defined at Vt. Stat. Ann. tit. 18, §4631a(8)). By participating in any promotion that is subject to the above limitations, Buyer warrants that it is not a current or potential healthcare customer or a federal, state and local government customer. Seller reserves the right to determine Buyer's status in the event of a dispute.

In the event that the promotion is available to Healthcare providers, Buyer acknowledges that the offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by Buyer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b) (3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

MISCELLANEOUS: (a) Buyer may not delegate any duties nor assign any rights or claims hereunder without Seller's prior written consent, and any such attempted delegation or assignment shall be void. (b) The rights and obligations of the parties hereunder shall be governed by and construed in accordance with the laws of the Commonwealth of Pennsylvania, without reference to its choice of law provisions. Each party hereby

irrevocably consents to the exclusive jurisdiction of the state and federal courts located in Allegheny County, Pennsylvania, USA, in any action arising out of or relating to this Agreement and waives any other venue to which it may be entitled by domicile or otherwise. (c) In the event of any legal proceeding between the Seller and Buyer relating to this Agreement, neither party may claim the right to a trial by jury, and both parties waive any right they may have under applicable law or otherwise to a right to a trial by jury. Any action arising under this Agreement must be brought within one (1) year from the date that the cause of action arose. (d) The application to this Agreement of the U.N. Convention on Contracts for the International Sale of Goods is hereby expressly excluded. (e) In the event that any one or more provisions contained herein shall be held by a court of competent jurisdiction to be invalid, illegal or unenforceable in any respect, the validity, legality and enforceability of the remaining provisions contained herein shall remain in full force and effect, unless the revision materially changes the bargain. (f) Seller's failure to enforce, or Seller's waiver of a breach of, any provision contained herein shall not constitute a waiver of any other breach or of such provision. (g) Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals. (h) Seller shall cooperate with Buyer to monitor invoicing accuracy and will conduct appropriate internal reviews upon request. Reviews shall be limited to once per year and shall cover the prior twelve month period. Reasonable credit shall be given to undercharges and overcharges. If Buyer engages external consultants to review invoice accuracy, Seller reserves the right to charge a reasonable fee for support provided to such party. (i) Buyer agrees that all pricing, discounts and technical information that Seller provides to Buyer are the confidential and proprietary information of Seller. Buyer agrees to (1) keep such information confidential and not disclose such information to any third party, and (2) use such information solely for Buyer's internal purposes and in connection with the Products supplied hereunder. Nothing herein shall restrict the use of information available to the general public. (j) Any notice or communication required or permitted hereunder shall be in writing and shall be deemed received when personally delivered or three (3) business days after being sent by certified mail, postage prepaid, to a party at the address specified herein or at such other address as either party may from time to time designate to the other. (k) Seller may, in its sole discretion, provide (1) applicable Product training to Buyer or its employees, or (2) samples of Products to Buyer for distribution to patients of Buyer. Buyer agrees that any such samples shall be distributed to patients for patient use or, if not so distributed, returned to Seller. Buyer shall not use such samples to provide care to patient and shall not bill patients or third party payers for the provision of such samples.

ACCEPTABLE PAYMENT METHODS: The Seller prefers to receive payment via ACH or other electronic interface methods that directly exchange funds between the Buyer's and Seller's bank accounts. The Seller also accepts checks mailed to one of its lockbox remittance locations. Although the Seller does accept credit card payments at the time of purchase, it does not accept credit card payments after the point of sale.

MEDICARE/MEDICAID REPORTING REQUIREMENTS: If Buyer is a recipient of Medicare/Medicaid funds, Buyer acknowledges that it has been informed of and agrees to fully and accurately account for, and report on its applicable cost report, the total value of any discount, rebate or other compensation paid hereunder in a way that complies with all applicable federal, state and local laws and regulations which establish "Safe Harbor" for discounts. Buyer shall make written request to Seller in the event Buyer requires additional information from Seller in order to meet its reporting requirements. Buyer acknowledges that agreement to such reporting requirement was a condition precedent to Seller's agreement to provide Products and that Seller would not have entered into this Agreement had Buyer not agreed to comply with such obligations.

PRODUCT DESCRIPTIONS: Product Descriptions are provided by the manufacturer of the product and are current as of the date of printing of this catalog. More current descriptions provided by the manufacturers are reflected at fisheredu.com.

LOOK TO FISHER SCIENCE EDUCATION FOR HASSLE-FREE SERVICE • PROMPT DELIVERY • GREAT PRICES • QUALITY PRODUCTS AND SELECTION

FISHER SCIENCE EDUCATION GUARANTEE

At Fisher Science Education, we are committed to providing you with the highest quality products and services possible. If, for any reason, you are not 100% satisfied with your purchase, simply contact our friendly customer service team within 30 days of the invoice date. We will exchange the product, repair or replace the product or issue you a credit for the full purchase price. Simply call us, toll-free, for a Return Authorization Number. Special order items, furniture, and closeouts cannot be exchanged or credited.

HOW TO ORDER

Fisher Science Education offers you many ways to order. So we can serve you quickly, it helps to have your account number available when placing your order. We gladly accept purchase orders issued on School or District letterhead. Registered users may also place orders online at fisheredu.com. For your convenience, we also accept VISA™, MasterCard™, and American Express™ credit cards. For credit card purchases, be ready to provide the card number, the type of card, exact name on the card, expiration date, full address and telephone number.

Mail: Fisher Science Education
4500 Turnberry Drive
Hanover Park, IL 60133

Phone: 1-800-955-1177 • 7 a.m. to 5:30 p.m.,
Central Time – Monday through Friday

Fax: 1-800-955-0740 • 24 hours a day, 7 days a week

Email: fse.custserv@thermofisher.com

Internet: fisheredu.com

Special note on international purchases: Please contact Global Export Customer Service if your order requires export and if the ultimate destination is outside of the United States.

Fisher Scientific Latin America and Global Export Customer Service Center:

3970 John's Creek Court, Suite 500, Suwanee, GA 30024 USA

Telephone (inside the U.S.): 1-800-955-5090

Telephone (outside the U.S.): 770-871-4725

Fax (outside the U.S.): 770-871-4726

PRICES AND PAYMENT

All pricing subject to change.

If you don't have an account and would like one, we'll be happy to set up a user account for you — simply call Customer Service at 1-800-955-1177.

To complete your credit card, check, or money order payments, please include applicable sales tax and 10% of the total price of all products (minimum of \$15) purchased to cover shipping, handling and insurance

Please provide payment (net) within 30 days from the date of the invoice.

If you have any questions about product selection, pricing, or invoicing, please call us at 1-800-955-1177.

SPECIAL PURCHASES

Fisher Science Education will sell BALANCES, CHEMICALS and BACTERIAL CULTURES to educational institutions ONLY. Please submit orders for these products on school letterhead or order them through your registered user account at fisheredu.com.

We will honor your requests for tax-exempt purchases if your account is coded tax-exempt (a signed tax-exempt certificate is on file with Fisher Science Education) or, if purchasing with a credit card, the card is an official tax-exempt card.

Some of the products we sell are time sensitive and must be delivered within a specified period of time. Others may be considered hazardous materials. For time-sensitive materials, please specify the desired delivery date when you want the products to arrive so that you will be available to receive them. There is a charge for time-sensitive and hazardous materials. (See: Shipping Your Order.)

SHIPPING YOUR ORDER

We want your purchase to arrive promptly and undamaged. To ensure this, we only deliver products to a street address, not to a Post Office Box — please remember this when placing your order. Other things to consider for product shipping include:

How should we ship your purchase? We will ship according to your special instructions, whenever possible. Otherwise, we will select the best shipment method, based on the nature of the material, destination, package size and weight.

Do you need extra-quick delivery? Overnight shipping is available for many products at an additional cost, but certain hazardous materials may be exempt from this service — please ask your Customer Service Representative when placing your order.

Is motor freight necessary? Some items must be shipped by motor freight (rather than by small parcel carriers) and will have different shipping costs. These items are noted in our catalog with a "Motor Freight" icon. Call your Customer Service Representative for assistance in calculating these expenses.

Is the material hazardous? There is a single charge of \$27.50 per order for hazardous materials, no matter how many hazardous items are on your order.

Is your purchase live or time-sensitive? Be sure to specify the desired delivery date for all live material orders. Some live products have minimum purchase requirements, require special handling fees or have shipping restrictions. Contact Customer Service for specific details.

ONCE YOU RECEIVE YOUR ORDER

Inspect your shipment! Open and inspect your order immediately upon arrival and notify us of any problems. Report claims for damage within the seven (7) day period after delivery to allow us to remedy problems with shipping companies.

Want to return a product? The dedicated staff at Fisher Science Education Customer Service is committed to resolving any questions or problems you may encounter. Here are the steps to take:

- Please notify us no later than seven (7) days after delivery of any damaged or unsatisfactory shipments
- Obtain a return authorization number from Fisher Science Education Customer Service prior to returning any product.
- Provide the Fisher Science Education Confirmation Number, found in the upper left hand corner of the packing slip or invoice, or your account number and the purchase order number.

NOTE: Products are subject to a fifteen percent (15%) restocking charge, unless they were not supplied in accordance with your order or are defective at the time of delivery.

Products in this catalog are intended for use by children age 13 or over.

LIMITED WARRANTY

Fisher Scientific Company L.L.C. ("Fisher") will make available the manufacturer's warranties, without recourse to Fisher.

However, except as otherwise set forth here, Fisher expressly disclaims all other warranties, whether expressed, implied or statutory, including the warranties of merchantability and fitness for a particular purpose.

In no event shall Fisher be liable for any special, indirect, incidental, consequential or exemplary damages, however arising, even if Fisher Scientific has been advised of the possibility of such damages, out of the sale or use of the products described in this catalog.

Your sole and exclusive remedy for any cause of action arising out of the sale or use of any product described in this catalog shall be repair or replacement of the product or refund of the purchase price at Fisher's option.



GOGGLE SANITIZING CABINET

Safe and easy to use.

Thoroughly sanitize up to 36 safety goggles and glasses in just 15 minutes with UV light!



UV Goggle Sanitizing Cabinet

Code	Description
S07817	UV Goggle Sanitizing Cabinet

fisheredu.com/eisco



Protector™ ClassMate™ Fume Hoods

360° visibility designed for instructional labs.

- High-performance hood
- Low profile airfoil with spill trough
- Angled sash allows closer, more comfortable viewing

fisheredu.com Search: ClassMate

LABCONCO



Safety

2022 Product Catalog

In the United States

Order online: fisheredu.com
Fax an order: 1-800-955-0740
Call customer service: 1-800-955-1177

Fisher Science Education

4500 Turnberry Drive
Hanover Park, IL 60133

