



Magne-Box™ Polycarbonate Cryo-Boxes

Patent pending magnetic closure design!

- Each box features magnets that securely connect base to hinged lid and ensure samples stay secure
- Polycarbonate design withstands temperatures from -196 to 121°C
- Boxes fit standard size freezer racks
- 25, 81, and 100-place boxes accommodate 1.5mL to 2mL cryogenic vials
- 81-tall boxes accommodate 4mL to 5mL cryovials



Ordering Information

Product	Dimensions	Color	Qty/PK	Cat. No.
25-place box	3 x 3 x 2.06"	Orange	8	03-391-290
81-place box	5.25 x 5.25 x 2.13"	Assorted	4	03-391-291
81-tall box	5.25 x 5.25 x 3.75"	Blue	4	03-391-292
100-place box	5.25 x 5.25 x 2.06"	Green	4	03-391-293

Find it at [fishersci.com](https://www.fishersci.com)

Contact us today:

In the United States

Order online: [fishersci.com](https://www.fishersci.com)
 Fax an order: 1-800-926-1166
 Call customer service: 1-800-766-7000

In Canada

Order online: [fishersci.ca](https://www.fishersci.ca)
 Fax an order: 1-800-463-2996
 Call customer service: 1-800-234-7437

