



Not all labware is created equal

Thermo Scientific Nalgene beakers and graduated cylinders are now available in convenient variety packs.

Get seven cylinders ranging in size from 10mL to 1L or seven beakers ranging from 50mL to 1L in each variety pack.

A convenient and cost effective way to get started setting up your new lab, assemble student labware kits for a new semester, or refresh an established lab. Break the habit of using glass labware and take advantage of the benefits of plastic: durability, light weight, and low cost.

Thermo Scientific™ Nalgene™ labware sets the standard for quality and superior performance for your valuable research:

- Safer alternative to glass – impact-resistant polypropylene copolymer material is shatter proof at temperatures above -40 °C
- Accurate to internationally recognized laboratory standards for precision work
- Made with high-purity materials to protect the integrity of your research solutions; suitable for food and beverage use
- Polypropylene copolymer material is compatible with a wide variety of laboratory reagents, including acids, bases and alcohols
- Beakers and cylinders are autoclavable. (*Note: Autoclaving will affect the accuracy of the graduations.*)

Great for general laboratory use, strong and durable



Includes one each 50ml, 100ml, 150ml, 250ml, 400ml, 600ml, and 1000ml beakers

- Complies with 21 CFR Part 177.1520 and ISO/DIN 7056 International Standards for Laboratory Plasticware Beakers
- No-drip spout design provides precision pouring control
- Ribbed and nests for easy stacking
- Light weight for easy handling – feel the difference

► **Beaker variety pack - Cat# 02-540-271**



Includes one each 10ml, 25ml, 50ml, 100ml, 250ml, 500ml, and 1000mL graduated cylinders

- Complies with 21 CFR Part 177.1520 and ASTM E1272 Class B requirements as well as all requirements of ISO Standard 6706
- Permanently molded-in graduations are inked for easy reading
- No meniscus simplifies reading accuracy
- Large blue octagonal base prevents rolling and tipping

► **Cylinder variety pack - Cat# 02-540-270**

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.



In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

PFCCGBKCYLPCK 0614

In Canada:

For customer service, call 1-800-234-7437
To fax an order, use 1-800-463-2996
To order online: www.fishersci.ca

BN0609143

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand