

**Chemicals** 

## Introducing the Elite 8 The Right Solvents for Your Applications

The Elite 8 Solvents from Thermo Fisher Scientific include the most commonly used stockroom, analytical laboratory, and production solvents. These solvents are available in a variety of grades and pack sizes, allowing you to find the best suited option for your application. No matter the application, we're dedicated to providing you with the quality you desire, the purity you require, and the certainty you deserve.

thermo scientific

> A458-1 1 L LOT 195943 Methanol

## Find the right grade of solvent to use with your analytical instruments

	Solvent Grade						
Solvent	Thermo Scientific™ UHPLC-MS	Fisher Chemical™ Optima™ LC-MS	Fisher Chemical GC Resolve	Fisher Chemical GC Headspace	Fisher Chemical Optima HPLC	Fisher Chemical Pesticide	Fisher Chemical HPLC
	Instrument						
	UHPLC-MS, UHPLC-UV	UHPLC-UV, LCMS	GC		GC, HPLC	GC, LC	HPLC
Acetone	-	-	A928-4	-	A929-4	A40-4	A949-4
Acetonitrile	A956-1	A955-4	-	-	A996-4	-	A998-4
Ethyl Acetate	-	-	-	-	E196-4	E1914	E195-4
Hexanes -60% n-hexane	-	-	H307-4	-	H303-4	H300-4	H302-4
Isopropanol	-	A461-4	-	-	A464-4	A519-4	A451-4
Methanol	A458-1	A456-4	A457-4	-	A454-4	-	A452-4
Methylene Chloride	-	_	D154-4	-	D151-4	D142-4	D143-4
Water	W8-1	W6-4	-	W10-1	W7-4	-	W5-4

Fisher Chemical GC Resolv Solvents are the premium option for pesticide analysis.

## Learn more at fishersci.com/chemicals and fishersci.ca/chemicals or fishersci.com/ThermoScientificChemicals and fishersci.ca/ThermoScientificChemicals

Distributed by Fisher Scientific. Contact us today:

In the United States Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

## In Canada

Order online: fishersci.ca Fax an order: 1-800-463-2996 Call customer service: 1-800-234-7437



© 2022 Thermo Fisher Scientific Inc. All rights reserved. Trademarks used are owned as indicated at fishersci.com/trademarks.