

# Complementary Solutions for Immunohistochemistry

Optimize your immunohistochemistry panel offerings! When it comes to performing IHC in your laboratory, it's vital that you have all the tools necessary for accurate and reproducible results. Whether research or clinical applications, Epredia has recommended products to compliment your assays in a variety of sizes to meet the specific needs of your laboratory.

## Research Ancillaries

Research Detection and Ancillaries			
Pretreatment	Volume	Catalog No.	Description
Dewax & HIER Buffer L	10 x 100 mL	TA-999-DHBL	All-in-one dewaxing and pH 6 citrate-based epitope retrieval buffer
Dewax & HIER Buffer M	10 x 100 mL	TA-999-DHBM	All-in-one dewaxing and pH 8 EDTA-based epitope retrieval buffer
Dewax & HIER Buffer H	10 x 100 mL	TA-999-DHBH	All-in-one dewaxing and pH 9 Tris-EDTA-based epitope retrieval buffer

Detection Kit	Volume	Catalog No.	Description
Quanto Detection Kit Mouse on Mouse (RUO)	60 mL	TL-060-QHDM	Peroxide Block, UV-Block, Amplifier, HRP Polymer, DAB Quanto
DAB Quanto	60 mL	TL-125-QHD	UltraVision Quanto Detection Kit HRP DAB

Ancillary Reagents	Volume	Catalog No.	Description
Mayer's Hematoxylin	60 mL	TA-060-MH	IHC Counterstain
Mayer's Hematoxylin	125 mL	TA-125-MH	IHC Counterstain
Antibody Diluent OP Quanto	125 mL	TA-125-ADQ	Antibody diluent for concentrates primary antibodies
Phosphate Buffered Saline and Tween 20	1000 mL	TA-999-PT	Rinse buffer for Autostainer
Tris Buffered Saline and Tween 20	1000 mL	TA-999-TT	Rinse buffer for Autostainer


## Clinical Ancillaries

### Clinical Detection and Ancillaries: IVD

Pretreatment	Volume	Catalog No.	Description
Dewax & HIER Buffer L	10 x 100 mL	TA-999-DHBL	All-in-one dewaxing and pH 6 citrate-based epitope retrieval buffer
Dewax & HIER Buffer M	10 x 100 mL	TA-999-DHBM	All-in-one dewaxing and pH 8 EDTA-based epitope retrieval buffer
Dewax & HIER Buffer H	10 x 100 mL	TA-999-DHBH	All-in-one dewaxing and pH 9 Tris-EDTA-based epitope retrieval buffer

Detection Kit	Volume	Catalog No.	Description
Quanto Detection Kit HRP DAB	60 mL	TL-060-QHD	Peroxide Block, UV-Block, Amplifier, HRP Polymer, DAB Quanto
Quanto Detection Kit HRP DAB	125 mL	TL-125-QHD	Peroxide Block, UV-Block, Amplifier, HRP Polymer, DAB Quanto

Ancillary Reagents	Volume	Catalog No.	Description
Mayer's Hematoxylin	60 mL	TA-060-MH	IHC Counterstain
Mayer's Hematoxylin	125 mL	TA-125-MH	IHC Counterstain
Antibody Diluent OP Quanto	125 mL	TA-125-ADQ	Antibody diluent for concentrates primary antibodies
Phosphate Buffered Saline and Tween 20	1000 mL	TA-999-PT	Rinse buffer for Autostainer
Tris Buffered Saline and Tween 20	1000 mL	TA-999-TT	Rinse buffer for Autostainer



## Clinical Primaries

Breast Panel: IVD			
Antibody	Clone	Catalog No.	Species
Estrogen Receptor	SP1	RM-9101	Rabbit monoclonal
Progesterone Receptor	SP2	RM-9102	Rabbit monoclonal
Her2/neu	SP3	RM-9103	Rabbit monoclonal
Ki67	SP6	RM-9106	Rabbit monoclonal
p53	Y5	RM-2103	Rabbit monoclonal
Keratin 7	OV-TL 12/30	MS-1352	Mouse monoclonal
Keratin 20	Ks20.8	MS-377	Mouse monoclonal

Prostate Panel: IVD			
Antibody	Clone	Catalog No.	Species
p504S	13H4	RM-9130	Rabbit monoclonal
PSA	EP1588Y	RM-2104	Rabbit monoclonal
PSAP	PASE/4LJ	MS-321	Mouse monoclonal
Ki67	SP6	RM-9106	Rabbit monoclonal
Keratin 5	EP1601Y	RM-2106	Rabbit monoclonal
Keratin 14	Polyclonal	RB-9020	Rabbit polyclonal
Keratin 20	Ks20.8	MS-377	Mouse monoclonal

Gastrointestinal Panel: IVD			
Antibody	Clone	Catalog No.	Species
Dog-1	SP31	RM-9132	Rabbit monoclonal
CD117	Polyclonal	RB-9038	Rabbit polyclonal
Helicobacter pylori	Polyclonal	RB-9070	Rabbit polyclonal
Ki67	SP6	RM-9106	Rabbit monoclonal
Pan Cytokeratin	AE1/AE3	MS-343	Mouse monoclonal
Keratin 7	OV-TL 12/30	MS-1352	Mouse monoclonal
Keratin 20	Ks20.8	MS-377	Mouse monoclonal

Distributed by Fisher Healthcare. Contact us today:

### In the United States

Order online: [fisherhealthcare.com](http://fisherhealthcare.com)

Call customer service: 1-800-640-0640

