




# Organic Synthesis

Products for Simple and Sophisticated Workflows

With more than 60 years of synthesis experience, TCI specializes in organic reagents and products that are hard to find and difficult to synthesize.



## Organic Building Blocks


From amylamine to xanthurenic acid, TCI has what you need to keep your research moving forward.

Cat. No.	Description	Quantity
A002425ML	Butyl Acetate, 99.0+%	25 mL
A002625ML	tert-Butyl Acetate 99.0+%	25 mL
A0082100G	Acetyl Chloride 98.0+%	100 g
A064825G	2-Aminoethanethiol 95.0+%	25 g
A014225ML	Butyl Acrylate (Stabilized with MEHQ), 99.0+%	25 mL
A014325ML	Ethyl Acrylate (Stabilized with MEHQ), 99.0+%	25 mL
A014525ML	Methyl Acrylate (Stabilized with MEHQ), 99.0+%	25 mL
A016125G	Adipic Acid, 99.0+%	25 g
A026925G	4-Aminobenzoic Acid, 99.0+%	25 g
A048025ML	p-Anisaldehyde, 99.0+%	25 mL

Visit [fishersci.com/TCI-chemical-building-blocks](https://fishersci.com/TCI-chemical-building-blocks) or [fishersci.ca/TCI-chemical-building-blocks](https://fishersci.ca/TCI-chemical-building-blocks) to learn more.

## Ligands for Research and Development

From general use to focused carbene ligands with no salt byproducts, find a variety of ligands for your lab.



Cat. No.	Description	Quantity
D013425G	1,4-Diazabicyclo[2.2.2]octane, 98.0+%	25 g
D044325ML	Dicyclopentadiene (Stabilized with BHT), 97.0+%	25 mL
E007725ML	Ethylenediamine Anhydrous, 98.0+%	25 mL
H009325G	Hexamethylenetetramine, 99.0+%	25 g
I000125G	Imidazole, 98.0+%	25 g
L006725ML	2,6-Lutidine, 98.0+%	25 mL
P005225ML	Acetylacetone, 99.0+%	25 mL
T014725ML	N,N,N',N'-Tetramethylethylenediamine, 98.0+%	25 mL
T047025ML	1,3,5-Trimethylbenzene, 97.0+%	25 mL
T29931G	Tris[(1-benzyl-1H-1,2,3-triazol-4-yl)methyl]amine, 97.0+%	1 g

Visit [fishersci.com/tci-chemical-ligands](https://fishersci.com/tci-chemical-ligands) or [fishersci.ca/tci-chemical-ligands](https://fishersci.ca/tci-chemical-ligands) to browse hundreds of options.

Visit [fishersci.com/organicsynthesis](https://fishersci.com/organicsynthesis) or [fishersci.ca/organicsynthesis](https://fishersci.ca/organicsynthesis) to learn more.

## Matrix Materials for MALDI

TCI offers various matrices for Matrix Assisted Laser Desorption Ionization (MALDI). To see which analytes are suitable for your applications, please refer to the selection table below.

Cat. No.	Synonym	Proteins	Peptides	Glycans	Lipids	Nucleic Acids	Synthetic Polymers	Other Methods
D2932	SA	✓						
C0002	CA	✓						
F1257	FA	✓						
N1103	NA	✓						
M0247	MBT		✓					
C1297	CMBT	✓	✓	✓				
C1768	α-CHCA	✓	✓	✓				PSD
D0561	2,4-DHAP	✓	✓	✓		✓		ISD
D1955	2,5-DHAP	✓	✓	✓		✓		
D1716	2,6-DHAP	✓	✓	✓		✓		
D2933	DHBA		✓	✓	✓			ISD
D5355	DHB/B		✓	✓	✓			
T1888	THAP		✓	✓	✓	✓		PSD
A0425	AMT			✓				
A0418	3-AQ			✓				
A0859	AHBA			✓		✓		
D5689	1,5-DAN		✓		✓			ISD
B1018	DMAN				✓			
P0421	PA					✓		
H0787	3-HPA					✓		
D1950	DABP					✓		
E1287	Esculetin					✓		
H1622	HABA		✓	✓			✓	
P1121	Norharman	✓	✓	✓	✓	✓	✓	
H1237	Harmaline	✓	✓	✓	✓	✓	✓	
I0025	IAA						✓	
T1166	DIT						✓	
T0168	TPB						✓	
B3635	DCTB						✓	
C0287	CSA						✓	
R0064	ATRA						✓	
A3299	9-ACA						✓	
N0987	9-NA						✓	
A0317	5-ASA		✓					ISD
A0358	5,1-ANL	✓	✓					ISD
H0635	3H4NBA		✓	✓				ISD, NBS
A3298	9-AA				✓			MALDI Imaging

Visit [fishersci.com/organicsynthesis](https://fishersci.com/organicsynthesis) or [fishersci.ca/organicsynthesis](https://fishersci.ca/organicsynthesis) to learn more.

Distributed by Fisher Scientific. Contact us today:

### In the United States

Order online: [fishersci.com](https://fishersci.com)

Call customer service: 1-800-766-7000

### In Canada

Order online: [fishersci.ca](https://fishersci.ca)

Call customer service: 1-800-234-7437

